

A novel targeted NGS-based assay for simultaneous detection of SNVs and structural variation in pharmacogenes

Introduction

Pharmacogenomics (PGx) can help optimizing drug therapy, yet current methods have limitations in terms of genomic coverage, costs, incidental findings and privacy concerns. Furthermore, several technologies struggle with accurate genotyping of complex but important genes like *CYP2D6* due to frequent copy number variations (CNVs) and hybrid gene formations.

To address the need for a simple, single-assay solution for detecting both single nucleotide variants (SNVs) and complex structural variations in pharmacogenes, we are developing two new digitalMLPA assays:

- 2015 Pharmacogenomics Panel 1: targets 29 actionable PGx genes, containing 205 probes for copy number determination and 217 SNV-specific probes, including all AMP Tier 1 and most Tier 2 variants.
- 2016 Pharmacogenomics Panel 2: contains all 2015 probes, plus additional probes for 52 PGx related genes (adding to a total of 436 probes for CN determination and 315 SNV-specific probes).

In this study, we evaluate the performance of these assays using reference materials to assess their ability to detect SNVs and structural variants in prominent pharmacogenes.

Performance evaluation of digitalMLPA pharmacogenomics assays

The digitalMLPA pharmacogenomics assays were tested on ~350 DNA samples with known genotype information based on GeT-RM publications, and Ensembl for variants without published reference materials.

The assays showed high concordance with known data (see Table 1)(see **Table 1**). Note that only a subset of tested genes is displayed. In various samples, the results did not match the GeT-RM result but did match the genotype information in Ensembl.

The number of probes that did not show concordance in all samples was limited and these are currently investigated to improve specificity.

For *CYP2D6*, discordance in 3 samples was due to involvement of hybrid alleles that were not mentioned in GeT-RM call.

Gene	Number of probes in 2016 Pharmacogenomics Panel 2		GeT-RM		Ensembl (n=152)
	# CNV probes	# SNV probes	# Samples	% Concordance	
<i>CYP2B6</i>	8	8	95	99	100
<i>CYP2C9</i>	5	12	115	99	100
<i>CYP2C19</i>	6	9	116	100	98
<i>CYP2D6</i>	25	22	120	97	99
<i>CYP3A4</i>	3	9	110	100	100
<i>CYP3A5</i>	3	3	95	100	100
<i>CYP4F2</i>	5	3	68	97	99
<i>DPYD</i>	28	18	59	100	100
<i>G6PD</i>	4	22	-	-	100
<i>HLA-A *31:01</i>	0	1	76	100	-
<i>HLA-B *15:02</i>	0	1	76	100	-
<i>HLA-B *57:01</i>	0	1	76	100	-
<i>HLA-B *58:01</i>	0	1	73	100	-
<i>NUOT15</i>	2	5	13	96	96
<i>RYR1</i>	0	23	-	-	100
<i>SLCO1B1</i>	10	8	95	92	100
<i>TPMT</i>	9	12	102	100	100
<i>UGT1A1</i>	6	6	78	100	100
<i>VKORC1</i>	3	4	22	100	100
Other genes	316	142	-	-	>99

Table 1. Concordance between output of the digitalMLPA 2015/2016 pharmacogenomics assays and published genotype data (GeT-RM and Ensembl), per gene. The number of investigated samples is indicated.

digitalMLPA: a high-throughput method for reliable SNV and CNV detection

digitalMLPA is a multiplex PCR-based method for detecting up to 1600 selected SNVs and structural variants in a single assay. It combines the sensitivity of CNV and SNV detection of the widely used MLPA technology with the broad scale of NGS for amplicon detection.

The 2015/2016 pharmacogenomics assays are compatible with:

- purified genomic DNA
- crude DNA extracts from whole blood without any purification
- crude DNA extracts from 3 mm punch of a DBS without any purification
- buccal swab samples without any purification

The workflow of digitalMLPA is illustrated in **Figure 1**.




Figure 1. Samples (containing 20-40 ng DNA) are first denatured, after which the left and right probe oligo bind to their target DNA during overnight hybridisation. This is followed by ligation of all hybridised probes. Successfully ligated probes are amplified by multiplex PCR. Amplicons of all reactions are combined and loaded on an Illumina sequencing platform (no library quantification is required). Data analysis is performed using the free analysis software Coffalyser digitalMLPA.

Tackling complexity: analysing structural rearrangements in *CYP2D6*




Figure 2. Schematic representation of the *CYP2D6* locus, its frequent structural variants and corresponding star alleles. Adapted from Turner et al., 2023. Clin Pharmacol Ther.

Gene	Probe Location	Probe Target	Associated Star Allele	*1/*1	*2/*5	*5/*5	*1x2/*2x2	*1/*13+*2 intron 1 switch	*1/*13+*2 intron 4 switch	*39/*13 exon 7-8 switch	*1/*83	*36+*10/ *36+*10	*1/*36+*10	*1/*68+*4
CYP2D7	Exon 1	Copy number	*13, *36, *68, *83 (CYP2D6)	1.9	2.2	1.8	1.8	2.0	1.8	1.8	1.8	2.1	2.0	2.1
	Intron 4			1.8	2.1	1.8	1.9	0.9	2.0	1.9	1.6	2.3	2.0	3.1
	Exon 9			1.8	2.0	1.7	2.0	1.0	0.9	0.9	1.7	4.3	3.1	2.8
	Exon 9			2.1	2.1	1.9	1.9	1.0	1.0	1.0	1.9	3.8	3.0	3.0
	Exon 9			1.9	2.0	1.8	2.0	1.0	0.9	2.0	3.9	2.9	2.9	
	4500 nt upstream	2.0		1.0	0.0	3.8	2.0	1.9	1.0	1.8	3.7	2.8	3.0	
	2500 nt upstream	2.1		1.0	0.0	3.8	2.1	1.8	0.9	1.9	3.9	3.0	3.0	
	2000 nt upstream	2.0		1.1	0.0	4.0	2.2	2.0	1.0	2.1	3.8	2.9	2.8	
	1500 nt upstream	2.0		1.1	0.0	3.9	2.1	1.9	0.9	2.0	4.0	3.1	3.0	
	250 nt upstream	2.1		1.1	0.1	3.8	2.0	1.9	1.0	1.9	3.9	2.9	2.9	
	Exon 1	2.0		1.0	0.0	3.9	2.1	1.9	1.0	2.0	3.9	3.0	3.2	
	Intron 1	2.0		1.0	0.0	4.0	2.0	2.1	0.9	2.0	3.9	3.0	2.8	
	Intron 2	2.1		1.1	0.0	4.0	3.0	1.9	1.0	2.0	4.0	3.1	2.1	
	Intron 2	1.9		1.1	0.0	3.8	2.9	1.9	0.9	1.9	3.7	2.9	1.9	
CYP2D6	Intron 2	Copy number	*N, *5, *13, *36, *68, *83	2.1	1.0	0.0	3.8	3.2	2.0	0.9	2.1	3.8	2.9	1.7
	Intron 2			2.0	1.1	0.0	3.9	3.1	1.9	0.9	2.1	3.8	3.0	2.0
	Intron 4			2.1	1.0	0.0	3.9	3.0	2.0	1.0	2.0	3.9	2.9	1.9
	Intron 4			2.0	1.0	0.0	3.9	3.1	2.9	0.9	2.0	3.8	2.9	2.0
	Intron 5			2.0	1.0	0.0	4.0	3.1	2.8	0.9	1.9	3.8	2.9	2.0
	Intron 6			2.0	1.0	0.0	4.0	3.1	2.9	0.9	2.0	4.0	3.0	2.0
	Exon 9			2.0	1.0	0.0	3.9	2.9	2.9	1.8	1.0	2.1	1.8	1.8
	Exon 9			1.8	0.9	0.0	3.7	2.9	2.9	1.9	1.0	1.8	1.8	2.0
	150 nt downstream			2.0	1.1	0.0	3.8	3.1	2.9	2.0	2.0	2.0	2.1	1.9
	450 nt downstream			2.2	1.1	0.0	4.3	3.0	3.0	2.1	2.0	2.0	2.0	2.2
	500 nt downstream			2.1	0.9	0.0	4.0	3.0	2.9	2.0	1.9	2.0	2.0	2.0
	Exon 1	rs5030862		*12 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
	Exon 1	rs1065852		*10 (AMP Tier 1)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	3.8	2.0	2.0
	Intron 1	rs201377835		*11	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Exon 2	rs28371706	*17 (AMP Tier 1)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0		
Exon 9	rs1135840	*2 (AMP Tier 1)	0.0	1.0	0.0	2.1	2.1	1.9	1.0	0.0	1.9	0.9	1.0	
Exon 3	rs1135822	*40 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 3	rs61736512	*29 (AMP Tier 1)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 3	rs5030655	*6 (AMP Tier 1)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 3	rs5030685	*14 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 3	rs5030685	*8 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Intron 3	rs3892097	*4 (AMP Tier 1)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 4	rs72549356	*40 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 5	rs35742686	*3 (AMP Tier 1)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 5	rs5030656	*9 (AMP Tier 1)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 6	rs79292917	*59 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 6	rs5030867	*7 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	
Exon 6	rs16947	*21*17 (AMP Tier 1)	0.0	1.1	0.0	2.1	2.1	1.8	0.0	0.0	0.0	0.0	0.0	0.0
Intron 6	rs28371725	*41 (AMP Tier 1)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Exon 7	rs58421388	*29 (AMP Tier 1)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Exon 7	rs72549347	*56 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Exon 7	rs72549346	*42 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Exon 9	rs267608319	*31 (AMP Tier 2)	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0

Figure 3. Heatmap of probes for CN determination (upper part) and SNVs (lower part) in *CYP2D6* and *CYP2D7*, tested on various samples with known structural variations using the 2015/2016 pharmacogenomics assays. Diplotypes of samples are displayed in upper row. Values represent the calculated number of copies.

CYP2D6 is a key gene in drug metabolism but challenging to characterise due to frequent CNVs and hybrid alleles (see **Figure 2**). digitalMLPA probes were designed to specifically discriminate between *CYP2D6* and *CYP2D7*. Various samples with known *CYP2D6* deletions and hybrid alleles were tested using the 2015/2016 assays. The correct copy number status as well as the approximate breakpoint in hybrid alleles was identified (**Figure 3**), indicating the value of these digitalMLPA assays to detect structural variants in *CYP2D6*.

For further information please contact
Paul van Vught at paul@mrcholland.com

Interested in our products?



Scan the QR code to be involved
and receive the latest news!

Download this poster here:



Conclusions

- The **digitalMLPA pharmacogenomics assays show high concordance** to known datasets when tested on ~350 DNA samples with known variants, **including all AMP Tier 1 and most Tier 2 targets**.
- These assays can be used for the detection of both selected SNVs and (complex) structural variants in pharmacogenes, **including *CYP2D6***, in a single reaction.
- As **digitalMLPA is a targeted method**, where probes are sequenced instead of DNA, this greatly reduces incidental findings and privacy concerns, making this assay a promising technique for the characterisation of pharmacogenes.