

Product Description SALSA® MLPA® Probemix P216-B2 Growth Hormone Deficiency mix-1

To be used with the MLPA General Protocol.

Version B2. Compared to version B1, four reference probes have been replaced and three probe lengths have been adjusted. For complete product history see page 7.

Catalogue numbers:

- **P216-025R:** SALSA MLPA Probemix P216 Growth Hormone Deficiency mix-1, 25 reactions.
- **P216-050R:** SALSA MLPA Probemix P216 Growth Hormone Deficiency mix-1, 50 reactions.
- **P216-100R:** SALSA MLPA Probemix P216 Growth Hormone Deficiency mix-1, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mlpa.com).

Certificate of Analysis: Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mlpa.com.

Precautions and warnings: For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mlpa.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information: The SALSA MLPA Probemix P216 Growth Hormone Deficiency mix-1 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *GH1*, *POU1F1*, *PROP1*, *GHRHR*, *LHX3*, *LHX4*, and *HESX1* genes, which are associated with Growth Hormone Deficiency (GHD).

GHD is characterised by short stature, delayed growth velocity, and delayed skeletal maturation. Idiopathic growth hormone deficiency is the most common cause of GHD in children.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK1347/>.

This SALSA MLPA Probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering: The exon numbering used in this P216-B2 Growth Hormone Deficiency mix-1 product description is the exon numbering from the RefSeq transcript NM_033343.3 (*LHX4*), NM_003865.2 (*HESX1*), NM_001122757.2 (*POU1F1*), NM_006261.4 (*PROP1*), NM_000823.3 (*GHRHR*), NM_178138.5 (*LHX3*), and NM_000515.4 (*GH1*), which is identical to the NG_008081.1, NG_008242.1, NG_008225.2, NG_015889.1, NG_021416.1, NG_008097.1, and NG_011676.1 sequences respectively. The exon numbering and NM_ sequence used have been retrieved on 06/2019. As changes to the NCBI database can occur after release of this product description, exon numbering may not be up-to-date.

Probemix content: The SALSA MLPA Probemix P216-B2 Growth Hormone Deficiency mix-1 contains 49 MLPA probes with amplification products between 130 and 505 nucleotides (nt). This includes four probes for the *GH1* gene, five probes for the *POU1F1* gene, three probes for the *PROP1* gene, 12 probes for the *GHRHR* gene, seven probes for the *LHX3* gene, six probes for the *LHX4* gene, and finally four probes for the *HESX1* gene (Tables 2a-g). In addition, eight reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mlpa.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mlpa.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal of 88 nt and 96 nt fragment indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique: The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mlpa.com).

MLPA technique validation: Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens: Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples: A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from unrelated individuals who are from families without a history of growth hormone deficiency. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol.

Positive control DNA samples: MRC-Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/home.html>) have a diverse collection of biological resources which may be used as a positive control DNA sample in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

Data analysis: Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mlpa.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results: The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the dosage quotient (DQ) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the DQ of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Dosage quotient
Normal	$0.80 < DQ < 1.20$
Homozygous deletion	$DQ = 0$
Heterozygous deletion	$0.40 < DQ < 0.65$
Heterozygous duplication	$1.30 < DQ < 1.65$
Heterozygous triplication/Homozygous duplication	$1.75 < DQ < 2.15$
Ambiguous copy number	All other values

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Incomplete DNA denaturation (e.g. due to salt contamination) can lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *LHX3* gene. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- When running MLPA products, the capillary electrophoresis protocol may need optimization. False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: lower injection voltage / injection time settings, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure:

- In most populations, the major cause of genetic defects in the *GH1*, *POU1F1*, *PROPI*, *GHRHR*, *LHX3*, *LHX4*, and *HESX1* genes are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P216 Growth Hormone Deficiency mix-1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNPs, point mutations, small indels) in the target sequence detected by a probe can cause false positive results. Mutations/SNPs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results: Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

Mutation databases: For all genes present in this P216 probemix a database is available: <https://databases.lovd.nl/shared/genes>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNPs and unusual results (e.g., a duplication of *GHRHR* exons 5 and 7 but not exon 6) to MRC-Holland: info@mlpa.com.

Table 1. SALSA MLPA Probemix P216-B2 Growth Hormone Deficiency mix-1

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a	
		Reference	Target gene
64-105	Control fragments – see table in probemix content section for more information		
130 *	Reference probe 20108-L27306	8p23	
137 †	GHRHR probe 07905-L31750		Exon 2
142 ‹	LHX3 probe 07230-L06880		Exon 6
148	POU1F1 probe 07241-L06891		Exon 6
154	GHRHR probe 07212-L06862		Exon 9
163	PROP1 probe 07904-L07642		Exon 1
170	LHX4 probe 07232-L08335		Exon 2
175 *	Reference probe 19291-L31751	4q25	
181	POU1F1 probe 07240-L08337		Exon 4
189	GH1 probe 17430-L21430		Exon 4
196	LHX4 probe 07235-L08508		Exon 5
202 ‡	PROP1 probe 21349-L31752		Exon 3
208 † Δ	Reference probe 16366-L31753	12q13	
214	HESX1 probe 07223-L20027		Exon 3
220 ‹ Δ	LHX3 probe 07225-L06875		Intron 1
226	POU1F1 probe 07238-L21577		Exon 2
232 ¶	GH1 probe 17015-SP0428-L21431		Upstream
240 ¶	HESX1 probe 17016-SP0429-L21432		Exon 4
247	POU1F1 probe 07239-L06889		Exon 3
251 ‹	LHX3 probe 07228-L20764		Exon 4
258 ‹	LHX3 probe 15634-L20686		Exon 1
265	GH1 probe 07218-L06868		Exon 3
271 ¶	GHRHR probe 17017-SP0430-L20065		Exon 3
278	LHX4 probe 07233-L20028		Exon 3
285 ‹	LHX3 probe 07907-L20029		Exon 2
292	Reference probe 08790-L11322	10q21	
301 ‹	LHX3 probe 07679-L06879		Exon 5
310	HESX1 probe 07221-L08511		Exon 1
319 *	Reference probe 17521-L21420	2q32	
328	GHRHR probe 07214-L06864		Exon 11
337	POU1F1 probe 07237-L06887		Exon 1
346	GHRHR probe 07207-L06857		Exon 4
355	LHX4 probe 07236-L06886		Exon 6
364	Reference probe 14946-L16679	6q22	
373	GHRHR probe 07213-L06863		Exon 10
382	PROP1 probe 07243-L06893		Exon 2
391	LHX4 probe 07234-L06884		Exon 4
400	GHRHR probe 07678-L06854		Exon 1
409	GHRHR probe 07208-L06858		Exon 5
419	GHRHR probe 07215-L06865		Exon 12
427	HESX1 probe 15091-L16864		Exon 2
436	GHRHR probe 07210-L06860		Exon 7
445	GHRHR probe 07216-L06866		Exon 13
454	LHX4 probe 12872-L06881		Exon 1
463	Reference probe 14308-L15978	15q13	
472	GH1 probe 17018-L20066		Exon 5
480 ‹	LHX3 probe 17019-L20067		Exon 3
490 ¶	GHRHR probe 17020-SP0431-L20068		Exon 6
505 *	Reference probe 14883-L27237	14q11	

a) See above section on exon numbering for more information.

* New in version B2.

‡ Changed in version B2. Minor alteration, no change in sequence detected.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

Table 2. P216-B2 probes arranged according to chromosomal location

Table 2a. *LHX4* gene

Length (nt)	SALSA MLPA probe	<i>LHX4</i> exon ^a	Ligation site NM_033343.3	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	<i>233-235 (Exon 1)</i>		
454	12872-L06881	Exon 1	195-196	GAGAGCGAGAGA-TCTCCGTAGACT	17.8 kb
170	07232-L08335	Exon 2	332-333	GCGCTGGCTGCA-ACCAGCACATCC	18.3 kb
278	07233-L20028	Exon 3	19 nt after exon 3	CATGGCCCCGCA-TGGTCCCCTCTC	4.8 kb
391	07234-L06884	Exon 4	704-705	AGGCTGGAGCTA-AGCGGCCCCGGA	0.5 kb
196	07235-L08508	Exon 5	865-866	AGAAGGGCCAAA-GAGAAACGCCTG	2.4 kb
355	07236-L06886	Exon 6	1058-1059	GGATTTATGGCA-ACGTGGGGGACG	
		<i>stop codon</i>	<i>1403-1405 (Exon 6)</i>		

Table 2b. *HESX1* gene

Length (nt)	SALSA MLPA probe	<i>HESX1</i> exon ^a	Ligation site NM_003865.2	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	<i>335-337 (Exon 1)</i>		
310	07221-L08511	Exon 1	322-321 reverse	CTCTCGTGGTCT-GCACAGAGCAAC	1.1 kb
427	15091-L16864	Exon 2	550-551	ATTCATTCCT-AGCGTGGTGGAT	0.5 kb
214	07223-L20027	Exon 3	748-749	ATCGATATTAGA-GAAGACTTAGCT	0.3 kb
240 Ж	17016-SP0429-L21432	Exon 4	912-913 and 949-950	CTAAACAAGTGA-37 nt spanning oligo-AAATATTAAGTG	
		<i>stop codon</i>	<i>890-892 (Exon 4)</i>		

Table 2c. *POU1F1* gene

Length (nt)	SALSA MLPA probe	<i>POU1F1</i> exon ^a	Ligation site NM_001122757.2	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	<i>126-128 (Exon 1)</i>		
337	07237-L06887	Exon 1	208-209	TCTGATAATGCA-TCACAGTGCTGC	3.0 kb
226	07238-L21577	Exon 2	391-392	TTATGGAAACCA-GCCATCAACCTA	8.9 kb
247	07239-L06889	Exon 3	487-488	TCCTATACACCA-GCCTCTTCTGGC	2.3 kb
181	07240-L08337	Exon 4	696-697	ATGGCTCTGAAT-TCAGTCAAACAA	2.2 kb
148	07241-L06891	Exon 6	942-943	AGATCATGAGGA-TGGCTGAAGAAC	
		<i>stop codon</i>	<i>1077-1079 (Exon 6)</i>		

Table 2d. *PROPI* gene

Length (nt)	SALSA MLPA probe	<i>PROPI</i> exon ^a	Ligation site NM_006261.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	<i>310-312 (Exon 1)</i>		
163	07904-L07642	Exon 1	266-265 reverse	TCTGACTTGAGA-TTCTCTGCTTC	1.6 kb
382	07243-L06893	Exon 2	2 nt before exon 2 reverse	AGCACTCGAGTC-TGAGAACGGAGA	1.7 kb
202	21349-L31752	Exon 3	999-1000	TGAGGTCAAACA-AGTACCACCAAG	
		<i>stop codon</i>	<i>988-990 (Exon 3)</i>		

Table 2e. *GHRHR* gene

Length (nt)	SALSA MLPA probe	<i>GHRHR</i> exon ^a	Ligation site NM_000823.3	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	49-51 (Exon 1)		
400	07678-L06854	Exon 1	98-99	GTTGAGCCCGTT-ACCGACCGTGAG	4.7 kb
137	07905-L31750	Exon 2	28 nt before exon 2 reverse	CAGGATGAGCCA-AGCCATTTGGGT	0.3 kb
271 Ж	17017-SP0430-L20065	Exon 3	236-237 and 260-261	GGATGGGCTGCT-24 nt spanning oligo-CGAGTGGGTCAC	0.8 kb
346	07207-L06857	Exon 4	340-341	GGGATTGTAATA-TCACTGGGCTGGT	1.3 kb
409	07208-L06858	Exon 5	485-486	AGCCCTCTTCGT-GGCCATCACCAT	0.8 kb
490 Ж	17020-SP0431-L20068	Exon 6	9 nt before exon 6 and 530-531	CATTCTCCCAT-27 nt spanning oligo-GAACTACGTCCA	2.1 kb
436	07210-L06860	Exon 7	709-710	TCAGCTGGCTGT-TGGCAGAAGCCG	1.0 kb
154	07212-L06862	Exon 9	914-915	CAAAGGGCCCAT-TGTCCTCTCGGT	0.8 kb
373	07213-L06863	Exon 10	990-991	CTGGAGCCAGCT-CAGGGCAGCCTC	0.6 kb
328	07214-L06864	Exon 11	1060-1061	TGATCCCACTCT-TTGGAAATCACT	0.9 kb
419	07215-L06865	Exon 12	1191-1192	TTCTCAACCAA-GAGGTGTGTGAT	1.9 kb
445	07216-L06866	Exon 13	1337-1336 reverse	GTGGACTCCAGT-GGCGTGATGAGG	
		<i>stop codon</i>	1318-1320 (Exon 13)		

Table 2f. *LHX3* gene

Length (nt)	SALSA MLPA probe	<i>LHX3</i> exon ^a	Ligation site NM_178138.5	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	118-120 (Exon 1)		
258 «	15634-L20686	Exon 1	4 nt after exon 1	GACTCGGGGTAA-GCCCCAGCAGGA	1.8 kb
220 «	07225-L06875	Intron 1	80-81 (NM_014564.4)	GGGGTGGCGTCG-CTGGGCCGGGA	2.4 kb
285 «	07907-L20029	Exon 2	274-275	GCCACTGGCACA-GCAAGTGCTCA	1.0 kb
480 «	17019-L20067	Exon 3	545-544, reverse	TTTCGTAGTCCG-CCTTGACACGA	0.7 kb
251 «	07228-L20764	Exon 4	634-635	TGGAGACGCTGA-AGAGCGTTACA	0.3 kb
301 «	07679-L06879	Exon 5	855-856	GACAGCGTTCAG-GAGGGCAGGAC	1.8 kb
142 «	07230-L06880	Exon 6	1783-1784	TTGGCCTTGCCT-GTCGAGGCAAGA	
		<i>stop codon</i>	1309-1311 (Exon 6)		

Table 2g. *GHI* gene

Length (nt)	SALSA MLPA probe	<i>LHX3</i> exon ^a	Ligation site NM_000515.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	77-79 (Exon 1)		
232 # Ж	17015-SP0428-L21431	Upstream	730 nt and 703 nt before exon 1	ATGGGAGGAAGA-27 nt spanning oligo-TTTCTGTTTCTT	1.5 kb
265 #	07218-L06868	Exon 3	313-312 reverse	GACTCTGAGAAA-CAGAGGGAGGTC	0.3 kb
189 #	17430-L21430	Exon 4	472-473	TACGGCGCTCT-GACAGCAACGTC	0.5 kb
472 #	17018-L20066	Exon 5	786-785 reverse	CTGGAGTGGCAA-CTTCCAGGGCCA	
		<i>stop codon</i>	728-730 (Exon 5)		

a) See above section on exon numbering for more information.

b) Only partial probe sequences are shown. Complete probe sequences are available at www.mlpa.com. Please notify us of any mistakes: info@mlpa.com.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

Related SALSA MLPA probemixes

P217 IGF1R	Contains probes for <i>IGF1R</i> and <i>IGFBP3</i> .
P262 GH1	Contains probes for <i>IGF1</i> , <i>GHR</i> , <i>JAK2</i> , and <i>STAT5B</i> .
P026 Sotos	Contains probes for <i>NSD1</i> and <i>NFIX</i> .
P018 SHOX	Contains probes for <i>SHOX</i> and several other probes in the PAR region.
P210 BTK	Contains probes for <i>BTK</i> .

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P216 Growth Hormone Deficiency mix-1

- Akcay A et al. (2013). Suprasellar mass mimicking a hypothalamic glioma in a patient with a complete PROP1 deletion. *Eur J Hum Genet*, 56(8), 445-451.
- Bertko E et al. (2017). Combined pituitary hormone deficiency due to gross deletions in the POU1F1 (PIT-1) and PROP1 genes. *J Hum Genet*, 62(8), 755.
- Mul D et al. (2012). A mosaic de novo duplication of 17q21–25 is associated with GH insensitivity, disturbed in vitro CD28-mediated signaling, and decreased STAT5B, PI3K, and NF-κB activation. *Eur J Endocrinol*, 166(4), 743-752.
- Petkovic V et al. (2013). Short stature in two siblings heterozygous for a novel bioinactive GH mutant (GH-P59S) suggesting that the mutant also affects secretion of the wild-type GH. *Eur J Endocrinol*, 168(3), K35-K43.
- Selice R et al. (2013). Prostate volume and growth during testosterone replacement therapy is related to visceral obesity in Klinefelter syndrome. *Eur J Endocrinol*, 169(6), 743-749.
- Shimada A et al. (2016). A novel mutation in OTX2 causes combined pituitary hormone deficiency, bilateral microphthalmia, and agenesis of the left internal carotid artery. *Horm Res Paediatr*, 86(1), 62-69.
- Takagi M et al. (2012). Gradual loss of ACTH due to a novel mutation in LHX4: comprehensive mutation screening in Japanese patients with congenital hypopituitarism. *PLoS One*, 7(9), e46008.
- Takagi M et al. (2016). A novel heterozygous intronic mutation in POU1F1 is associated with combined pituitary hormone deficiency. *Endocr. J.* EJ16-0361.

P216 Product history

Version	Modification
B2	Four reference probes have been replaced and three probe lengths have been adjusted.
B1	Reference probes are included and several target probes are exchanged.
A1	First release.

Implemented changes in the product description

Version B2-01 — 14 August 2019 (02P)

- Product description rewritten and adapted to a new template.
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).
- Ligation sites of the probes targeting the *LHX3* gene updated according to new version of the NM_ reference sequence.
- Warning added to Table 2 for probe specificity relying on a single nucleotide difference between target gene and related gene or pseudogene.

Version 13 – 31 August 2015 (55)

- Product description adapted to a new lot (lot number, changes in Table 1 and Table 2, new picture

included).

- Exon numbering of the *LHX3* gene has been updated according to NM_178138.4.
- Various minor textual changes.

Version 12 (54) – 26 February 2015

- New sample picture included in product description.
- Data analysis method has been modified.
- Updated link for "Database of Genomic Variants".

Version 11 (51)

- Product description adapted to a new product version (version number changed, lot number added, changes in Table 1 and Table 2, new picture included).
- New references added on page 1.
- Warning added below Table 2 that the NM_ reference sequence has changed.
- Warning added: ⚠ The significance of exon 1 deletions is not clear as this exon is non-coding and alternative transcript variants using other transcription start sites are known.
- Various minor textual changes on page 1.
- Various minor layout changes.
- Several warnings added in Table 1 and 2.
- Data analysis method has been modified.
- Ligation sites of the probes targeting the *POU1F1* gene updated according to new version of the NM_ reference sequence.

Version 10 (48)

- Electropherogram pictures using the new MLPA buffer (introduced in December 2012) added.

More information: www.mlpa.com; www.mlpa.eu

	MRC-Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mlpa.com (information & technical questions); order@mlpa.com (orders)
Phone	+31 888 657 200