

Product Description SALSA® MLPA® Probemix P015-F2 MECP2

To be used with the MLPA General Protocol.

Version F2. Compared to version F1, two reference probes have been replaced and the 118 nt Y fragment has been removed. For complete product history see page 11.

Catalogue numbers:

- **P015-025R:** SALSA MLPA probemix P015 MECP2, 25 reactions.
- **P015-050R:** SALSA MLPA probemix P015 MECP2, 50 reactions.
- **P015-100R:** SALSA MLPA probemix P015 MECP2, 100 reactions.

To be used in combination with a SALSA® MLPA® reagent kit, available for various number of reactions. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman capillary sequencers, respectively (see www.mlpa.com).

Certificate of Analysis: Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mlpa.com.

Precautions and warnings: For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mlpa.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

Intended use: The SALSA MLPA probemix P015 MECP2 is an in vitro diagnostic (IVD)¹ or a research use only (RUO) assay for the detection of deletions or duplications in the human *MECP2* gene, in order to confirm a potential cause and clinical diagnosis for classic and atypical Rett syndrome, and for *MECP2* duplication syndrome. This assay can also be used for the detection of deletions or duplications in the human *NTNG1* and *CDKL5* genes, in order to confirm a potential cause and clinical diagnosis for atypical Rett syndrome, and for the detection of deletions or duplications in the human *ARX* and *CDKL5* genes, in order to confirm a potential cause and clinical diagnosis for X-linked intellectual disability syndrome. This product can also be used for molecular genetic testing of at-risk family members. This assay is for use with human DNA extracted from peripheral blood.

Deletions or duplications detected with the P015 MECP2 probemix must be confirmed by another technique. In particular, deletions or duplications detected by only a single probe always require validation by another method. Most defects in *MECP2*, *CDKL5*, *ARX* and *NTNG1* genes are point mutations, none of which will be detected by MLPA. It is therefore recommended to use this SALSA MLPA probemix in combination with sequence analysis of the aforementioned genes. Not all exons of *CDKL5*, *ARX* and *NTNG1* genes are covered. The P189 is available for deletion or duplication analysis of other *CDKL5*, *ARX* and *NTNG1* exons. This assay is not intended to be used as a standalone assay for clinical decisions. The results of this test must be interpreted by a clinical molecular geneticist or equivalent.

¹Please note that this probemix is for In Vitro Diagnostic use (IVD) in the countries specified at the end of this product description. In all other countries, the product is for Research Use Only (RUO).

Clinical background: Rett syndrome is a dominant X-linked brain disorder characterized by microcephaly, loss of achieved psychomotor abilities, intellectual disability, and autistic behaviours. The prevalence of Rett syndrome is approximately 1 in 10.000 live female births. Rett syndrome is usually caused by mutations in the *MECP2* gene, although patients with mutations in two other genes, *CDKL5* and *FOXP1*, may also exhibit a Rett-like phenotype. While loss-of-function mutations in *MECP2* result in Rett syndrome, gain-of-function mutations are associated with *MECP2* duplication syndrome which occurs almost exclusively in males. *MECP2* duplication syndrome and Rett syndrome share overlapping clinical phenotypes including intellectual disability, motor deficits, epilepsy, hypotonia, and progressive spasticity.

Atypical Rett syndrome is the early onset seizure variant of Rett syndrome and is characterized by seizures in the first months of life with subsequent development of Rett syndrome features. It has an estimated

prevalence of 1 in 45.000. *MECP2* alterations in Rett syndrome patients comprise mainly point mutations (approximately 80% in classical and 40% in atypical cases), while deletions or duplications are less common (8% in classical and 3% in atypical cases). Approximately 6.5-10% of patients with atypical Rett syndrome have large deletions in *CDKL5* (Rett syndrome database RettBASE).

One report described a patient with atypical Rett syndrome who presented with early onset of epileptic seizures (not infantile spasms) and a *de novo* translocation in intron 6 of *NTNG1* gene (Borg et al. 2005). Please note that MLPA will not detect balanced translocations.

Mutations in *CDKL5* mainly cause X-linked early infantile epileptic encephalopathy, also called X-linked infantile spasm syndrome. *ARX* gene mutations account for 9.5% of all cases of X-linked intellectual disability syndromes, including X-linked lissencephaly with abnormal genitalia, X-linked infantile spasm syndrome, and Partington syndrome (neurological disorder with movement problems). Since these genes are involved in multiple syndromes and are covered in multiple probemixes, Table 1 is provided to give an overview.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK1497/> and <https://www.ncbi.nlm.nih.gov/books/NBK1284/>.

Table 1: Overview on the genes and probemixes related to (atypical) Rett syndrome.

Condition	Genes	Probemix and coverage	Remarks
Classic Rett syndrome	<i>MECP2</i>	P015*: Each exon P374: Exon 2, 3, 4 P245*: Exon 1, 3, 4	MECP2 probe exon 2 in P374 is different from P015. MECP2 probes in P245 are identical to P015 and P374 (different length).
MECP2 duplication syndrome	<i>MECP2</i>	P015*: Each exon P374: Exon 2, 3, 4 P245*: Exon 1, 3, 4	MECP2 probe exon 2 in P374 is different from P015. MECP2 probes in P245 are identical to P015 and P374 (different length).
Atypical Rett syndrome	<i>MECP2</i>	P015*: Each exon P374: Exon 2, 3, 4 P245*: Exon 1, 3, 4	MECP2 probe exon 2 in P374 is different from P015. MECP2 probes in P245 are identical to P015 and P374 (different length).
	<i>CDKL5</i>	P189*: Each exon P015*: Exons 3, 6, 9, 10 P259: Exon 2, 17, 20	CDKL5 probes in P015 and in P259 are identical to P189 (different length).
	<i>NTNG1</i>	P189*: Each exon P015*: Exons 2, 3, 5, 6	NTNG1 probes in P015 are identical to P189 except for exon 2 (different length).
FOXG1 syndrome	<i>FOXG1</i>	P189*: Each exon P075: Each exon + upstream P395: Each exon + upstream	FOXG1 probes in P189 are identical to P075 and P395 (different length). One FOXG1 exon 1 probe is different between P075 and P395.
X-linked intellectual disability syndrome	<i>CDKL5</i>	P189*: Each exon P015*: Exons 3, 6, 9, 10 P259: Exon 2, 17, 20	CDKL5 probes in P015 and in P259 are identical to P189 (different length).
	<i>ARX</i>	P189*: Each exon P015*: Exons 1, 5 P106: Exons 1, 2, 4	ARX probes in P015 and P106 are identical to P189 (different length).

*IVD. All other for Research Use Only (RUO).

Gene structure:

The *MECP2* gene (4 exons) spans ~76 kb of genomic DNA on chromosome Xq28. The status of *MECP2* LRG_764 is pending, but available at <http://www.lrg-sequence.org>, and identical to GenBank NG_007107.2. The *CDKL5* gene (21 exons) spans ~228 kb of genomic DNA on chromosome Xp22.13. The *ARX* gene (5 exons) spans ~12 kb of genomic DNA on chromosome Xp21.3. The *NTNG1* gene (6 exons) spans ~345 kb of genomic DNA on chromosome 1p13.3.

Transcript variants:

The *MECP2* LRG_764 is identical to RefSeq transcript NM_004992.3 (10241 nt, coding sequence 227-1687): <https://www.ncbi.nlm.nih.gov/gene/4204>. *MECP2* MLPA probes are indicated according to NM_004992.3 in Table 3.

Transcript variant NM_003159.2 (3434 nt, coding sequence 254-3346) for *CDKL5* is described: <https://www.ncbi.nlm.nih.gov/gene/6792>. This sequence is a reference standard in the RefSeqGene project. *CDKL5* MLPA probes are indicated according to NM_003159.2 in Table 3.

Transcript variant NM_139058.2 (2885 nt, coding sequence 212-1900) for *ARX* is described: <https://www.ncbi.nlm.nih.gov/gene/170302>. This sequence is a reference standard in the RefSeqGene project. *ARX* MLPA probes are indicated according to NM_139058.2 in Table 3.

Transcript variant NM_014917.3 (6199 nt, coding sequence 808-2124) for *NTNG1* is described: <https://www.ncbi.nlm.nih.gov/gene/22854>. This sequence is a reference standard in the RefSeqGene project. *NTNG1* MLPA probes are indicated according to NM_014917.3 in Table 3.

Exon numbering: The exon numbering used in this P015-F2 *MECP2* product description is the exon numbering from the RefSeq transcript NM_004992.3 for *MECP2*, which is identical the LRG_764 sequence. For *CDKL5* the exon numbering used is from the RefSeq transcript variant NM_003159.2, for *ARX* from the RefSeq transcript NM_139058.2, and for *NTNG1* from the RefSeq transcript NM_014917.3. The exon numbering and NM sequences used are from 07/2018, but can be changed by NCBI after the release of the product description.

Probemix content: The P015-F2 probemix contains 46 probes with amplification products between 130 nt and 467 nt. The P015-F2 probemix contains 17 probes for the *MECP2* gene, covering each exon. Furthermore, several probes are present for genes in close proximity to *MECP2*. One of these probes is located within the pseudo autosomal region 2 (*PAR2*). The P015-F2 probemix contains four *CDKL5* probes, two *ARX* probes, and four *NTNG1* probes. More probes for the *CDKL5*, *ARX* and *NTNG1* genes are present in the P189 *CDKL5/ARX/FOXG1* probemix. In addition, 10 reference probes are included in this probemix, detecting several different autosomal chromosomal locations. The identity of the genes detected by the reference probes is available online (www.mlpa.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity Fragments (Q-fragments), two DNA Denaturation Fragments (D-fragments), one benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mlpa.com.

Length (nt)	Name
64-70-76-82	Q-fragments (Only visible with <100 ng sample DNA)
88-96	D-fragments (Low signal of 88 or 96 fragment indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique: The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mlpa.com).

MLPA technique validation: Internal validation of the MLPA technique using 16 DNA samples from healthy individuals of the same sex is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation <0.10 for all probes over the experiment.

Required specimens: Extracted DNA from peripheral blood, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples: All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from unrelated individuals who are from families without a history of Rett syndrome. It is recommended to use reference and patient samples of the same sex to facilitate interpretation. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol.

Positive control DNA samples: MRC-Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Biobank (<https://catalog.coriell.org>) and DSMZ (<https://www.dsmz.de/home.html>) have a diverse collection of biological resources which may be used as a positive control DNA sample in your MLPA experiments. Female sample ID number NA23676 from the Coriell Institute has been tested at MRC-Holland and can be used as positive control sample to detect a heterozygous whole gene duplication of *MECP2*. The quality of cell lines can change, therefore samples should be validated before use.

Performance characteristics: The frequency of *MECP2* deletions or duplications in classic Rett syndrome is 8% and 3% in atypical cases. Approximately 6.5-10% of patients with atypical Rett syndrome have large deletions in *CDKL5* (Rett syndrome database RettBASE). The frequency of *ARX* deletions or duplications is ~1% in X-linked intellectual disability syndrome cases <https://www.omim.org/entry/300382>. No deletions or duplications in the *NTNG1* gene have been described so far. The analytical sensitivity and specificity for the detection of deletions or duplications in the *MECP2*, *CDKL5*, *ARX*, and *NTNG1* genes is very high and can be considered >99% (based on a 2000-2016 literature review).

Analytical performance can be compromised by: SNPs or other polymorphisms (e.g. indels) in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.

Data analysis: Coffalyser.Net software must be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mlpa.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results: The expected results for the probes detecting autosomal sequences are allele copy numbers of 2 (normal), 1 (heterozygous deletion), or 3 (heterozygous duplication). The same results can be expected for the X-chromosome-specific probes in female samples. For the X-chromosome-specific probes in male samples, expected copy numbers are 1 (normal), 0 (deletion) or 2 (duplication).

The standard deviation of all probes in the reference samples should be <0.10. When this criteria is fulfilled, the following cut-off values for the dosage quotient (DQ) of the probes can be used to interpret MLPA results when **reference samples of the same sex** have been used:

Copy Number status		Dosage quotient
Autosomal sequences X-chromosome sequences in female	X-chromosome sequences in male	
Normal	Normal	$0.80 < DQ < 1.20$
Homozygous deletion	Deletion	$DQ = 0$
Heterozygous deletion	-	$0.40 < DQ < 0.65$
Heterozygous duplication	-	$1.30 < DQ < 1.65$
Heterozygous triplication/ Homozygous duplication	Duplication	$1.75 < DQ < 2.15$
Ambiguous copy number	Ambiguous copy number	All other values

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Incomplete DNA denaturation (e.g. due to salt contamination) can lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *MECP2*, *CDKL5*, *ARX*, *NTNG1* genes. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might in some cases not result in inactivation of that gene copy.
- Copy number changes detected by reference probes are unlikely to have any relation to the condition tested for.
- Deletion of a probe's recognition sequence on the X-chromosome will lead to a complete absence of the corresponding probe amplification product in males, whereas female heterozygotes are recognizable by a 35-50% reduction in relative peak height.

Limitations of the procedure:

- In most populations, the major cause of genetic defects in the *MECP2*, *CDKL5*, *ARX*, and *NTNG1* genes are small (point) mutations, none of which will be detected by using SALSA MLPA probemix P015-MECP2.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNPs, point mutations, small indels) in the target sequence detected by a probe can cause false positive results. Mutations/SNPs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.
- Not all exons of the *CDKL5*, *ARX* and *NGNT1* genes are covered in the P015 MECP2 probemix.

Confirmation of results: Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by one or more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

RettBASE mutation database: <http://mecp2.chw.edu.au/mecp2/index.php>. We strongly encourage users to deposit positive results in the RettBASE mutation database. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNPs and unusual results (e.g., a duplication of *MECP2* exons 1 and 3 but not exon 2) to MRC-Holland: info@mlpa.com.

Table 2. SALSA MLPA P015-F2 MECP2 probemix

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^(a)				
		Reference	MECP2 Xq28	CDKL5 Xp22.13	ARX Xp21.3	NTNG1 1p13
64-105	Control fragments – see table in probemix content section for more information					
130	Reference probe 00797-L21056	5q31				
136 « ↖	GDI1 probe 16875-L19669		Xq28			
143	CDKL5 probe 06466-L06567			Exon 10		
149	NTNG1 probe 18447-L24143				Exon 2	
154 ~	MECP2 probe 18442-L24243		Exon 4b			
161 «	ARX probe 18440-L24338			Exon 1		
166	Reference probe 08222-L24146	10q26				
172 « ↖	SLC6A8 probe 01879-L24147		Xq28			
178 « ↖	IRAK1 probe 10835-L24148		Xq28			
182 «	MECP2 probe 10836-L24339		Exon 1			
190	CDKL5 probe 06458-L23617			Exon 3		
196	MECP2 probe 10839-L23618		Exon 3			
202 «	MECP2 probe 03409-L16570		Exon 1			
208	Reference probe 09865-L08705	13q32				
214 «	MECP2 probe 03770-L13387		Upstream			
221 ↖	F8 probe 06288-L05892		Xq28			
229 ~	MECP2 probe 18441-L12494		Exon 4b			
235 «	ARX probe 06455-L21229			Exon 5		
241 ↖	L1CAM probe 07051-L14553		Xq28			
247 *	Reference probe 05959-L05376	7p11				
254 «	MECP2 probe 03768-L23833		Exon 2			
260 « ±	MECP2 probe 01769-L23834		Exon 4b			
266 *	Reference probe 10728-L22588	6p12				
274 ~	MECP2 probe 01768-L13824		Exon 4b			
280	Reference probe 16434-L24149	18q21				
286	CDKL5 probe 06462-L24150			Exon 6		
292 ~	MECP2 probe 18444-L24151		Exon 4b			
299	Reference probe 16621-L24152	1q32				
305 «	MECP2 probe 02002-L24153		Exon 1			
312	NTNG1 probe 06487-L24154			Exon 5		
319 ↖	DKC1 probe 18761-L24340		Xq28			
338	Reference probe 09776-L24156	15q21				
346 ~	MECP2 probe 01347-L24157		Exon 4b			
356	MECP2 probe 10841-L24158		Exon 4b			
365	MECP2 probe 01348-L24159		Exon 3			
373 «	MECP2 probe 01349-L24160		Exon 2			
383 ↖	IDH3G probe 01887-L23933		Xq28			
391	MECP2 probe 14737-L24161		Exon 3			
400	CDKL5 probe 06465-L05991			Exon 9		
409	Reference probe 17462-L21218	12p13				
418 «	MECP2 probe 18446-L23620		Exon 4b			
427	NTNG1 probe 06488-L23934			Exon 6		
438	NTNG1 probe 06483-L24162			Exon 3		
448 « ↖	FLNA probe 04138-L24163		Xq28			
457 ↖	VAMP7 probe 01094-L24164		Xq28- PAR2			
467	Reference probe 02674-L24165	11q22				

(a) The exon numbering used in this P015-F2 product description is identical for *MECP2* to the exon numbering from the LRG_764 sequence and to the RefSeq transcript NM_004992.3, for *CDKL5* to the RefSeq transcript variant NM_003159.2, for *ARX* to the RefSeq transcript NM_139058.2, and for *NTNG1* to the RefSeq transcript NM_014917.3. Exon numbering used here may differ from literature. Please notify us of any mistakes.

* New in version F2 (from lot F2-0118 onwards).

‡ SNP rs267608346 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

- « Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.
- ~ Probe located within, or close to, a hotspot region of the *MECP2* gene in which small deletions occur frequently (Laccone et al. 2004, Huppke et al. 2005). In addition to small deletions, many SNPs are also present in this region, making this probe more prone to false positive deletion results.
- ▾ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

Table 3. P015 probes arranged according to chromosomal location

Table 3a. *MECP2* / Xq28 region (from telomere to centromere)

Length (nt)	SALSA MLPA probe	Gene/ Exon ^(a)	Ligation site ^(b)	Partial sequence ^(c) (24 nt adjacent to ligation site)	Distance to next probe
<i>q-telomere</i>					<i>~90.0 kb</i>
457 ▸	01094-L24164	VAMP7 Exon 8 (PAR2 region)		TGTGGGAAAAGT-GTTTCCATTCTG	957.0 kb
221 ▸	06288-L05892	F8		TTCAGGGAGTCT-GGCCAAGGAAAA	218.1 kb
319 ▸	18761-L24340	DKC1		CATGATGTGCTT-GATGCTCAGTGG	331.9 kb
136 « ▸	16875-L19669	GDI1		CCTGACCATGGA-CGAGGAATACGA	79.7 kb
448 « ▸	04138-L24163	FLNA		TGACGGCACGTA-TACAGTGGCCTA	221.7 kb
MECP2 NM_004992.3					
		<i>start codon</i>	<i>227-229 (ex 2)</i>		
214 «	03770-L13387	Upstream	1093 nt before ex 1	GCAAGAATGTTA-GTTTGCTGTCTG	0.7 kb
305 «	02002-L24153	Exon 1	380 nt before ex 1, reverse	GGGACGCCTGTT-TGCGCTGCTCTG	0.1 kb
202 «	03409-L16570	Exon 1	318 nt before ex 1	CATTAATCCTTA-ACATTCAAATTC	0.4 kb
182 «	10836-L24339	Exon 1	52-53	GGGCTGTGGTAA-AAGCCGTCGGGA	5.3 kb
254 «	03768-L23833	Exon 2	44 nt before ex 2	GAAAAAGGTCGT-GCAGCTCAATGG	0.1 kb
373 «	01349-L24160	Exon 2	188-189	GACTCCCCAGAA-TACACCTTGCTT	59.6 kb
196	10839-L23618	Exon 3	47 nt before ex 3	ACTTGTTCTGCA-GACTGGCATGTT	0.2 kb
365	01348-L24159	Exon 3	385-386	GCCCACCACTCT-GCTGAGCCCGCA	0.2 kb
391	14737-L24161	Exon 3	582-583	CTCTGCTGGGAA-GTATGATGTGTA	0.7 kb
356	10841-L24158	Exon 4b	32 nt before ex 4b	AGAGCCTCTAAT-TGTTCTTGTGT	0.3 kb
229 ~	18441-L12494	Exon 4b	896-897	TCCTTGTCAGA-TGCCTTTTCAA	0.3 kb
346 ~	01347-L24157	Exon 4b	1243-1244	CTGAAGACCTGT-AAGAGCCCTGGG	0.2 kb
154 ~	18442-L24243	Exon 4b	1479-1480	AGAGGAGAAGAT-GCCCAGAGGAGG	0.1 kb
274 ~	01768-L13824	Exon 4b	1622-1623	TTTCATCCTCCA-TGCCAAGGCCAA	0.1 kb
292 ~	18444-L24151	Exon 4b	1682-1683	CCGAGAGAGTTA-GTGTACTTTACA	2.2 kb
418 «	18446-L23620	Exon 4b	3871-3872	TGCTGCCATGAA-CTGTCAAGTGTG	3.2 kb
260 « ±	01769-L23834	Exon 4b	7036-7037	CAGTAACACATA-GACTGTGCGCAT	5.9 kb
		<i>stop codon</i>	<i>1685-1687 (ex 4b)</i>		
178 « ▸	10835-L24148	IRAK1		TTTATGAAGCTT-TTCCAGGCTCCC	154.4 kb
241 ▸	07051-L14553	L1CAM		CAGCGGGTGAAA-ACTACAGTGTCTG	70.3 kb
383 ▸	01887-L23933	IDH3G		TCCCCGAACTT-CGCACCCCGTCTG	99.8 kb
172 « ▸ #	01879-L24147	SLC6A8		ACCCCGCTGGTC-TGCATGGTAAGG	-

Table 3b. *CDKL5* (chromosome X)

Length (nt)	SALSA MLPA probe	CDKL5 Exon ^(a)	Ligation site ^(b) NM_003159.2	Partial sequence ^(c) (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	<i>254-256 (ex 2)</i>		
190	06458-L23617	Exon 3	42 nt before ex 3	GAGCTTTGTAGT-TTGTATGCGTGC	69.1 kb
286	06462-L24150	Exon 6	573-574	GCCAAATGGAGT-TCCACCTGAGAA	8.2 kb
400	06465-L05991	Exon 9	913-914	ATTGACCAACTT-TTTACTATTTCAG	7.3 kb
143	06466-L06567	Exon 10	1026-1027	TCCTCAGTCCTT-GGAAAGAAGATA	-
		<i>stop codon</i>	<i>3344-3346 (ex 21)</i>		

Table 3c. *ARX* (chromosome X)

Length (nt)	SALSA MLPA probe	ARX Exon ^(a)	Ligation site ^(b) NM_139058.2	Partial sequence ^(c) (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	212-214 (ex 1)		
161 «	18440-L24338	Exon 1	98-99	AGATCGCAATAA-TATCCGTTATAA	11.2 kb
235 «	06455-L21229	Exon 5	1964-1965	CAGCACCCTCA-AGACCAAATGGA	-
		<i>stop codon</i>	1898-1900 (ex 5)		

Table 3d. *NTNG1* (chromosome 1)

Length (nt)	SALSA MLPA probe	NTNG1 Exon ^(a)	Ligation site ^(b) NM_014917.3	Partial sequence ^(c) (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	808-810 (ex 2)		
149	18447-L24143	Exon 2	578-579	ACAAGTATGTTA-GGCTTCCACCAA	176.5 kb
438	06483-L24162	Exon 3	1606-1607	GGATAAGGCTGT-TAAGACCAGCCG	83.1 kb
312	06487-L24154	Exon 5	227 nt after ex 5	TATGACTTTTCT-GACTACTCTTAA	72.7 kb
427	06488-L23934	Exon 6	1907-1908	GAATGTCTGCGA-CAACGAGCTCCT	-
		<i>stop codon</i>	2122-2124 (ex 6)		

(a) The exon numbering used in this P015-F2 MECP2 product description is identical for *MECP2* to the exon numbering of the LRG_764 sequence and to the RefSeq transcript NM_004992.3, for *CDKL5* to the RefSeq transcript variant NM_003159.2, for *ARX* to the RefSeq transcript NM_139058.2, and for *NTNG1* to the RefSeq transcript NM_014917.3.

(b) Ligation sites of the P015 MECP2 MLPA probes are indicated according to RefSeq sequence NM_004992.3 for *MECP2* containing 4 exons, for *CDKL5* to the RefSeq sequence NM_003159.2 containing 21 exons, for *ARX* to the RefSeq sequence NM_139058.2 containing 5 exons, and for *NTNG1* to the RefSeq sequence NM_014917.3 containing 6 exons.

(c) Only partial probe sequences are shown. Complete probe sequences are available at www.mlpa.com. Please notify us of any mistakes: info@mlpa.com.

~ Probe located within, or close to, a hotspot region of the *MECP2* gene in which small deletions occur frequently (Laccone et al. 2004, Huppke et al. 2005). In addition to small deletions, many SNPs are also present in this region, making this probe more prone to false positive deletion results.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

‡ SNP rs267608346 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

∩ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

Related SALSA® MLPA® probemixes

- P189 CDKL5/ARX/FOXG1: Contains probes for the *CDKL5*, *NTNG1*, *ARX* and *FOXG1* genes.
- P075 TCF4-FOXG1: Contains probes for the *TCF4* and *FOXG1* genes (Pitt-Hopkins syndrome).
- P336 UBE3A: For the detection of small rearrangements in the *UBE3A* gene, involved in Angelman syndrome.
- P395 MEFC2-FOXG1: Contains probes for the *MEFC2* and *FOXG1* genes (mental retardation).
- P106 MRX: Contains probes for various genes that cause X-linked mental retardation when defect.
- ME028 Prader Willi / Angelman syndrome: For the detection of both copy number as well as methylation changes of the Prader Willi / Angelman region.
- P245 Microdeletion Syndromes-1A: Includes probes for *MECP2*.
- P374: Microdeletion Syndromes 8: Includes probes for *MECP2*.
- P259: RPS6KA3 (Coffin-Lowry syndrome): Includes probes for *CDKL5*.

References

- Borg I et al. (2005). Disruption of Netrin G1 by a balanced chromosome translocation in a girl with Rett syndrome. *Eur J Hum Genet.* 13(8):921-7
- Huppke P et al. (2005). Molecular diagnosis of Rett syndrome. *J Child Neurol.* 20(9):732-6

- Laccone F et al. (2004). Large deletions of the MECP2 gene detected by gene dosage analysis in patients with Rett syndrome. *Hum Mutat.* 23(3):234-44
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801

Selected publications using SALSA® MLPA® Probemix P015 MECP2

- Aradhya S et al. (2012) Exon-level array CGH in a large clinical cohort demonstrates increased sensitivity of diagnostic testing for Mendelian disorders. *Genet Med.* 14:594-603
- Archer HL et al. (2006) Gross rearrangements of the MECP2 gene are found in both classical and atypical Rett syndrome patients. *J Med Genet.* 43:451-456
- Ariani F et al. (2004) Real-time quantitative PCR as a routine method for screening large rearrangements in Rett syndrome: Report of one case of MECP2 deletion and one case of MECP2 duplication. *Hum Mutat.* 24:172-177
- Bartholdi D et al. (2006) Clinical profiles of four patients with Rett syndrome carrying a novel exon 1 mutation or genomic rearrangement in the MECP2 gene. *Clin Genet.* 69:319-326
- Bebbington A et al. (2012) The phenotype associated with a large deletion on MECP2. *Eur J Hum Genet* 20:921-927
- Bijlsma E et al. (2012) Xq28 duplications including MECP2 in five females: Expanding the phenotype to severe mental retardation. *Eur J Med Genet.* 55:404-413
- Colak D et al. (2011) Genomic and transcriptomic analyses distinguish classic Rett and Rett-like syndrome and reveals shared altered pathways. *Genomics.* 97:19-28
- Erlandson A et al. (2003) Multiplex ligation-dependent probe amplification (MLPA) detects large deletions in the MECP2 gene of Swedish Rett syndrome patients. *Genet Test.* 7:329-332
- Fernández RM et al. (2010) Novel association of severe neonatal encephalopathy and Hirschsprung disease in a male with a duplication at the Xq28 region. *BMC Med Genet.* 11:137
- Fieremans N et al. (2014) De novo MECP2 duplications in two females with intellectual disability and unfavorable complete skewed X-inactivation. *Hum Genet.* 133:1359-1367
- Hara M et al. (2015) De novo SHANK3 mutation causes Rett syndrome-like phenotype in a female patient. *Am J Med Genet A.* 167:1593-1596
- Hardwick SA et al. (2007) Delineation of large deletions of the MECP2 gene in Rett syndrome patients, including a familial case with a male proband. *Eur J Hum Genet.* 15:1218-1229
- Jezela-Stanek A et al. (2011) Cryptic x; autosome translocation in a boy—Delineation of the phenotype. *Pediatr Neurol.* 44:221-224
- Kalman LV et al. (2014) Development of a genomic DNA reference material panel for Rett syndrome (MECP2-related disorders) genetic testing. *J Mol Diagn.* 16:273-279
- Kerkhof J et al. (2017). Clinical Validation of Copy Number Variant Detection from Targeted Next-Generation Sequencing Panels. *J Mol Diagn.* 19:905-920.
- Khajuria R et al. (2009) Rapid detection of deletions in hotspot C-terminal segment region of MECP2 by routine PCR method: report of two classical Rett syndrome patients of Indian origin. *Genet Test Mol Biomarkers.* 13:277-280
- Li X et al. (2017). Clinical and molecular genetic characterization of familial MECP2 duplication syndrome in a Chinese family. *BMC Med Genet.* 18:131
- Makrythanasis P et al. (2010) De novo duplication of MECP2 in a girl with mental retardation and no obvious dysmorphic features. *Clin Genet.* 78:175-180
- Mayo S et al. (2011) De novo interstitial triplication of MECP2 in a girl with neurodevelopmental disorder and random X chromosome inactivation. *Cytogenet Genome Res.* 135:93-101
- Parrini E et al. (2017) Diagnostic Targeted Resequencing in 349 Patients with Drug-Resistant Pediatric Epilepsies Identifies Causative Mutations in 30 Different Genes. *Hum Mutat.* 38:216-225
- Petel-Galil Y et al. (2006) Comprehensive diagnosis of Rett's syndrome relying on genetic, epigenetic and expression evidence of deficiency of the methyl-CpG-binding protein 2 gene: study of a cohort of Israeli patients. *J Med Genet.* 43:e56
- Pieras JI et al. (2011) Somatic mosaicism for Y120X mutation in the MECP2 gene causes atypical Rett syndrome in a male. *Brain Dev.* 33:608-611

- Pratt DW et al. (2013) Genotyping FOXP1 Mutations in Patients with Clinical Evidence of the FOXP1 Syndrome. *Mol Syndromol.* 3:284-287
- Raizis AM et al. (2009) Spectrum of MECP2 mutations in New Zealand Rett syndrome patients. *N Z Med J.* 122:21-28
- Ravn K et al. (2005) Large genomic rearrangements in MECP2. *Hum Mutat.* 25:324
- Sanmann JN et al. (2012) Characterization of six novel patients with MECP2 duplications due to unbalanced rearrangements of the X chromosome. *Am J Med Genet A.* 158:1285-1291
- Saunders CJ et al. (2010) Allele drop-out in the MECP2 gene due to G-quadruplex and i-motif sequences when using polymerase chain reaction-based diagnosis for Rett syndrome. *Genet Test Mol Biomarkers.* 14:241-247
- Scala E et al. (2007) MECP2 deletions and genotype-phenotype correlation in Rett syndrome. *Am J Med Genet A.* 143A:2775-2784
- Todorov T et al. (2010) Clinical and molecular data on mental retardation in Bulgaria. *Balkan J Med Genet.* 13:11-25
- Trump N et al. (2016) Improving diagnosis and broadening the phenotypes in early-onset seizure and severe developmental delay disorders through gene panel analysis. *J Med Genet.* 0:1-8
- Utine G et al. (2011) Searching for copy number changes in nonsyndromic X-linked intellectual disability. *Mol syndromol.* 2:64-71
- Vieira JP et al. (2015) Variant Rett syndrome in a girl with a pericentric X-chromosome inversion leading to epigenetic changes and overexpression of the MECP2 gene. *Int J Dev Neurosci.* 46:82-87
- Westenberger A et al. (2013) X-linked Dystonia-Parkinsonism manifesting in a female patient due to atypical turner syndrome. *Mov Disord.* 28:675-678
- Wooderchak-Donahue W et al. (2011) Verification of multiplex ligation-dependent probe amplification probes in the absence of positive samples. *Genet Test Mol Biomarkers.* 15:793-799
- Yon DK et al. (2017) A sibship with duplication of Xq28 inherited from the mother; genomic characterization and clinical outcomes. *BMC Med Genet.* 18:30
- Zahorakova D et al. (2016) MECP2 mutations in Czech patients with Rett syndrome and Rett-like phenotypes: novel mutations, genotype-phenotype correlations and validation of high-resolution melting analysis for mutation scanning. *J Hum Genet.* 1-9
- Zahorakova D et al. (2007) Mutation analysis of the MECP2 gene in patients of Slavic origin with Rett syndrome: novel mutations and polymorphisms. *J Hum Genet.* 52:342-348
- Zhang Q et al. (2017) Familial cases and male cases with MECP2 mutations. *Am J Med Genet.* 9999:1-7
- Zhang Q et al. (2015) MECP2 duplication syndrome in a Chinese family. *BMC Med Genet.* 16:112
- Zhang X et al. (2012) Molecular characteristics of Chinese patients with Rett syndrome. *Eur J Med Genet.* 55:677-681
- Zhu HY et al. (2011) [Molecular genetic study of MECP2 gene for a patient with typical Rett syndrome]. *Zhonghua Yi Xue Yi Chuan Xue Za Zhi* 28:625-629

P015 MECP2 Product history	
<i>Version</i>	<i>Modification</i>
F2	Two reference probes have been replaced, the 118 nt Y fragment has been removed.
F1	Four MECP2 probes (exons 3 and 4) and two additional reference probes have been included. One ARX, one NTNG1 and several reference probes have been replaced. One CDKL5 and one IRAK1 probe have been removed. The 88 and 96 nt denaturation control probes have been replaced (QDX2).
E1	One MECP2 exon 3 probe has been removed.
D2	Several new MECP2 probes have been included and extra probes up- and downstream of MECP2 were added. In addition, some probes for the CDKL5, ARX and NTNG1 genes are included.
D1	Test Lot.
C2	Extra control fragments at 88, 96, 100 and 105 nt have been added.
B1	One extra probe for MECP2 exon 1, three extra MECP2 exon 4 probes and extra reference probes have been included.
A1	First release.

Implemented changes in the product description

Version F2-03 – 18 May 2020 (04)

- Israel added as country with IVD status.

Version F2-02 – 30 July 2018 (04)

- Warning added to Table 3 for probe specificity relying on a single nucleotide difference between target gene and related gene or pseudogene.
- Selected publications were updated using P015.

Version F2-01 – 01 March 2018 (04)

- Product description restructured and adapted to a new template.
- Product description adapted to a new product version (version number changed, changes in Table 2 and Table 3).
- Table 1 is included to give an overview of the included genes in different probemixes related to (atypical) Rett syndrome.

Version 15 - 27 January 2017 (55)

- Probe number (182nt MECP2 probe 10836-L24339) adjusted as it was incorrect displayed in previous product description version.

Version 14 - 15 November 2015 (55)

- Lot number added, new pictures included.
- CDKL5 exon numbering adjusted.
- Various minor textual changes.

Version 13 (53)

- Various minor layout changes.
- Updated link for "Database of genomic variants".
- "Peak area" replaced with "peak height".

Version 12 (52)

- Reference to previous lot adjusted.

Version 11 (49)

- Product description adapted to a new version (lot number added, small changes in Table 1 and Table 2, small textual changes, new pictures included).

More information: www.mlpa.com; www.mlpa.eu

	MRC-Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mlpa.com (information & technical questions); order@mlpa.com (orders)
Phone	+31 888 657 200

IVD	EUROPE*  ISRAEL
RUO	ALL OTHER COUNTRIES

*comprising EU (candidate) member states and members of the European Free Trade Association (EFTA).
The product is for RUO in all other European countries.