

# Product Description

## SALSA® MLPA® Probemix P368-B1 DCLRE1C

To be used with the MLPA General Protocol.

### Version B1

For complete product history see page 7.

### Catalogue numbers:

- **P368-025R:** SALSA MLPA Probemix P368 DCLRE1C, 25 reactions.
- **P368-050R:** SALSA MLPA Probemix P368 DCLRE1C, 50 reactions.
- **P368-100R:** SALSA MLPA Probemix P368 DCLRE1C, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see [www.mrcholland.com](http://www.mrcholland.com)).

### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at [www.mrcholland.com](http://www.mrcholland.com).

### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: [www.mrcholland.com](http://www.mrcholland.com). It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

### General information

The SALSA MLPA Probemix P368 DCLRE1C is a **research use only (RUO)** assay for the detection of deletions or duplications in the *DCLRE1C*, *MEIG1*, and *DCLRE1CP1* genes, which are associated with radiosensitive Severe Combined Immunodeficiency (SCID) or radiosensitive Omenn syndrome.

The nuclease ARTEMIS (encoded by *DCLRE1C*, DNA crosslink repair 1C) is an essential factor in the process of V(D)J recombination and an important component of the nonhomologous end joining (NHEJ) DNA double-strand break (DSB) repair pathway. Patients with mutations in the *DCLRE1C* gene suffer from radiosensitive Severe Combined Immunodeficiency (SCID) or radiosensitive Omenn syndrome.

Pannicke et al. (2010) found that by far the most frequent mutations (59%) were gross deletions of exons 1-3 or 1-4, due to homologous recombination of the wild-type *DCLRE1C* gene with the *DCLRE1CP1* gene, located 62.5 kilobases (kb) upstream of *DCLRE1C*. *MEIG1* located between the gene and pseudogene is lost by these gross deletions.

**This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.**

### Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM\_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

### Exon numbering

The *DCLRE1C*, *MEIG1*, and *DCLRE1CP1* exon numbering used in this P368-B1 DCLRE1C product description is the exon numbering from the LRG\_54, NM\_001080836.3, and NR\_144469.1 sequences. The exon numbering of the NM\_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM\_ sequence and exon numbering may not be up-to-date.

### Probemix content

The SALSA MLPA Probemix P368-B1 DCLRE1C contains 33 MLPA probes with amplification products between 136 and 420 nucleotides (nt). This includes 18 probes for the *DCLRE1C* gene, one probe for each exon of the gene and two probes for exon 1 and exon 14, and five probes for the 10p13 region upstream of *DCLRE1C*, comprising three probes for *MEIG1* and two probes for the *DCLRE1CP1* gene. In addition, ten reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online ([www.mrcholland.com](http://www.mrcholland.com)).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at [www.mrcholland.com](http://www.mrcholland.com).

| Length (nt) | Name   |
|-------------|--|
| 64-70-76-82 | Q-fragments (only visible with <100 ng sample DNA)         |
| 88-96       | D-fragments (low signal indicates incomplete denaturation) |
| 92          | Benchmark fragment   |
| 100         | X-fragment (X chromosome specific)                         |
| 105         | Y-fragment (Y chromosome specific)                         |

### MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

### Required specimens

Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

### Reference samples

A sufficient number ( $\geq 3$ ) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of radiosensitive SCID or Omenn syndrome. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

## Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at [www.mrcholland.com](http://www.mrcholland.com). Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

## Interpretation of results

The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$  and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

| Copy number status                               | Final ratio (FR)   |
|--|--------------------|
| Normal   | $0.80 < FR < 1.20$ |
| Homozygous deletion                              | $FR = 0$           |
| Heterozygous deletion                            | $0.40 < FR < 0.65$ |
| Heterozygous duplication                         | $1.30 < FR < 1.65$ |
| Heterozygous triplication/homozygous duplication | $1.75 < FR < 2.15$ |
| Ambiguous copy number                            | All other values   |

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun

the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

### Limitations of the procedure

- In most populations, the major cause of genetic defects in the *DCLRE1C*, *MEIG1*, and *DCLRE1CP1* genes are small (point) mutations, none of which will not be detected by using SALSA MLPA Probemix P368 DCLRE1C.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

### Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

### LOVD mutation database

<https://databases.lovd.nl/shared/genes>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *DCLRE1C* exons 5 and 7 but not exon 6) to MRC Holland: [info@mrcholland.com](mailto:info@mrcholland.com).

**Table 1. SALSA MLPA Probemix P368-B1 DCLRE1C**

| Length (nt) | SALSA MLPA probe   | Chromosomal position (hg18) <sup>a</sup> |                 |                 |
|-------------|--|--|-----------------|-----------------|
|             |  | Reference                                | DCLRE1C         | Other           |
| 64-105      | Control fragments – see table in probemix content section for more information |  |                 |                 |
| 136         | Reference probe 09285-L09516   | 11q                                      |                 |                 |
| 142         | <b>DCLRE1C probe</b> 15595-L17451  |  | <b>Exon 8</b>   |                 |
| 148 ~       | <b>MEIG1 probe</b> 15596-L18176  |  |                 | <b>Exon 1</b>   |
| 154         | <b>DCLRE1C probe</b> 15597-L17453  |  | <b>Exon 6</b>   |                 |
| 166         | <b>DCLRE1C probe</b> 15598-L17454  |  | <b>Exon 14</b>  |                 |
| 172 ~       | <b>MEIG1 probe</b> 15599-L17455  |  |                 | <b>Exon 2</b>   |
| 178         | <b>DCLRE1C probe</b> 15600-L17456  |  | <b>Exon 5</b>   |                 |
| 185         | Reference probe 13048-L14231   | 15q                                      |                 |                 |
| 191         | <b>DCLRE1C probe</b> 15601-L17457  |  | <b>Exon 9</b>   |                 |
| 197         | <b>DCLRE1C probe</b> 15602-L18178  |  | <b>Intron 3</b> |                 |
| 203 ~       | <b>MEIG1 probe</b> 15603-L18177  |  |                 | <b>Exon 3</b>   |
| 211 X       | <b>DCLRE1C probe</b> 15604-SP0262-L17460                                       |  | <b>Exon 7</b>   |                 |
| 220         | Reference probe 09570-L10024   | 22q                                      |                 |                 |
| 230         | <b>DCLRE1C probe</b> 15605-L17461  |  | <b>Exon 2</b>   |                 |
| 238         | <b>DCLRE1C probe</b> 15606-L17462  |  | <b>Exon 13</b>  |                 |
| 245         | Reference probe 14028-L15626   | 7q                                       |                 |                 |
| 256         | <b>DCLRE1C probe</b> 16313-L18693  |  | <b>Exon 4</b>   |                 |
| 265         | <b>DCLRE1C probe</b> 15608-L17464  |  | <b>Exon 14</b>  |                 |
| 274         | Reference probe 11696-L17540   | 17q                                      |                 |                 |
| 283         | <b>DCLRE1C probe</b> 15610-L17466  |  | <b>Exon 10</b>  |                 |
| 292         | Reference probe 13325-L14751   | 18q                                      |                 |                 |
| 301         | <b>DCLRE1C probe</b> 15611-L17467  |  | <b>Exon 1</b>   |                 |
| 310         | <b>DCLRE1C probe</b> 15612-L17468  |  | <b>Exon 3</b>   |                 |
| 319         | <b>DCLRE1C probe</b> 15613-L17469  |  | <b>Exon 12</b>  |                 |
| 337         | Reference probe 07824-L07578   | 1q                                       |                 |                 |
| 346         | <b>DCLRE1C probe</b> 15615-L17471  |  | <b>Exon 11</b>  |                 |
| 352         | <b>DCLRE1C probe</b> 15616-L18180  |  | <b>Exon 1</b>   |                 |
| 361         | <b>DCLRE1C probe</b> 15614-L17470  |  | <b>Intron 2</b> |                 |
| 373 ~       | <b>DCLRE1CP probe</b> 16160-L18345   |  |                 | <b>Exon 3</b>   |
| 381         | Reference probe 12579-L11003   | 6p                                       |                 |                 |
| 400 ~       | <b>DCLRE1CP probe</b> 16159-L18344   |  |                 | <b>Intron 2</b> |
| 409         | Reference probe 09615-L09910   | 20p                                      |                 |                 |
| 420         | Reference probe 08839-L08899   | 2p                                       |                 |                 |

<sup>a</sup> See section Exon numbering on page 1 for more information.

X This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

~ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

**Table 2. P368 probes arranged according to chromosomal location**

| Length (nt) | SALSA MLPA probe    | Exon <sup>a</sup>  | Ligation site   | Partial sequence <sup>b</sup> (24 nt adjacent to ligation site) | Distance to next probe |
|-------------|---------------------|--------------------|---|---|------------------------|
|             |                     | <b>DCLRE1C</b>     | <b>NM_001033855.3</b>                                     |   |                        |
|             |                     | <i>stop codon</i>  | 2163-2165 (Exon 14)                                       |   |                        |
| 265         | 15608-L17464        | Exon 14            | 1514-1515   | CTGGGCTCTGTA-CTTCACCTGCAA                                       | 0.2 kb                 |
| 166         | 15598-L17454        | Exon 14            | 1306-1307   | CAAAGTTCATA-CCCGGAACTTT   | 10.5 kb                |
| 238         | 15606-L17462        | Exon 13            | 1179-1180   | CTTCCCAAAGTA-CGGAGCCAAAGT                                       | 3.2 kb                 |
| 319         | 15613-L17469        | Exon 12            | 1096-1097   | TCCTGTGAACGC-ATATCCAAATGT                                       | 3.8 kb                 |
| 346         | 15615-L17471        | Exon 11            | 56 nt after exon 11                                       | TGCAGCCTTGGT-AGTTGACAGGTA                                       | 1.3 kb                 |
| 283         | 15610-L17466        | Exon 10            | 903-904   | AATTACCCTGTG-GAATTACTTCCA                                       | 4.8 kb                 |
| 191         | 15601-L17457        | Exon 9             | 773-774   | TAGGTTTCATGTG-AATAAGCTAGAC                                      | 1.6 kb                 |
| 142 #       | 15595-L17451        | Exon 8             | 636-637   | AGGAGTGTAA-GTGGAGTCTTAG   | 0.1 kb                 |
| 211 Ж #     | 15604-SP0262-L17460 | Exon 7             | 20 and 50 nt after exon 7                                 | GCCTGGAGGAAC-30 nt spanning oligo-TAGTAGGAAGTT                  | 1.0 kb                 |
| 154 #       | 15597-L17453        | Exon 6             | 106 nt before exon 6                                      | ACCATTGTATTG-TAATGGATATGT                                       | 0.9 kb                 |
| 178         | 15600-L17456        | Exon 5             | 425-426   | CTCTTACCAGCT-GGTCACTGTCCG                                       | 3.1 kb                 |
| 256 #       | 16313-L18693        | Exon 4             | 163 nt after exon 4                                       | AAAAATCTAGTC-ATTGAGTGGGGT                                       | 2.7 kb                 |
| 197         | 15602-L18178        | Intron 3           | 2430 nt before exon 4 (NM_001033858.3: 553-554)           | TGGAAGGCAGAA-TGGCCCACTAGG                                       | 2.9 kb                 |
| 310         | 15612-L17468        | Exon 3             | 275-276   | TACTGTTCACCT-GTGACTAAGGAG                                       | 2.3 kb                 |
| 361         | 15614-L17470        | Intron 2           | 1524 nt after exon 2 (NM_001033858.3: 13 nt after exon 3) | TCTGTCTGATTG-CAGAGCCTATGT                                       | 1.6 kb                 |
| 230         | 15605-L17461        | Exon 2             | 50 nt before exon 2                                       | CCCTTGGGCTAA-GGAATCCTCTGG                                       | 4.8 kb                 |
| 352         | 15616-L18180        | Exon 1             | 192-193   | CCCACTGCCACA-AAGGTGAGTGAG                                       | 0.1 kb                 |
| 301         | 15611-L17467        | Exon 1             | 127-128   | GTATCCAACAT-CTCCATAGACCG  | 5.5 kb                 |
|             |                     | <i>start codon</i> | 87-89 (Exon 1)  |   |                        |
|             |                     | <b>MEIG1</b>       | <b>NM_001080836.3</b>                                     |   |                        |
|             |                     | <i>start codon</i> | 200-202 (Exon 2)  |   |                        |
| 148 ~       | 15596-L18176        | Exon 1             | 79-80   | GAGTAGAGAACG-CAAGCACCCACG                                       | 7.1 kb                 |
| 172 ~       | 15599-L17455        | Exon 2             | 308-309   | GGGATGAAACCG-AATATAGACAAG                                       | 6.0 kb                 |
| 203 ~       | 15603-L18177        | Exon 3             | 362-363   | CGGAGACAGGAT-ATGTGAAGAAAC                                       | 45.6 kb                |
|             |                     | <i>stop codon</i>  | 464-466 (Exon 3)  |   |                        |
|             |                     | <b>DCLRE1CP</b>    | <b>NR_144469.1</b>  |   |                        |
| 373 ~#      | 16160-L18345        | Exon 3             | 163-164   | ATAATGGAACCT-TCTTGACACAG  | 2.2 kb                 |
| 400 ~#      | 16159-L18344        | Intron 2           | 1509 nt after exon 2                                      | AATAGTGTGATA-GTGTGTGGCCTT                                       |                        |

<sup>a</sup> See section Exon numbering on page 1 for more information.

<sup>b</sup> Only partial probe sequences are shown. Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com). Please notify us of any mistakes: [info@mrcholland.com](mailto:info@mrcholland.com).

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

~ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

# This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.



## References

- Pannicke U et al. (2010). The most frequent DCLRE1C (ARTEMIS) mutations are based on homologous recombination events. *Hum mutat* 31.2: 197-207.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

## Selected publications using SALSA MLPA Probemix P368 DCLRE1C

- Vignesh P et al. (2021). Clinical, Immunological, and Molecular Features of Severe Combined Immune Deficiency: A Multi-Institutional Experience From India. *Front Immunol*, 11, 3747.

| P368 product history |  |
|----------------------|--|
| Version              | Modification   |
| B1                   | New QDX2 fragments and inclusion of <i>DCLRE1C</i> exon 6 probe. |
| A1                   | First release.   |

| Implemented changes in the product description   |
|--|
| <p>Version B1-02 – 13 December 2021 (04P)</p> <ul style="list-style-type: none"> <li>- Product description rewritten and adapted to a new template.</li> <li>- Ligation sites of the probes targeting the <i>DCLRE1C</i> and <i>MEIG1</i> genes updated according to new version of the NM_ reference sequences.</li> </ul> <p>Version B1-01 – 29 May 2018 (01P)</p> <ul style="list-style-type: none"> <li>- Product description restructured and adapted to a new template.</li> <li>- Product description adapted to a new product version (version number changed, lot number added, changes in Table 1 and Table 2, new picture included).</li> <li>- Warning added to Table 2 for probe specificity relying on a single nucleotide difference between target gene and related gene or pseudogene.</li> </ul> <p>Version 07 – 20 November 2015 (55)</p> <ul style="list-style-type: none"> <li>- Product description adapted to a new lot (lot number added, small changes in Table 1 and Table 2, new picture included).</li> <li>- Exon numbering of the <i>DCLRE1C</i> gene has been changed in Table 1 and 2.</li> <li>- Ligation sites of the probes targeting the <i>DCLRE1C</i> gene updated according to new version of the NM_ reference sequence.</li> </ul> <p>Version 06 (48)</p> <ul style="list-style-type: none"> <li>- Figure(s) based on the use of old MLPA buffer (replaced in December 2012) removed</li> </ul> |

| More information: <a href="http://www.mrcholland.com">www.mrcholland.com</a> ; <a href="http://www.mrcholland.eu">www.mrcholland.eu</a> |   |
|---|---|
|    | MRC Holland bv; Willem Schoutenstraat 1<br>1057 DL, Amsterdam, The Netherlands  |
| E-mail  | <a href="mailto:info@mrcholland.com">info@mrcholland.com</a> (information & technical questions)<br><a href="mailto:order@mrcholland.com">order@mrcholland.com</a> (orders) |
| Phone   | +31 888 657 200   |