

## Product Description

### SALSA® MLPA® Probemix P071-B2 LMNB1-PLP1-NOTCH3

To be used with the MLPA General Protocol.

#### Version B2

For complete product history see page 7.

#### Catalogue numbers:

- **P071-025R:** SALSA MLPA Probemix P071 LMNB1-PLP1-NOTCH3, 25 reactions.
- **P071-050R:** SALSA MLPA Probemix P071 LMNB1-PLP1-NOTCH3, 50 reactions.
- **P071-100R:** SALSA MLPA Probemix P071 LMNB1-PLP1-NOTCH3, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see [www.mrcholland.com](http://www.mrcholland.com)).

#### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at [www.mrcholland.com](http://www.mrcholland.com).

#### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: [www.mrcholland.com](http://www.mrcholland.com). It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

#### General information

The SALSA MLPA Probemix P071 LMNB1-PLP1-NOTCH3 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *LMNB1*, *PLP1* and *NOTCH3* genes, which are associated with autosomal dominant leukodystrophy with autonomic disease (ADLD), Pelizaeus-Merzbacher disease (PMD), and cerebral autosomal dominant arteriopathy with subcortical infarcts and leukoencephalopathy (CADASIL), respectively.

Autosomal dominant leukodystrophy with autonomic disease (ADLD) is a slowly progressive neurological disorder that is characterized by widespread loss of myelin in the central nervous system. Symptoms develop in the fourth or fifth decade of life and include problems with the autonomous nervous system, followed by movement difficulties. ADLD is caused by mutations in the *LMNB1* gene on chromosome 5q23.2, with most cases resulting from a duplication of the *LMNB1* gene.

Pelizaeus-Merzbacher disease (PMD) is a rare X-linked neurological disorder that is caused by demyelination in the central nervous system. The age of onset is typically during infancy or early childhood. PMD is caused by mutations in the *PLP1* gene on chromosome Xq22.2. Duplication of the complete *PLP1* gene accounts for 50-70% of PMD cases, whereas point mutations and deletions are involved in the majority of the remaining cases. In addition to PMD, defects in the *PLP1* gene can also result in spastic paraplegia type 2 (SPG2).

Cerebral autosomal dominant arteriopathy with subcortical infarcts and leukoencephalopathy (CADASIL) is a condition characterized by migraine, strokes, mood disorders, and cognitive decline that can result in dementia. CADASIL has a mid-adult onset and is caused by mutations in the *NOTCH3* gene on chromosome 19p13.12. Defects in the *NOTCH3* gene can also lead to lateral meningocele syndrome.

**This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.**

### Gene structure and transcript variants

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM\_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

### Exon numbering

The *LMNB1*, *PLP1* and *NOTCH3* exon numbering used in this P071-B2 LMNB1-PLP1-NOTCH3 product description is the exon numbering from the NG\_008360.2, NG\_016452.2 and NG\_009819.1 sequence. As changes to the databases can occur after release of this product description, the NM\_ sequence and exon numbering may not be up-to-date.

### Probemix content

The SALSA MLPA Probemix P071-B2 LMNB1-PLP1-NOTCH3 contains 34 MLPA probes with amplification products between 122 and 415 nucleotides. This includes 12 probes for the *LMNB1* gene, one probe for each exon and two probes for exon 1. It also includes seven probes for the *PLP1* gene, one probe for each exon with the exception of exon 1, and three probes for the *NOTCH3* gene, targeting exons 2, 7 and 28. Three flanking probes are included that target regions centromeric and telomeric of the *PLP1* gene. In addition, nine reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online ([www.mrcholland.com](http://www.mrcholland.com)).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at [www.mrcholland.com](http://www.mrcholland.com).

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

### MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

### Required specimens

Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

### Reference samples

A sufficient number ( $\geq 3$ ) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of autosomal dominant leukodystrophy with autonomic disease (ADLD), Pelizaeus-

Merzbacher disease (PMD), and cerebral autosomal dominant arteriopathy with subcortical infarcts and leukoencephalopathy (CADASIL). It is recommended to use samples of the same sex to facilitate interpretation. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

### Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at [www.mrcholland.com](http://www.mrcholland.com). Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

### Interpretation of results

The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$  and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results when **reference samples of the same sex** have been used:

Copy number status		Final ratio (FR)
Autosomal sequences and X chromosome sequences in females	X chromosome sequences in males	
Normal	Normal	$0.80 < FR < 1.20$
Homozygous deletion	Deletion	$FR = 0$
Heterozygous deletion		$0.40 < FR < 0.65$
Heterozygous duplication		$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	Duplication	$1.75 < FR < 2.15$
Ambiguous copy number		All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.

- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

### Limitations of the procedure

- In most populations, the major cause of genetic defects in the *LMNB1*, *PLP1* and *NOTCH3* genes are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P071-B2 LMNB1-PLP1-NOTCH3.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

### Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

### *LMNB1*, *PLP1* and *NOTCH3* mutation databases:

<https://databases.lovd.nl/shared/genes/LMNB1>, <https://databases.lovd.nl/shared/genes/PLP1>, <https://databases.lovd.nl/shared/genes/NOTCH3>. We strongly encourage users to deposit positive results in the Malacards Human Disease Database. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *LMNB1* exons 5 and 7 but not exon 6) to MRC Holland: [info@mrcholland.com](mailto:info@mrcholland.com).

**Table 1. SALSA MLPA Probemix P071-B2 LMNB1-PLP1-NOTCH3**

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) <sup>a</sup>			
		Reference	LMNB1	NOTCH3	PLP1
64-105	Control fragments – see table in probemix content section for more information				
122	Reference probe 02844-L02274	18q			
129 ~	AR probe 01690-L00423				Xq12
136	Reference probe 07292-L06929	6q			
142	<b>PLP1 probe</b> 09208-L11144				<b>Exon 5</b>
151	<b>PLP1 probe</b> 01274-L11199				<b>Exon 2</b>
157	<b>PLP1 probe</b> 09209-L00828				<b>Exon 6</b>
165	Reference probe 08357-L08210	17q			
171	Reference probe 16243-L18501	11p			
177	<b>LMNB1 probe</b> 09221-L09539		<b>Exon 5</b>		
184	<b>LMNB1 probe</b> 09224-L09542		<b>Exon 8</b>		
197	<b>LMNB1 probe</b> 09220-L10868		<b>Exon 4</b>		
206	<b>NOTCH3 probe</b> 09228-L10869			<b>Exon 2</b>	
214	<b>LMNB1 probe</b> 09223-L09541		<b>Exon 7</b>		
220	<b>NOTCH3 probe</b> 09232-L11200			<b>Exon 28</b>	
226	<b>LMNB1 probe</b> 09217-L09535		<b>Exon 1</b>		
232	<b>PLP1 probe</b> 01276-L31608				<b>Exon 3</b>
239 ~	FRMD7 probe 08248-L22041				Xq26
246	<b>LMNB1 probe</b> 09222-L09540		<b>Exon 6</b>		
256	<b>LMNB1 probe</b> 09225-L09543		<b>Exon 9</b>		
264	<b>PLP1 probe</b> 01278-L11204				<b>Exon 4</b>
268	<b>LMNB1 probe</b> 09218-L11205		<b>Exon 2</b>		
283	Reference probe 09101-L09160	4q			
292	<b>LMNB1 probe</b> 09226-L09544		<b>Exon 10</b>		
301	<b>LMNB1 probe</b> 09216-L09534		<b>Exon 1</b>		
310	<b>NOTCH3 probe</b> 09229-L09547			<b>Exon 7</b>	
317	<b>LMNB1 probe</b> 09219-L09537		<b>Exon 3</b>		
328	<b>LMNB1 probe</b> 09227-L09545		<b>Exon 11</b>		
337	Reference probe 09937-L12248	8q			
346	Reference probe 07207-L06857	7p			
355	Reference probe 07236-L06886	1q			
364	<b>PLP1 probe</b> 09185-L00830				<b>Exon 8</b>
373	<b>PLP1 probe</b> 01282-L00965				<b>Exon 7</b>
385 ~	PAK3 probe 02908-L03178				Xq23
415	Reference probe 07594-L29165	15q			

<sup>a</sup> See section Exon numbering on page 2 for more information.

– Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

**Table 2. P071-B2 probes arranged according to chromosomal location****Table 2a. LMNB1**

Length (nt)	SALSA MLPA probe	LMNB1 exon <sup>a</sup>	Ligation site NM_005573.4	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
		Start Codon	374-376 (Exon 1)		
301	09216-L09534	Exon 1	35-36	CAAAGTGCTGCG-AGCAGGAGACGG	0.7 kb
226	09217-L09535	Exon 1	687-688	CAAGCTGCAGAT-CGAGCTGGGCAA	27.1 kb
268	09218-L11205	Exon 2	836-837	CTGCATTGGTG-ACAAAAAAGTT	0.7 kb
317	09219-L09537	Exon 3	925-926	AAAAAACAGTTA-GCAGATGAACT	4.6 kb
197	09220-L10868	Exon 4	1058-1059	CGCGCTTGGTAG-AGGTGGATTCTG	1.6 kb
177	09221-L09539	Exon 5	1225-1226	ATGAATACTTCT-ACTGTCAACAGT	7.2 kb
246	09222-L09540	Exon 6	1403-1404	CAGACAAAGAGA-GAGAGATGGCGG	2.0 kb
214	09223-L09541	Exon 7	1657-1658	GAGGCGAGTAGT-AGTGTAGCATC	1.8 kb
184	09224-L09542	Exon 8	1831-1832	AAATATACCTCA-AGATATGTGCTG	3.2 kb
256	09225-L09543	Exon 9	1963-1964	GTGAAGTTATA-TTGAAAAATTCT	6.7 kb
292	09226-L09544	Exon 10	2046-2047	GGAGGAGGAGGA-AGAAGCAGCTGG	4.1 kb
328	09227-L09545	Exon 11	2695-2696	AGCACTCTGGAT-GATGGATTCCAC	
		Stop Codon	2132-2134 (Exon 11)		

**Table 2b. PLP1**

Length (nt)	SALSA MLPA probe	PLP1 exon <sup>a</sup>	Ligation site NM_001128834.3	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
129 ~	01690-L00423	AR		CATGCAACTCCT-TCAGCAACAGCA	36.2 Mb
		Start Codon	286-288 (Exon 2)		
	No probe	Exon 1			
151	01274-L11199	Exon 2	243-244	AGCCACAAAGCA-GACTAGCCAGCC	8.7 kb
232	01276-L31608	Exon 3	360-361	ACTGGATTGTGT-TTCTTTGGGGTG	1.0 kb
264	01278-L11204	Exon 4	684-685	CATCAAGCTCAT-TCTTTGGAGCGG	1.2 kb
142	09208-L11144	Exon 5	811-812	TGCCTGTGTACA-TTACTTCAACA	0.6 kb
157	09209-L00828	Exon 6	939-940	TTCCCTGGCAAG-GTTTGTGGCTCC	0.9 kb
373	01282-L00965	Exon 7	1013-1014	GTTTATTGCTGC-ATTTGTGGGGGC	1.2 kb
364	09185-L00830	Exon 8	1077-1078	GCCACTTACAAC-TTGCCGTCCTT	7.3 Mb
		Stop Codon	1117-1119 (Exon 8)		
385 ~	02908-L03178	PAK3		CGGGATTCTTCA-GCACTCAACCAC	20.8 Mb
239 ~	08248-L22041	FRMD7		ATTCTGCAGCCA-TTCTGGAAATAA	

**Table 2c. NOTCH3**

Length (nt)	SALSA MLPA probe	NOTCH3 exon <sup>a</sup>	Ligation site NM_000435.3	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
		Start Codon	91-93 (Exon 1)		
206	09228-L10869	Exon 2	252-253	AATGGAGGTCGT-TGCACCCAGCTG	8.1 kb
310	09229-L09547	Exon 7	1155-1156	GATGACGCCTGT-GTCAGCAACCCC	19.3 kb
220	09232-L11200	Exon 28	5257-5258	CAGACTGGATGG-ACACAGAGTGCC	
		Stop Codon	7054-7056 (Exon 33)		

<sup>a</sup> See section

Exon numbering on page 2 for more information.

<sup>b</sup> Only partial probe sequences are shown. Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com). Please notify us of any mistakes: [info@mrcholland.com](mailto:info@mrcholland.com).

~ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.



SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com).

## Related SALSA MLPA probemixes

P022 PLP1 Contains probes for the *PLP1* gene, involved in Pelizaeus-Merzbacher disease.

## References

- Gong Y et al. (2001). LDL receptor-related protein 5 (LRP5) affects bone accrual and eye development. *Cell*, 107(4). 513-523.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

## Selected publications using SALSA MLPA Probemix P071 LMNB1-PLP1-NOTCH3

- Brussino A et al. (2010). A family with autosomal dominant leukodystrophy linked to 5q23.2-q23.3 without lamin B1 mutations. *Eur J Neurol.* 17:541-549.
- Dai Y et al. (2017). An LMNB1 duplication caused adult-onset autosomal dominant leukodystrophy in Chinese family: clinical manifestations, neuroradiology and genetic diagnosis. *Front Mol Neurosci.* 10:215.
- Sandoval-Rodriguez V et al. (2017). Autosomal dominant leukodystrophy presenting as Alzheimer's-type dementia. *Mult Scler Relat Disord.* 17:230-233.
- Westenberger A et al. (2013). X-linked Dystonia-Parkinsonism manifesting in a female patient due to atypical turner syndrome. *Mov Disord.* 28:675-678.

P071 product history	
Version	Modification
B2	One reference probe and one flanking probe have been removed, and one target probe has been changed in length but not in sequence detected.
B1	Three reference probes have been replaced and one NOTCH3 probe has been removed.
A2	The 118 nt Y-chromosome probe has been removed and the control fragments have been replaced by QDX2 control fragments.
A1	First release.

**Implemented changes in the product description**


Version B2-02 – 22 September 2022 (04P)

- Product description rewritten and adapted to a new template.
- Ligation sites of the probes targeting the *PLP1* gene updated according to new version of the NM\_ reference sequence.
- Salt sensitivity warnings have been removed from table 1 and 2.

Version B2-01 – 07 August 2019 (02P)

- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).
- Product description rewritten and adapted to a new template.
- Ligation sites of the probes targeting the *LMNB1* and *NOTCH3* genes updated according to new version of the NM\_ reference sequence.

**More information: [www.mrcholland.com](http://www.mrcholland.com); [www.mrcholland.eu](http://www.mrcholland.eu)**

	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	<a href="mailto:info@mrcholland.com">info@mrcholland.com</a> (information & technical questions) <a href="mailto:order@mrcholland.com">order@mrcholland.com</a> (orders)
Phone	+31 888 657 200