

Product Description

SALSA® MLPA® Probemix P203-B2 PKLR

To be used with the MLPA General Protocol.

Version B2

For complete product history see page 6.

Catalogue numbers:

- **P203-025R:** SALSA MLPA Probemix P203 PKLR, 25 reactions.
- **P203-050R:** SALSA MLPA Probemix P203 PKLR, 50 reactions.
- **P203-100R:** SALSA MLPA Probemix P203 PKLR, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P203 PKLR is a **research use only (RUO)** assay for the detection of deletions or duplications in the *PKLR* gene, which is associated with Pyruvate kinase (PK) deficiency.

PK deficiency is a common cause of hereditary non-spherocytic haemolytic anaemia. It is an autosomal recessive disorder commonly caused by mutations in the gene encoding erythrocyte and liver-type pyruvate kinase (*PKLR*), a key enzyme of the glycolytic pathway. The clinical symptoms of the disease can be variable, ranging from chronic non-spherocytic haemolytic anaemia to neonatal jaundice requiring erythrocyte transfusions.

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering

The *PKLR* exon numbering used in this P203-B2 PKLR product description is the exon numbering from the LRG_1136 sequence. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P203-B2 PKLR contains 21 MLPA probes with amplification products between 177 and 409 nucleotides (nt). This includes twelve probes for the *PKLR* gene, one probe for each exon and two probes for exon 1. In addition, nine reference probes are included that detect autosomal chromosomal

locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of Pyruvate kinase deficiency. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	$FR = 0$
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *PKLR* gene are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P203 PKLR.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

***PKLR* mutation database**

<https://databases.lovd.nl/shared/genes/PKLR>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *PKLR* exons 5 and 7 but not exon 6) to MRC Holland: info@mrcholland.com.

Table 1. SALSA MLPA Probemix P203-B2 PKLR

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a	
		Reference	PKLR
64-105	Control fragments – see table in probemix content section for more information		
177	Reference probe 04359-L03779	7p	
184	PKLR probe 16749-L19376		Exon 8
193	PKLR probe 16750-L19377		Exon 11
202	Reference probe 02498-L01929	17q	
214	Reference probe 10698-L11280	6p	
236	PKLR probe 06512-L20596		Exon 7
244	PKLR probe 06507-L20694		Exon 1
253	Reference probe 06236-L01311	21q	
261	PKLR probe 06508-L20695		Exon 2
274	Reference probe 08053-L07834	5p	
282 ✕	PKLR probe 16664-SP0378-L19230		Exon 1
294	PKLR probe 06510-L22536		Exon 5
320	Reference probe 14346-L16015	2q	
333	PKLR probe 06509-L20595		Exon 3
344	PKLR probe 06514-L23128		Exon 9
352	PKLR probe 16665-L19231		Exon 4
361	Reference probe 08674-L08686	9q	
373	PKLR probe 06511-L22279		Exon 6
384	PKLR probe 06515-L22280		Exon 10
391	Reference probe 01795-L01358	13q	
409	Reference probe 17224-L20551	11p	

^a See section Exon numbering on page 1 for more information.

✕ This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2. PKLR probes arranged according to chromosomal location

Length (nt)	SALSA MLPA probe	PKLR exon ^a	Ligation site NM_000298.6	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	44-46 (Exon 1)		
244	06507-L20694	Exon 1	62-63	TCCAGGAGAACA-TATCATCCCTGC	0.4 kb
282 ± ✕	16664-SP0378-L19230	Exon 1	296 nt and 324 nt after exon 1 (2-3 and 30-31 in NM_181871.3)	CAGTGGTGACAC-28 nt spanning oligo-GCCCAGCTGGGC	0.8 kb
261	06508-L20695	Exon 2	276-277	CCTACTGGACAT-TGACTCCGAGCC	4.4 kb
333	06509-L20595	Exon 3	335-336	TAGGGCCAGCAT-CTCGCTCCGTGG	0.2 kb
352	16665-L19231	Exon 4	479-480	CAGGTTCCCCAC-TCAGCTACCGGC	0.4 kb
294	06510-L22536	Exon 5	713-714	ACGACGGGCTCA-TCTCCCTAGTGG	0.6 kb
373	06511-L22279	Exon 6	981-982	GATCATCAGCAA-AATTGAGAACCA	0.2 kb
236	06512-L20596	Exon 7	1096-1097	AAGGTTTTCTG-GCTCAGAAGATG	0.9 kb
184	16749-L19376	Exon 8	1300-1301	GAAGCGGTGAAG-ATGCAGCATGCG	0.2 kb
344	06514-L23128	Exon 9	1365-1366	GTTTGAGGAGCT-ACGTCGGGCAGC	1.5 kb
384	06515-L22280	Exon 10	1608-1609	TGAACCTCCAGA-AGCCATCTGGGC	1.2 kb
193	16750-L19377	Exon 11	1706-1707	TGGTGATTGTGG-TGACAGGCTGGC	
		<i>stop codon</i>	1766-1768 (Exon 11)		

^a See section Exon numbering on page 1 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

[‡] The 282 nt probe is present in transcript variant 2 (NM_181871.3).

[✕] This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P203 PKLR

- Laas C et al. (2021). Improving the laboratory diagnosis of pyruvate kinase deficiency. *Br J Haematol*, 193(5), 994-1000.
- van Oirschot BA et al. (2014). Novel type of red blood cell pyruvate kinase hyperactivity predicts a remote regulatory locus involved in PKLR gene expression. *Am J Hematol*, 89(4), 380-384.

P203 product history	
Version	Modification
B2	Four reference probes have been replaced.
B1	Two probes for <i>PKLR</i> have been added and two replaced. All reference probes have been replaced and QDX2 fragments have been added.
A1	First release.

Implemented changes in the product description	
Version B2-02 – 22 September 2022 (04P)	
- Product description rewritten and adapted to a new template.	
Version B2-01 – 25 July 2019 (02P)	
- Product description rewritten and adapted to a new template.	
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).	
- Ligation sites of the probes targeting the <i>PKLR</i> gene updated according to new version of the NM_ reference sequence.	
- Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.	

More information: www.mrcholland.com ; www.mrcholland.eu	
	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mrcholland.com (information & technical questions) order@mrcholland.com (orders)
Phone	+31 888 657 200