

Product Description

SALSA® MLPA® Probemix P202-C1 IKZF1-ERG

To be used with the MLPA General Protocol.

Version C1

For complete product history see page 11.

Catalogue numbers:

- **P202-025R:** SALSA MLPA Probemix P202 IKZF1-ERG, 25 reactions.
- **P202-050R:** SALSA MLPA Probemix P202 IKZF1-ERG, 50 reactions.
- **P202-100R:** SALSA MLPA Probemix P202 IKZF1-ERG, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P202 IKZF1-ERG is a **research use only (RUO)** assay for the detection of deletions or duplications in the *IKZF1* (7p12.2) and *ERG* (21q22.2) genes, which are frequently altered in acute lymphoblastic leukemia (ALL). In addition, this probemix can be used to detect copy number aberrations of the *CDKN2A/2B* genes (9p21.3) and the 14q32.33 chromosomal regions, which are frequently altered in ALL.

Partial or complete deletions of the *IKZF1* (also known as *IKAROS*) gene are detected in ALL, especially in cases that also carry the *BCR-ABL1* gene fusion (Philadelphia chromosome). In ALL, *IKZF1* deletions have been associated with relapse and poor clinical outcome (Mullighan et al. 2009, Martinelli et al. 2009, and Iacobucci et al. 2009). Partial or complete gene deletions of *IKZF1* are detected in ~80% of paediatric and 60-90% of adult *BCR-ABL1* positive ALL cases. Partial gene deletions of *IKZF1* frequently affect exons 4-7, but smaller intragenic deletions, down to single exon deletions, have been reported and have been suggested to be associated with unfavourable prognosis in paediatric B-cell precursor (BCP) ALL (Boer et al. 2016).

Short intragenic deletions of *ERG* have been described in BCP-ALL patients and have been shown to be associated with good outcome and, moreover, *ERG* deletion is suggested to define a subgroup of superior outcome among patients with *IKZF1* deletions (Clappier et al. 2014 and Zaliouva et al. 2014).

In chronic-phase chronic myeloid leukemia (CML), *IKZF1* copy number changes are rare. However, in CML-blast crisis, deletions of *IKZF1* are more frequent (25-66%) and might therefore have a role in CML transformation from chronic phase to blast crisis (Alpár et al. 2012).

In common variable immunodeficiency (CVID) disorder, characterized by late-onset hypogammaglobulinemia and a poor antibody response to infectious and vaccine antigens, families with germline heterozygous *IKZF1* deletions have been detected (Kuehn et al. 2016).

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Matched Annotation from NCBI and EMBL-EBI (MANE): <https://www.ncbi.nlm.nih.gov/refseq/MANE/> and <http://tark.ensembl.org/>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering

The *IKZF1*, *ERG*, *CDKN2A* and *CDKN2B* exon numbering used in this P202-C1 IKZF1-ERG product description is the exon numbering from the MANE sequence. The *ERG* and *CDKN2A* exon numbering has been changed; the exon numbering (LRG) used in previous versions of this product description can be found in between brackets in Table 2. **From description version 04 onwards, we have adopted the MANE exon numbering. Please be aware that the MANE and LRG exon numbering do not always correspond, and MANE exon numbering used here may differ from literature.** As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P202-C1 IKZF1-ERG contains 59 MLPA probes with amplification products between 118 and 504 nucleotides (nt). This includes 21 probes for the *IKZF1* gene, 13 probes for the *ERG* gene, three probes for the *CDKN2A/2B* genes, four probes for the 14q32.33 chromosomal region, and for both *IKZF1* and *ERG* a telomeric and a centromeric flanking probe. In addition, 14 reference probes are included that detect relatively copy number stable regions in acute lymphoblastic leukemia. Complete probe sequences are available online (www.mrcholland.com) and the identity of the genes detected by the reference probes is available in Table 3.

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola (2012).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different healthy individuals without a history of leukemia. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID numbers from the Coriell Institute (described in the table below) have been tested with this P202-C1 probemix at MRC Holland and can be used as positive control sample. The quality of cell lines can change; therefore samples should be validated before use.

Coriell sample	Chromosomal position (hg18) of copy number aberration*	Altered target genes in P202-C1	Expected copy number alteration
NA10925	7p12.2	<i>ZBPB</i> , <i>IKZF1</i> and <i>FIGNL1</i>	Heterozygous deletion
NA07081	7p12.2	<i>ZBPB</i> , <i>IKZF1</i> and <i>FIGNL1</i>	Heterozygous duplication
NA01750	9p21.3	<i>CDKN2A</i> and <i>CDKN2B</i>	Heterozygous duplication
NA08123	14q32.33	<i>CEP170B</i> , <i>MTA1</i> , <i>CRIP2</i> and <i>IGHD</i>	Heterozygous deletion
NA09868	21q22.13-q22.2	<i>KCNJ6</i> , <i>ERG</i> and <i>ETS2</i>	Heterozygous deletion

* Indicated chromosomal bands accommodate genes targeted by MLPA probes, however, the whole extent of copy number alteration (CNA) present in this cell line cannot be determined by this P202-C1 IKZF1-ERG probemix.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 . When this criterion is fulfilled, the following cut-off values for the final ratio (FR) of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	$FR = 0$
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

Please note that these above mentioned final ratios are only valid for germline testing. Final ratios are affected both by percentage of tumour cells and by possible subclonality.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in subclonal cases.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *IKZF1* gene. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

P202 specific note

In samples from tumour tissues, reference probes are more prone to have deviating copy number results as compared to blood derived germline samples. When regions targeted by reference probes are affected by copy number alterations, it can help to turn the slope correction off in Coffalyser.Net analysis to get the correct copy number interpretation on the target region.

Limitations of the procedure

- Most genetic alterations in cancer are small (point) mutations. If present, these type of mutations in *IKZF1*, *CDKN2A/2B* and *ERG* will not be detected by using SALSA MLPA Probemix P202 IKZF1-ERG.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.
- MLPA analysis on tumour samples provides information on the average situation in the cells from which the DNA sample was purified. Gains or losses of genomic regions or genes may not be detected if the percentage of tumour cells is low. In addition, subclonality of the aberration affects the final ratio of the corresponding probe. Furthermore, there is always a possibility that one or more reference probes *do* show a copy number alteration in a patient sample, especially in samples with more chaotic karyotypes.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

COSMIC mutation database

<http://cancer.sanger.ac.uk/cosmic>. We strongly encourage users to deposit positive results in the COSMIC. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report false positive results due to SNVs and unusual results (e.g., a deletion of *IKZF1* exons 3 and 5 but not exon 4) to MRC Holland: info@mrcholland.com.

Table 1. SALSA MLPA Probemix P202-C1 IKZF1-ERG

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a			
		Reference	IKZF1	ERG	Other regions
64-105	Control fragments – see table in probemix content section for more information				
118	Reference probe 20126-L20708	11p11			
124 ~	KCNJ6 probe 21195-L31204			21q22.13	
129 «	IKZF1 probe 21892-L30651		Exon 1		
136	Reference probe 13224-L31079	1p21			
142	IKZF1 probe 21511-L30640		Exon 5		
148 «	IKZF1 probe 17215-L31078		Exon 1		
154	IKZF1 probe 15407-L17580		Exon 8		
160 ~	ZBPB probe 15408-L17236		7p12.2		
166	Reference probe 07394-L20506	12q13			
172	IKZF1 probe 15410-L31063		Exon 7		
177	ERG probe 18151-L22661			upstream	
184	ERG probe 18152-L22662			Exon 10	
190	IKZF1 probe 15409-L30639		Exon 6		
196	Reference probe 07815-L30637	3p22			
202	IKZF1 probe 15424-L17583		Exon 4		
207 Δ	IGHD probe 16524-L28592				14q32.33
213 «	IKZF1 probe 14056-L20508		Exon 2		
220	Reference probe 08940-L20509	11p15			
226 ~	FIGNL1 probe 20418-L28594		7p12.2		
231	IKZF1 probe 15416-L17244		Exon 3		
236	IKZF1 probe 21893-L30652		Exon 3		
244	Reference probe 18056-L22446	16q23			
250	ERG probe 18155-L22665			Exon 7	
256	IKZF1 probe 15426-L17587		Exon 6		
262	CDKN2A probe 15675-L18954				9p21.3
269 «	IKZF1 probe 13877-L15918		Exon 1		
275	Reference probe 16270-L30644	20q11			
283	ERG probe 21887-L31064			upstream	
288	IKZF1 probe 17109-L20256		Exon 8		
296 ~	ETS2 probe 09515-L31066			21q22.2	
301	ERG probe 21894-L31065			Exon 3	
308 «	IKZF1 probe 21895-L30654		Exon 2		
315	ERG probe 21885-L31210			Exon 4	
322	CDKN2A probe 21890-L22800				9p21.3
329	Reference probe 03918-L30417	15q21			
336	ERG probe 02833-L30411			Exon 2	
343	IKZF1 probe 13869-L30647		Exon 7		
352	ERG probe 20881-L22666			Exon 1	
358	IKZF1 probe 21896-L30655		Exon 5		
366 «	CEP170B probe 21897-L30656				14q32.33
373	ERG probe 18157-L31067			Exon 5	
382	ERG probe 18158-L22668			Exon 8	
389	Reference probe 08835-L30638	2p13			
395 « ∅	IKZF1 probe 21883-L17250		Intron 1		
402 ∅	IKZF1 probe 21898-L31074		Intron 3		
409	CDKN2B probe 03814-L03851				9p21.3
417	Reference probe 03073-L15904	5p15			
423	IKZF1 probe 21886-L31075		Exon 4		
427 «	MTA1 probe 14071-L31076				14q32.33
436	Reference probe 13809-L31211	5q14			
443	Reference probe 12790-L31073	2q13			

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a			
		Reference	<i>IKZF1</i>	<i>ERG</i>	Other regions
449 «	CRIP2 probe 21901-L31116				14q32.33
459 ð	IKZF1 probe 21903-L31254		Intron 3		
468	Reference probe 13538-L31070	19p13			
475	ERG probe 18159-L22669			Exon 9	
481	ERG probe 18160-L22670			upstream	
490	IKZF1 probe 21904-L31069		upstream		
497	ERG probe 19022-L25058			Exon 6	
504	Reference probe 15203-L22928	3p12			

^a See section Exon numbering on page 2 for more information.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

– Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

∅ This probe is targeting an alternative exon of *IKZF1* present in NM_001291837.2 transcript variant 14.

ð This probe is targeting an alternative exon of *IKZF1* present in NM_001291845.2 transcript variant 15.

∩ This probe is targeting an alternative exon of *IKZF1* present in NM_001291846.2 transcript variant 16.

Δ More variable. This probe may be more variable, as a high number of variations is identified in healthy control samples, see further details in DGV database. Aberrant results should be treated with caution.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2a. *IKZF1* probes arranged according to chromosomal location

Length (nt)	SALSA MLPA probe	Gene, exon ^a	Chromosomal band (hg18) / Ligation site	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
<i>IKZF1</i> gene, 7p12.2					
Indicated ligation sites are in NM_006060.6, unless otherwise noted.					
160 ~	15408-L17236	<i>ZBPB</i>	7p12.2	CCACACGTGTTA-TGTGTAACGCAA	193.0 kb
490	21904-L31069	<i>IKZF1</i> , up	30 kb upstream of ex 1	TCATGTTCACAA-AATCTTGGGCAT	29.3 kb
129 «	21892-L30651	<i>IKZF1</i> , ex 1	539 nt before ex 1 reverse; NM_001291839.2: 63 nt after exon 1 reverse	GAAAACCTTTCGCA-ATCGCCCGGGC	0.2 kb
148 «	17215-L31078	<i>IKZF1</i> , ex 1	316 nt before ex 1	GCAGGTCGAGCA-GGGACCGCCAGC	0.5 kb
269 «	13877-L15918	<i>IKZF1</i> , ex 1	187-188	TCTTGGCCCCAA-AGCGCGACGCAC	4.0 kb
395 « ∅	21883-L17250	<i>IKZF1</i> , intr 1	3.9 kb after ex 1; NM_001291837.2; 172-173	TGAAAAGGCAG-CTCTCACTTGGC	10.1 kb
308 «	21895-L30654	<i>IKZF1</i> , ex 2	13 nt before ex 2 reverse	TGAGAAAGAGAG-GAAGGGATTTTA	0.1 kb
213 «	14056-L20508	<i>IKZF1</i> , ex 2	253-254	AGACATGTCCCA-AGTTTCAGGTGA	8.6 kb
231	15416-L17244	<i>IKZF1</i> , ex 3	288-289	CTGTAAGCGATA-CTCCAGATGAGG	0.1 kb
236	21893-L30652	<i>IKZF1</i> , ex 3	362-361 reverse	CTGTCACTCTTG-GAGCTTGTCTGT	65.4 kb
459 ð	21903-L31254	<i>IKZF1</i> , intr 3	8.4 kb before ex 4; NM_001291845.2: 498-499	TCCTATCATGTA-AATATCGTACGT	1.4 kb
402 ∩	21898-L31074	<i>IKZF1</i> , intr 3	7.0 kb before ex 4; NM_001291846.2; 469-470	ATCTTCTCACAC-AAGCGGCTACTT	7.0 kb
202	15424-L17583	<i>IKZF1</i> , ex 4	0 nt before ex 4	TGTTTCTTTCAG-CCAGTAATGTTA	0.2 kb
423	21886-L31075	<i>IKZF1</i> , ex 4	590-591	GATATCTGTGGG-ATCATTTCATC	5.8 kb
142	21511-L30640	<i>IKZF1</i> , ex 5	680-681	TGCGGGCCCTCA-TTCACCCAGAAG	0.1 kb
358	21896-L30655	<i>IKZF1</i> , ex 5	797-796 reverse	GAGTGCCTCCTC-AGGTGGCCAGTG	4.7 kb
256	15426-L17587	<i>IKZF1</i> , ex 6	812-813	TTTTCTGCAGTT-GGTAACCTCAC	0.1 kb
190	15409-L30639	<i>IKZF1</i> , ex 6	897-898	GCTGCCACAACCT-ACTTGAAAGCA	4.3 kb
172	15410-L31063	<i>IKZF1</i> , ex 7	26 nt before ex 7	AAGCCTTTCTAA-ACTGGCCTCTCT	0.1 kb
343	13869-L30647	<i>IKZF1</i> , ex 7	994-995	CAAGATAGGATC-AGAGAGATCTCT	8.5 kb
154	15407-L17580	<i>IKZF1</i> , ex 8	1441-1442	CAACGAGGAGCA-GCGCAGCGGTCT	0.8 kb
288	17109-L20256	<i>IKZF1</i> , ex 8	2271-2272	GGTGTGCCGCCA-CCCAAGTGCCAA	44.9 kb

Length (nt)	SALSA MLPA probe	Gene, exon ^a	Chromosomal band (hg18) / Ligation site	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
226 ~	20418-L28594	<i>FIGNL1</i>	7p12.2	AAAGCCACCATA-AAGGAAATAGTT	-

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

~ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

∅ This probe is targeting an alternative exon of *IKZF1* present in NM_001291837.2 transcript variant 14.

∂ This probe is targeting an alternative exon of *IKZF1* present in NM_001291845.2 transcript variant 15.

∪ This probe is targeting an alternative exon of *IKZF1* present in NM_001291846.2 transcript variant 16.

Clinical and/or diagnostic significance of copy number alterations in the above mentioned three alternative exons is not yet established.

Table 2b. *CDKN2A/2B* probes arranged according to chromosomal location

Length (nt)	SALSA MLPA probe	Gene, exon ^a	Ligation site	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
<i>CDKN2A</i> and <i>CDKN2B</i> genes, 9p21.3					
Indicated ligation sites for <i>CDKN2A</i> are in NM_000077.5 and NM_058195.4. Indicated ligation site for <i>CDKN2B</i> is for NM_004936.4.					
262	15675-L18954	<i>CDKN2A</i> , ex 3 (4)	NM_000077.5; 182 nt after ex 3; NM_058195.4; 182 nt after ex 3	TGAAATGCGGTT-AAAATGATGAAT	7.4 kb
322	21890-L22800	<i>CDKN2A</i> , ex 1 (2)	NM_000077.5; 71 nt before ex 1; NM_058195.4; 3.7 kb before ex 2	GCACCGGAGGAA-GAAAGAGGAGGG	33.9 kb
409	03814-L03851	<i>CDKN2B</i> , ex 1	NM_004936.4; 462-463	CCTGGAAGCCGG-CGCGGATCCCAA	-

Table 2c. 14q32.33 probes arranged according to chromosomal location

Length (nt)	SALSA MLPA probe	Gene	Chromosomal band (hg18)	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
14q32.33 region					
Chromosome 14, which contains immunoglobulin heavy locus (<i>IGH</i>), is frequently trisomic in high hyperdiploid B-cell ALL. This can lead to a higher number of <i>IGH</i> rearrangements than in cases with disomy 14. Four probes are included in this probemix targeting the 14q32.33 chromosomal region. Information about 14q32.33 copy number is suggested to be used in combination with <i>IGH</i> rearrangements for minimal residual disease detection in BCP-ALL (Csinady et al. 2009).					
366 «	21897-L30656	<i>CEP170B</i>	14q32.33	CCCCTGAACCTCT-CCAGGGCATCTT	579.0 kb
427 «	14071-L31076	<i>MTA1</i>	14q32.33	CACACAGTCTTA-CCAGTGGTATTC	9.0 kb
449 «	21901-L31116	<i>CRIP2</i>	14q32.33	CACCCTGCAGCC-ACTGCCATTTC	372.2 kb
207 Δ	16524-L28592	<i>IGHD</i>	14q32.33	TCCGTGACTGTC-ACCTGGTACATG	-

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Δ More variable. This probe may be more variable, as high number of variation is identified in healthy control samples, see further details in DGV database. Aberrant results should be treated with caution.

Table 2d. *ERG* probes arranged according to chromosomal location

Length (nt)	SALSA MLPA probe	Gene, exon ^a	Chromosomal band (hg18) / Ligation site	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
<i>ERG</i> gene, 21q22.2					
Indicated ligation sites are in NM_182918.4, unless otherwise noted.					
124 ~	21195-L31204	<i>KCNJ6</i>	21q22.13	AGCTCCTACATC-ACCAGTGAGATC	757.6 kb
184	18152-L22662	<i>ERG</i> , ex 10 (12)	1462-1463	TACTGGAATTCA-CCAACTGGGGGT	7.5 kb
475	18159-L22669	<i>ERG</i> , ex 9 (11)	991-992	ATTCTTGGACCA-ACAAGTAGCCGC	0.7 kb
382	18158-L22668	<i>ERG</i> , ex 8 (10)	931-932	CCATCTCCTTCC-ACAGTGCCCAA	0.7 kb
250	18155-L22665	<i>ERG</i> , ex 7 (9)	905-906	CCACGCCCCAGT-CGAAAGGTACAG	8.2 kb
497	19022-L25058	<i>ERG</i> , ex 6 (8)	791-792	CTTTTATTTTCC-CAAATACTTCAG	2.0 kb
373	18157-L31067	<i>ERG</i> , ex 5 (7)	702-703	TCCTCTCCACA-TTTGACTTCAGA	0.9 kb
315	21885-L31210	<i>ERG</i> , ex 4 (6)	632-633	ACTTCCAGAGGC-TCACCCAGCT	19.9 kb

Length (nt)	SALSA MLPA probe	Gene, exon ^a	Chromosomal band (hg18) / Ligation site	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
301	21894-L31065	<i>ERG</i> , ex 3 (5)	408-407 reverse	TGTAGCTGCCGT-AGTTCATCCCAA	22.1 kb
336	02833-L30411	<i>ERG</i> , ex 2 (4)	152-153	ACCAGTCGTTGT-TTGTGTGCCT	52.9 kb
352	20881-L22666	<i>ERG</i> , ex 1 (intr 3)	23 nt before exon 1	TGGCTGACTTCA-TTTCCAGACTT	77.2 kb
283	21887-L31064	<i>ERG</i> , up (3)	77 kb before exon 1; NM_001136154.1; 270-271	GCTTACTGAAGG-ACATGATTGAGA	9.2 kb
481	18160-L22670	<i>ERG</i> , up (2)	86 kb before exon 1; NM_001136154.1; 169-170	CGTGTGACCAA-AAGCAAGACAAA	76.8 kb
177	18151-L22661	<i>ERG</i> , up (1)	163 kb before exon 1; NM_001136154.1; 76-77	CGCTCCGGGACG-GTCGTGACGGCC	151.4 kb
296 -	09515-L31066	<i>ETS2</i>	21q22.2	AAATGAAGAGCA-AACACTGCAAGA	-

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

^a See section Exon numbering on page 2 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 3. Reference probes arranged according to chromosomal location

Length (nt)	SALSA MLPA probe	Gene	Chromosomal band (hg18)	Partial sequence (24 nt adjacent to ligation site)	Distance to next probe
136	13224-L31079	<i>COL11A1</i>	1p21	CAGATGGTGTCA-GAGGTCTCAAGG	-
389	08835-L30638	<i>DYSF</i>	2p13	GACTGAGAGCAA-AATCCCAGCACG	37.1 Mb
443	12790-L31073	<i>EDAR</i>	2q13	AGAATCAAGGCT-TTTGTGATATGT	-
196	07815-L30637	<i>SCN5A</i>	3p22	TGGTTCGAGACA-TTCATCATCTTC	43.1 Mb
504	15203-L22928	<i>GBE1</i>	3p12	GACCTAGAGGGA-CTCATGATCTTT	-
417	03073-L15904	<i>CTNND2</i>	5p15	CATCAGCCTCAG-AGAAGACGAGTT	78.3 Mb
436	13809-L31211	<i>ADGRV1</i>	5q14	ATGCGAGACGAA-CAGTCTGCAGTC	-
220	08940-L20509	<i>SLC6A5</i>	11p15	TTGCCTCTCAGG-TGTGGAAAGATG	26.7 Mb
118	20126-L20708	<i>MYBPC3</i>	11p11	ACGTCTCTGACA-CCACGGTCTCCC	-
166	07394-L20506	<i>COL2A1</i>	12q13	TCACCTCCTTCT-TGCTCACAGGGT	-
329	03918-L30417	<i>FBN1</i>	15q21	CCTACAGATGTG-AATGCTTCCCTG	-
244	18056-L22446	<i>PLCG2</i>	16q23	GATCCAGCAGTA-CTTCCCATCCAA	-
468	13538-L31070	<i>CACNA1A</i>	19p13	TGTGCAGTCCTT-CAAGGTGAGTCC	-
275	16270-L30644	<i>SAMHD1</i>	20q11	TGACGACATGGA-AGCCTATACTAA	-

Complete probe sequences are available at www.mrcholland.com.

Related SALSA MLPA probemixes

- **P335 ALL-*IKZF1***: Contains probes for *IKZF1*, *CDKN2A/2B*, *EBF1*, PAR1-region, *PAX5*, *ETV6*, *BTG1* and *RB1*.
- **P327 *iAMP21-ERG***: Contains probes for *RUNX1*, *ERG* and *iAMP21* detection in ALL.
- **P383 T-ALL**: Contains probes for the *STIL-TAL1*, *LEF1*, *CASP8AP2*, *MYB*, *EZH2*, *CDKN2A/B*, *MTAP*, *MLL3*, *NUP214-ABL1*, *PTEN*, *LMO1*, *LMO2*, *NF1*, *SUZ12*, *PTPN2* and *PHF6* genes involved in T-ALL.
- **P329 *CRLF2-CSF2RA-IL3RA***: Contains probes for the *CRLF2*, *CSF2RA* and *IL3RA* genes involved in ALL.
- **P377 Hematologic Malignancies**: Contains probes for screening DNA samples on the most common copy number changes associated with ALL, AML, CLL, CML, MDS and lymphomas.
- **P419 *CDKN2A/2B-CDK4***: Contains several probes for each exon of *CDKN2A/2B* and at least one probe for each exon of *CDK4*.
- **P037 *CCL-1*, P038 *CLL-2* and P040 *CLL***: Contain probes for various genes involved in CLL.
- **ME024 9p21 *CDKN2A/2B* region**: Contains probes for detection of methylation and/or copy number status of the chromosomal region 9p21.3 (*CDKN2A/2B*, *CDKN2B-AS1*, *MTAP*, *MIR31* and *PAX5*).

References

- Alpar D et al. (2012). MLPA is a powerful tool for detecting lymphoblastic transformation in chronic myeloid leukemia and revealing the clonal origin of relapse in pediatric acute lymphoblastic leukemia. *Cancer Genet.* 205:465-9.
- Boer JM et al. (2016). Prognostic value of rare IKZF1 deletion in childhood B-cell precursor acute lymphoblastic leukemia: an international collaborative study. *Leukemia.* 30:32-8.
- Clappier E et al. (2014). An intragenic ERG deletion is a marker of an oncogenic subtype of B-cell precursor acute lymphoblastic leukemia with a favorable outcome despite frequent IKZF1 deletions. *Leukemia.* 28:70-7.
- Csinady E et al. (2009). Chromosome 14 copy number-dependent IGH gene rearrangement patterns in high hyperdiploid childhood B-cell precursor ALL: implications for leukemia biology and minimal residual disease analysis. *Leukemia.* 23:870-6.
- Hömig-Hölzel C and Savola S (2012). Multiplex ligation-dependent probe amplification (MLPA) in tumor diagnostics and prognostics. *Diagn Mol Pathol.* 21:189-206.
- Iacobucci I et al. (2009). Identification and molecular characterization of recurrent genomic deletions on 7p12 in the IKZF1 gene in a large cohort of BCR-ABL1-positive acute lymphoblastic leukemia patients: on behalf of Gruppo Italiano Malattie Ematologiche dell'Adulto Acute Leukemia Working Party (GIMEMA AL WP). *Blood.* 114:2159-67.
- Kuehn HS et al. (2016). Loss of B cells in patients with heterozygous mutations in IKAROS. *N Engl J Med.* 374:1032-43.
- Martinelli G et al. (2009). IKZF1 (Ikaros) deletions in BCR-ABL1-positive acute lymphoblastic leukemia are associated with short disease-free survival and high rate of cumulative incidence of relapse: a GIMEMA AL WP report. *J Clin Oncol.* 27:5202-7.
- Mullighan CG et al. (2009). Deletion of IKZF1 and prognosis in acute lymphoblastic leukemia. *N Engl J Med.* 360:470-80.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.
- Zaliova M et al. (2014). ERG deletion is associated with CD2 and attenuates the negative impact of IKZF1 deletion in childhood acute lymphoblastic leukemia. *Leukemia.* 28:182-5.

Selected publications using SALSA MLPA Probemix P202 IKZF1-ERG

- Alpar D et al. (2012) MLPA is a powerful tool for detecting lymphoblastic transformation in chronic myeloid leukemia and revealing the clonal origin of relapse in pediatric acute lymphoblastic leukemia. *Cancer Genet.* 205:465-9.
- Boer JM et al. (2015). Expression profiling of adult acute lymphoblastic leukemia identifies a BCR-ABL1-like subgroup characterized by high non-response and relapse rates. *Haematologica.* 100(7):e261-4.
- Boer JM et al. (2016). Prognostic value of rare IKZF1 deletion in childhood B-cell precursor acute lymphoblastic leukemia: an international collaborative study. *Leukemia.* 30:32-8.
- Buitenkamp TD et al. (2012). Outcome in children with Down's syndrome and acute lymphoblastic leukemia: role of IKZF1 deletions and CRLF2 aberrations. *Leukemia.* 26:2204-11.
- Clappier E et al. (2014). An intragenic ERG deletion is a marker of an oncogenic subtype of B-cell precursor acute lymphoblastic leukemia with a favorable outcome despite frequent IKZF1 deletions. *Leukemia.* 28:70-7.
- Clappier E et al. (2015). IKZF1 deletion is an independent prognostic marker in childhood B-cell precursor acute lymphoblastic leukemia, and distinguishes patients benefiting from pulses during maintenance therapy: results of the EORTC Children's Leukemia Group study 58951. *Leukemia.* 29:2154-61.
- Dorge P et al. (2013). IKZF1 deletion is an independent predictor of outcome in pediatric acute lymphoblastic leukemia treated according to the ALL-BFM 2000 protocol. *Haematologica.* 98:428-32.

- Kim M et al. (2015). Impact of IKZF1 deletions on long-term outcomes of allo-SCT following imatinib-based chemotherapy in adult Philadelphia chromosome-positive ALL. *Bone Marrow Transplant*. 50:354-62.
- Kuehn HS et al. (2016). Loss of B cells in patients with heterozygous mutations in IKAROS. *N Engl J Med*. 374:1032-43.
- Lopes BA et al. (2019). IKZF1 Deletions with COBL Breakpoints Are Not Driven by RAG-Mediated Recombination Events in Acute Lymphoblastic Leukemia. *Transl Oncol*. 12:726-32.
- Maciel ALT et al. (2022). IKZF1 deletions associate with CRLF2 overexpression leading to a poor prognosis in B-cell precursor acute lymphoblastic leukaemia. *Transl Oncol*. 15:101291.
- Öfverholm T et al. (2013). Impact of IKZF1 deletions and PAX5 amplifications in pediatric B-cell precursor ALL treated according to NOPHO protocols. *Leukemia*. 27:1936-9.
- Palmi C et al. (2013). What is the relevance of Ikaros gene deletions as a prognostic marker in pediatric Philadelphia-negative B-cell precursor acute lymphoblastic leukemia? *Haematologica*. 98:1226-31.
- Ribera J et al. (2019). Molecular profiling refines minimal residual disease-based prognostic assessment in adults with Philadelphia chromosome-negative B-cell precursor acute lymphoblastic leukemia. *Genes Chromosomes Cancer*. 58:815-9.
- Schwab CJ et al. (2013). Genes commonly deleted in childhood B-cell precursor acute lymphoblastic leukemia: association with cytogenetics and clinical features. *Haematologica*. 98:1081-8.
- Steeghs EMP et al. (2019). Copy number alterations in B-cell development genes, drug resistance, and clinical outcome in pediatric B-cell precursor acute lymphoblastic leukemia. *Sci Rep*. 9:4634.
- van der Veer A et al. (2013). Independent prognostic value of BCR-ABL1-like signature and IKZF1 deletion, but not high CRLF2 expression, in children with B-cell precursor ALL. *Blood*. 122:2622-9.
- van der Veer A et al. (2014). IKZF1 status as a prognostic feature in BCR-ABL1-positive childhood ALL. *Blood*. 123:1691-8.

P202 product history	
Version	Modification
C1	Probemix content completely revised. Probes have been added for the <i>ERG</i> gene and flanking regions, and removed for the <i>IKZF2</i> and <i>IKZF3</i> genes. In addition, several probes for <i>IKZF1</i> and several reference probes have been replaced.
B2	A few probes have a change in length but no change in sequences detected.
B1	Two new <i>IKZF1</i> probes, three new <i>IKZF2</i> probes, one new <i>IKZF3</i> probe and four probes for 14q32.33 added. New QDX2 fragments added.
A1	First release.

Implemented changes in the product description
<p>Version C1-04 – 18 October 2022 (04P)</p> <ul style="list-style-type: none"> - Product description rewritten and adapted to a new template. - Various minor textual or layout changes. - Exon numbering of the <i>ERG</i> and <i>CDKN2A</i> genes has been changed according to MANE. - Ligation sites of the probes targeting the <i>CDKN2A/2B</i> and <i>ERG</i> genes updated according to new version of the NM_ reference sequence. - Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products. <p>Version C1-03 – 21 April 2021 (01P)</p> <ul style="list-style-type: none"> - Warning added to Tables 1 & 2c about possible higher variability of the IGHD probe (16524-L28592) at 207 nt. <p>Version C1-02 – 05 February 2020 (01P)</p> <ul style="list-style-type: none"> - Gene name GPR98 has been changed to ADGRV1 in Table 2e, following HUGO Gene Nomenclature. - Various minor layout changes. - New references added in Selected publications on page 9.

Version C1-01 – 16 January 2019 (01P)

- Product description adapted to a new product version (version number changed, changes in Table 1 and Tables 2) and to a new template.
- For uniformity, the chromosomal positions and bands in this document are now all based on hg18 (NCBI36).

Version 10 – 01 February 2018 (T08)

- Publication information of the references updated on page 2.
- Exon and ligation site information corrected for IKZF1 probes (17215-L20539, 15422-L17250 and 15424-L17583).
- Shortened version of the chromosomal bands introduced for reference probes in Table 1 and 2.
- Typos corrected.

Version 09 – 15 September 2017 (T08)

- Warning added in Table 1, 296 nt probe 04606-L17670, 373 nt probe 13657-L15111, 427 nt probe 14071-L18749, and 445 nt probe 16357-L15668.

More information: www.mrcholland.com; www.mrcholland.eu

	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mrcholland.com (information & technical questions) order@mrcholland.com (orders)
Phone	+31 888 657 200