

# Product Description

## SALSA® MLPA® Probemix P335-C2 ALL-IKZF1

To be used with the MLPA General Protocol.

### Version C2

As compared to version C1, lengths of several probes are adjusted but no changes in the sequences detected. For complete product history see page 15.

### Catalogue numbers:

- **P335-025R:** SALSA MLPA Probemix P335 ALL-IKZF1, 25 reactions.
- **P335-050R:** SALSA MLPA Probemix P335 ALL-IKZF1, 50 reactions.
- **P335-100R:** SALSA MLPA Probemix P335 ALL-IKZF1, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see [www.mrcholland.com](http://www.mrcholland.com)).

### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at [www.mrcholland.com](http://www.mrcholland.com).

### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: [www.mrcholland.com](http://www.mrcholland.com). It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

### Intended purpose

The SALSA MLPA Probemix P335 ALL-IKZF1 is an in vitro diagnostic (IVD)<sup>1</sup> or a research use only (RUO) semi-quantitative assay<sup>2</sup> for the detection of deletions of the *IKZF1* gene for stratification of patients with acute lymphoblastic leukemia (ALL) into prognostic subgroups. The SALSA MLPA Probemix P335 ALL-IKZF1 is a RUO assay<sup>2</sup> for the detection of deletions or duplications in B-cell differentiation and cell cycle control genes (*EBF1*, *CDKN2A/B*, *PAX5*, *ETV6*, *BTG1* and *RB1*) and in the *PAR1* region. This assay is for use on genomic DNA isolated from human peripheral whole blood and bone marrow specimens.

Copy number alterations (CNAs) detected with P335 ALL-IKZF1 should be confirmed with a different technique. In particular, CNAs detected by only a single probe always require confirmation by another method. In the majority of patients, defects in the *IKZF1* gene are deletions, but point mutations can occur which will not be detected by MLPA. It is therefore recommended to use this assay in combination with sequence analysis of the *IKZF1* gene.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes or population screening.

<sup>1</sup> Please note that this probemix is for in vitro diagnostic (IVD) use in the countries specified at the end of this product description. In all other countries, the product is for research use only (RUO).

<sup>2</sup> To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

Note that the clinical relevance of some genes in the P335 ALL-IKZF1 probemix is not yet fully established. Therefore, the **CE mark for diagnostic use only applies to the *IKZF1* gene**, and all other genes in the probemix are meant to be used in a research setting only.

## Clinical background

The overall incidence rate of ALL amounts to 1.6 in 100,000 per year (Malard and Mohty 2020). The peak incidence lies in childhood at an age of less than 5 years; thereafter, the incidence rate declines continually until the age of 50 years. After that, incidence slightly rises a second time. In patients over 50 years it rises a second time and reaches another peak at the age of over 80 years). There is a slight predominance of males (1.2:1). B-cell ALL accounts for 75% of all cases of ALL and T-cell ALL accounts for the remaining 25% of cases (<https://www.lls.org/leukemia/acute-lymphoblastic-leukemia/diagnosis/all-subtypes>).

Partial or complete deletions of the *IKZF1* (IKAROS family zinc finger 1) gene are frequently detected in ALL cases (Mullighan et al. 2008), especially in those patients who also carry the *BCR-ABL1* gene fusion (Philadelphia chromosome). *IKZF1* deletions can be identified in approximately 70% of the children with Philadelphia chromosome-positive (Ph+) ALL (2-4% of all paediatric ALL cases), in 10-15% of Philadelphia chromosome-negative (Ph-) paediatric ALL, and in 40% of adult B-cell precursor ALL cases (Bernt and Hunger 2014; Lejman et al. 2022; van der Sligte et al. 2015). Deletion of *IKZF1* is associated with a poor prognosis in B-ALL patients (Mullighan et al. 2009a) and a higher chance of relapse (Kuiper et al. 2010).

Several other (partial) gene deletions and duplications, such as in *PAX5*, *ETV6*, *RB1*, *BTG1*, *EBF1* and *CDKN2A/2B*, have also been described in ALL patients. Prognostic profiles combining these aberrations, such as the *IKZF1*<sup>plus</sup> profile (*IKZF1* deletions that co-occur with deletions in *CDKN2A*, *CDKN2B*, *PAX* or *PAR1* in the absence of *ERG* deletion; Stanulla et al. 2018), have been described in recent years. See Table 2 for more information on all genes and regions covered in this probemix.

## Gene structure

The *IKZF1* gene, located on chromosome 7p12.2, spans ~86 kb of genomic DNA and contains 8 exons. The *IKZF1* MANE select transcript sequence is available at <https://www.ncbi.nlm.nih.gov/refseq/MANE/> and is identical to NM\_006060.6.

## Transcript variants

For *IKZF1*, the longest isoform, also known as Ik-1 as described by Ezzat et al. (2003) is encoded by transcript variant NM\_006060.6 (6255 nt; coding sequence 222-1781; <http://www.ncbi.nlm.nih.gov/gene/10320>). This is the MANE Select transcript sequence and also a reference standard in the NCBI RefSeq project. The ATG translation start site is located in exon 2 and the stop codon is located in exon 8. In addition to this major transcript, multiple other minor transcript variants have been described.

## Exon numbering

The *IKZF1* exon numbering used in this P335-C2 ALL-*IKZF1* product description is the exon numbering from the MANE Select transcript sequence (<https://www.ncbi.nlm.nih.gov/refseq/MANE/>). From product description version C2-02 onwards, we have adopted the MANE exon numbering. Please be aware that the MANE and LRG exon numbering do not always correspond, and MANE exon numbering used here may differ from literature. As changes to the databases can occur after the release of this product description, the NM\_ sequence and exon numbering may not be up-to-date. Exon numbering used here may differ from literature.

## Probemix content

The SALSA MLPA Probemix P335-C2 ALL-*IKZF1* contains 57 MLPA probes with amplification products between 120 and 504 nucleotides (nt). This includes one probe for each of the eight exons of the *IKZF1* gene (7p12.2). Furthermore, this probemix also contains seven probes for *PAX5* (9p13.2), six probes for *ETV6* (12p13.2), five probes for *RB1* (13q14.2), four probes for *BTG1* and the *BTG1* downstream region (12q21.33), four probes for *EBF1* (5q33.3), three probes for *CDKN2A/CDKN2B* (9p21.3) and five probes for the Xp22.33 region (PAR1 region; *SHOX* area, *CRLF2*, *CSF2RA*, *IL3RA* and *P2RY8* genes). In addition, one probe at Yp11.31 (*ZFY*) and one probe at 9p24.1 (*JAK2*) have been included to help determine the extent of a deletion/duplication detected in patient samples. See page 6 and 7 for more information about interpretation of results of the *ZFY* probe. Finally, 13 reference probes are included that target relatively copy number stable regions in ALL. Complete probe sequences and the identity of the genes detected by the reference probes are available in Table 3 and online ([www.mrcholland.com](http://www.mrcholland.com)).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at [www.mrcholland.com](http://www.mrcholland.com).

| Length (nt) | Name   |
|-------------|--|
| 64-70-76-82 | Q-fragments (only visible with <100 ng sample DNA)         |
| 88-96       | D-fragments (low signal indicates incomplete denaturation) |
| 92          | Benchmark fragment   |
| 100         | X-fragment (X chromosome specific)                         |
| 105         | Y-fragment (Y chromosome specific)                         |

### MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola (2012).

### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from male, healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

### Required specimens

Extracted DNA from peripheral whole blood or bone marrow, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol. The minimum percentage of tumour cells required for reliable analysis is 30% (Al Zaabi et al. 2010, Coll-Mulet et al. 2008). We advise using tumour samples with at least 50% tumour cell content. Therefore, tumour samples should be evaluated by a pathologist before the extraction of DNA.

### Reference samples

A sufficient number ( $\geq 3$ ) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. For example, when testing DNA samples derived from bone marrow, use reference samples from bone marrow in the same experiment for optimal data normalisation. Reference samples should be derived from healthy **male** individuals who are from families without a history of ALL. For more information on the need of male reference samples refer to section P335 specific note below. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

The following cell lines available from Coriell and DSMZ have been tested with P335-C2 probemix and can be used as positive DNA samples:

| Sample name          | Source            | Chromosomal position of CNA* | Altered target (and flanking) genes in P335-C2         | Expected copy number alteration                  |
|----------------------|-------------------|------------------------------|--|--|
| NA01353              | Coriell Institute | Xp22.33-PAR1 region          | <i>SHOX, CRLF2, CSF2RA</i> (part of PAR1 region)       | Heterozygous deletion                            |
| NA01750              | Coriell Institute | 9p21.3-p24.1                 | <i>JAK2, CDKN2A, CDKN2B</i>                            | Heterozygous duplication                         |
| NA04371              | Coriell Institute | 5q33.3                       | <i>EBF1</i>  | Heterozygous duplication                         |
|                      |                   | Xp22.33-PAR1 region          | <i>CSF2RA</i>  | Heterozygous duplication                         |
| NA05067              | Coriell Institute | 9p13.2-p24.1                 | <i>JAK2, CDKN2A, CDKN2B, PAX5</i> (all exons)          | Heterozygous duplication                         |
| NA07081              | Coriell Institute | 7p12.2                       | <i>IKZF1</i> (all exons)                               | Heterozygous duplication                         |
| NA07981              | Coriell Institute | 12p13.2                      | <i>ETV6</i> (all exons)                                | Heterozygous triplication/homozygous duplication |
| NA09403              | Coriell Institute | Xp22.33-PAR1 region          | <i>SHOX, CRLF2, CSF2RA, IL3RA, P2RY8</i> (PAR1 region) | Heterozygous deletion                            |
| NA10925              | Coriell Institute | 7p12.2                       | <i>IKZF1</i>   | Heterozygous deletion                            |
|                      |                   | Xp22.33-PAR1 region          | <i>CRLF2, CSF2RA</i> (part of PAR1 region)             | Heterozygous duplication                         |
| NA12606              | Coriell Institute | 13q14.2                      | <i>RB1</i>   | Heterozygous duplication                         |
| NA12722 <sup>◇</sup> | Coriell Institute | 9p21.3-p24.1                 | <i>JAK2, CDKN2A, CDKN2B</i>                            | Heterozygous duplication                         |
| NA14164              | Coriell Institute | 13q14.2                      | <i>RB1</i>   | Heterozygous deletion                            |
| BV-173 <sup>◇</sup>  | DSMZ              | 7p12.2                       | <i>IKZF1</i> (exons 1-7)                               | Heterozygous deletion                            |
|                      |                   | 9p21.3                       | <i>CDKN2A</i> (exon 4)                                 | Heterozygous deletion                            |
|                      |                   | 9p21.3                       | <i>CDKN2A</i> (exon 2), <i>CDKN2B</i> (exon 2)         | Homozygous deletion                              |
|                      |                   | Xp22.33-PAR1 region          | PAR1 region  | Heterozygous deletion                            |

\* Indicated chromosomal bands accommodate genes targeted by MLPA probes, however, the whole extent of copy number alteration (CNA) present in this cell line cannot be determined by this P335-C2 ALL-*IKZF1* probemix.

<sup>◇</sup> In this indicated cell line sample some of the reference probes are also affected by CNAs.

### Performance characteristics

*IKZF1* deletions can be identified in approximately 70% of children with *BCR-ABL1* positive (Philadelphia chromosome, Ph+) ALL (2% of all pediatric ALL cases), in 10-15% of Philadelphia chromosome negative (Ph-) pediatric ALL cases and in 40% of adult B-ALL cases.

Genomic deletions in *IKZF1* are either whole gene deletions (25-50% of all *IKZF1* deletions; Mullighan et al. 2009a; Palmi et al. 2013) or intragenic deletions, with the most frequent intragenic deletion (exons 4-7) comprising 30-55% of all deletions (Kastner et al. 2013, Iacobucci et al. 2009). The analytical sensitivity and specificity for the detection of deletions or duplications in the *IKZF1* gene is very high and can be considered >99%.

Analytical performance can be compromised by: SNVs or other polymorphisms in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.

### Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at [www.mrcholland.com](http://www.mrcholland.com). Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

### Interpretation of results

The expected results for *IKZF1* gene specific MLPA probes are allele copy number of 2 (normal), 1 (heterozygous deletion) and 0 (homozygous deletion). Duplication of the *IKZF1* gene is not expected in ALL patient samples, except in the context of larger duplications within or of chromosome arm 7p (observed in less than 5% of ALL cases according to the Progenetix.org database). Duplication of *IKZF1* is not expected to impact patient stratification.

For all other (pseudo)autosomal genes in the P335 ALL-*IKZF1* probemix, the expected results are allele copy number of 2 (normal), 1 (heterozygous deletion), 0 (homozygous deletion), 3 (heterozygous gain) or  $\geq 4$  (heterozygous triplication/homozygous gain). More information on expected deletions or duplications per gene can be found in Table 2.

The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$ . When this criterion is fulfilled, the following cut-off values for the Final ratio (FR) of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

| Copy number status                        | Final ratio (FR)   |
|---|--------------------|
| Normal                                    | $0.80 < FR < 1.20$ |
| Homozygous deletion                       | FR = 0             |
| Heterozygous deletion                     | $0.40 < FR < 0.65$ |
| Heterozygous gain                         | $1.30 < FR < 1.65$ |
| Heterozygous triplication/homozygous gain | $1.75 < FR < 2.15$ |
| Ambiguous                                 | All other values   |

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software (Calculations, cut-offs and interpretation remain unchanged). Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

**Please note that these above mentioned final ratios are affected both by percentage of tumour cells and by possible subclonality.**

In case of a heterozygous deletion that is subclonal and/or the sample contains a lower percentage of tumour cells, the FR might not be in the expected range of 0.40-0.65 (as indicated in the table above). For example, in case of ~50% tumour cell content or a copy number alteration present in ~50% of the tumour cells, the FR will be ~0.75. However, the same FR (0.75) will also be found for a sample with a tumour cell percentage of 25% or a subclone comprising 25% of all tumour cells that harbours a homozygous deletion. The MLPA technique cannot discriminate between these two scenarios.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in subclonal cases.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination in the DNA sample) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples
































with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.

- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

### P335 specific note:

In this probemix, the 120 nt Y-specific target probe (ZFY) is included as a flanking probe to help determine the extent of a deletion in the PAR1 region and to detect loss of the whole Y chromosome in male samples. To ensure that the comparative analysis can be completed for all samples, **only male reference samples must be used**. The comparative analysis will fail if female reference samples are used.

- In male samples, the 120 nt Y-specific target probe will be counted as an “additional” probe when using Coffalyser.Net software (see Figure 1). In healthy male individuals a final ratio of ~1 will be obtained for the Y-specific target probe after the comparative analysis. If the Y chromosome is lost in male samples, rather than deletion of the PAR1 region only, this can also be confirmed by visual examination of the 105 nt Y-fragment peak in the electropherogram (fragment analysis).
- In female samples, the 120 nt Y-specific target probe will be reported as “absent”, i.e. will have a final ratio of 0, after the comparative analysis.

| sample name        | sample type | bin smpl                 | FRSS  | FMRS  | probes | △ | DNA   | DD  | X   | Y   |
|--------------------|-------------|--------------------------|---|---|--------|---|---|---|---|---|
| Reference sample 1 | reference   | <input type="checkbox"/> |  |  | 57/57  | ✓ |  |  |  |  |
| Reference sample 2 | reference   | <input type="checkbox"/> |  |  | 57/57  | ✓ |  |  |  |  |
| Reference sample 3 | reference   | <input type="checkbox"/> |  |  | 57/57  | ✓ |  |  |  |  |
| Sample             | sample      | <input type="checkbox"/> |  |  | 56/56  | ✓ |  |  |  |  |
| noDNA              | no DNA      | <input type="checkbox"/> |  |   | 0/0    | ✓ |  |  |  |  |

**Figure 1.** Example fragment analysis overview in Coffalyser.Net with three male reference samples, a female test sample and a no DNA reaction. Note that the number of detected/expected probes differs between male (57/57) and female (56/56) samples.

### Limitations of the procedure

- In the majority of patients, defects in the *IKZF1* gene are deletions, but point mutations can occur which will not be detected using SALSA MLPA Probemix P335 ALL-*IKZF1*.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe

signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

- MLPA analysis on tumour samples provides information on the *average* situation in the cells from which the DNA sample was purified. Gains or losses of genomic regions or genes may not be detected if the percentage of tumour cells is low. In addition, subclonality of the aberration affects the final ratio of the corresponding probe. Furthermore, there is always a possibility that one or more reference probes *do* show a copy number alteration in a patient sample.

### Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

### IKZF1 mutation databases

<https://databases.lovd.nl/shared/genes/IKZF1>, <http://www.hgmd.cf.ac.uk/ac/gene.php?gene=IKZF1>, and <http://cancer.sanger.ac.uk/cosmic/>.

We strongly encourage users to deposit positive results in one of these databases. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report false positive results due to SNVs and unusual results (e.g., a deletion of *IKZF1* exons 6 and 8 but not exon 7) to MRC Holland: [info@mrcholland.com](mailto:info@mrcholland.com).

**Table 1. SALSA MLPA Probemix P335-C2 ALL-IKZF1**

| Length (nt) | SALSA MLPA probe   | Chromosomal position (hg18) |                     | Location (hg18) in kb |
|-------------|--|-----------------------------|---------------------|-----------------------|
|             |  | Reference                   | Target region       |                       |
| 64-105      | Control fragments – see table in probemix content section for more information |                             |                     |                       |
| 120 ~       | <b>ZFY probe</b> S0135-L27684  |                             | Yp11.31             | Y-002.889             |
| 124         | Reference probe 15370-L13762   | 7q11                        |                     | 07-075.448            |
| 130         | Reference probe 13867-L15385   | 16p13                       |                     | 16-008.765            |
| 136         | <b>CRLF2 probe</b> 13889-L15427  |                             | Xp22.33-PAR1 region | X-001.281             |
| 142         | <b>IKZF1 probe</b> 21511-L30001  |                             | 7p12.2              | 07-050.418            |
| 149 ¥       | <b>SHOX-area probe</b> 05648-L32510  |                             | Xp22.33-PAR1 region | X-000.771             |
| 157 ¥ ~     | <b>JAK2 probe</b> 07452-L32509   |                             | 9p24.1              | 09-005.113            |
| 161 ¥       | <b>PAX5 probe</b> 12501-L32513   |                             | 9p13.2              | 09-037.024            |
| 166         | <b>CSF2RA probe</b> 13892-L16221   |                             | Xp22.33-PAR1 region | X-001.374             |
| 172         | <b>PAX5 probe</b> 14647-L15394   |                             | 9p13.2              | 09-037.011            |
| 178 ~       | <b>BTG1-area probe</b> 18021-L22630  |                             | 12q21.33            | 12-090.658            |
| 184         | <b>P2RY8 probe</b> 17837-L15740  |                             | Xp22.33-PAR1 region | X-001.545             |
| 192 ¥       | Reference probe 06941-L32511   | 11q12                       |                     | 11-061.484            |
| 196         | <b>ETV6 probe</b> 17838-L22035   |                             | 12p13.2             | 12-011.936            |
| 202         | <b>PAX5 probe</b> 17839-L22036   |                             | 9p13.2              | 09-036.957            |
| 208 «       | <b>IKZF1 probe</b> 14056-L15654  |                             | 7p12.2              | 07-050.329            |
| 214 ¥       | Reference probe 13265-L32718   | 1p21                        |                     | 01-103.204            |
| 220         | <b>RB1 probe</b> 01782-L01346  |                             | 13q14.2             | 13-047.821            |
| 226         | <b>EBF1 probe</b> 12509-L13559   |                             | 5q33.3              | 05-158.459            |
| 232         | <b>BTG1 probe</b> 21378-L30126   |                             | 12q21.33            | 12-091.063            |
| 239         | <b>CDKN2B probe</b> 16059-L30167   |                             | 9p21.3              | 09-021.996            |
| 244         | <b>ETV6 probe</b> 13874-L17160   |                             | 12p13.2             | 12-011.883            |
| 252         | <b>CDKN2A probe</b> 10333-L30127   |                             | 9p21.3              | 09-021.965            |
| 258         | Reference probe 04534-L22019   | 2q24                        |                     | 02-166.606            |
| 264         | <b>IKZF1 probe</b> 13873-L15917  |                             | 7p12.2              | 07-050.412            |
| 269 «       | <b>IKZF1 probe</b> 13877-L15918  |                             | 7p12.2              | 07-050.315            |
| 274         | <b>PAX5 probe</b> 17840-L22037   |                             | 9p13.2              | 09-036.872            |
| 282         | <b>PAX5 probe</b> 13870-L15920   |                             | 9p13.2              | 09-036.993            |
| 288         | <b>IKZF1 probe</b> 17109-L20256  |                             | 7p12.2              | 07-050.436            |
| 295 +       | Reference probe 10435-L22110   | 9q34                        |                     | 09-136.850            |
| 301 «       | <b>ETV6 probe</b> 14058-L15656   |                             | 12p13.2             | 12-011.694            |
| 309         | <b>CDKN2A probe</b> 17814-L22631   |                             | 9p21.3              | 09-021.958            |
| 315         | <b>RB1 probe</b> 01789-L22025  |                             | 13q14.2             | 13-047.851            |
| 324         | Reference probe 03918-L20270   | 15q21                       |                     | 15-046.585            |
| 330         | <b>BTG1 probe</b> 12553-L22632   |                             | 12q21.33            | 12-091.062            |
| 337         | <b>PAX5 probe</b> 17841-L22038   |                             | 9p13.2              | 09-036.913            |
| 344 ¥       | <b>IKZF1 probe</b> 13869-L32515  |                             | 7p12.2              | 07-050.427            |
| 351         | <b>IL3RA probe</b> 13907-L22294  |                             | Xp22.33-PAR1 region | X-001.416             |
| 358         | <b>RB1 probe</b> 01792-L22295  |                             | 13q14.2             | 13-047.928            |
| 364         | Reference probe 14675-L16327   | 3p25                        |                     | 03-010.142            |
| 372         | <b>EBF1 probe</b> 14059-L30509   |                             | 5q33.3              | 05-158.137            |
| 379         | <b>IKZF1 probe</b> 15427-L22113  |                             | 7p12.2              | 07-050.338            |
| 389 ¥ «     | <b>ETV6 probe</b> 14060-L32514   |                             | 12p13.2             | 12-011.694            |
| 394         | <b>PAX5 probe</b> 17842-L22633   |                             | 9p13.2              | 09-036.830            |
| 401         | <b>ETV6 probe</b> 13875-L22014   |                             | 12p13.2             | 12-011.797            |
| 409 ~       | <b>BTG1-area probe</b> 18022-L22363  |                             | 12q21.33            | 12-091.006            |
| 418         | Reference probe 10063-L30260   | 8q22                        |                     | 08-100.251            |
| 427         | <b>RB1 probe</b> 01797-L16909  |                             | 13q14.2             | 13-047.945            |
| 436         | <b>EBF1 probe</b> 13868-L22053   |                             | 5q33.3              | 05-158.072            |
| 445         | <b>RB1 probe</b> 01799-L01362  |                             | 13q14.2             | 13-047.949            |
| 454         | Reference probe 18691-L02476   | 5p15                        |                     | 05-009.094            |



| Length (nt) | SALSA MLPA probe                | Chromosomal position (hg18) |               | Location (hg18) in kb |
|-------------|---------------------------------|-----------------------------|---------------|-----------------------|
|             |                                 | Reference                   | Target region |                       |
| 463 ¥       | <b>EBF1 probe</b> 23047-L32673  |                             | 5q33.3        | 05-158.058            |
| 470         | <b>IKZF1 probe</b> 14061-L22112 |                             | 7p12.2        | 07-050.422            |
| 478         | Reference probe 14846-L22111    | 3q11                        |               | 03-099.783            |
| 485         | <b>ETV6 probe</b> 13871-L22009  |                             | 12p13.2       | 12-011.914            |
| 494         | Reference probe 15203-L16978    | 3p12                        |               | 03-081.775            |
| 504         | Reference probe 09870-L19465    | 2p15                        |               | 02-061.126            |

¥ Changed in version C2. Minor alteration, no change in sequence detected.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

± SNP rs556167410 could influence the 295 nt probe signal.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

**Table 2. P335-C2 probes arranged according to chromosomal location**

| Length (nt)  | SALSA MLPA probe | Gene / Exon <sup>a</sup> | Location (hg18)/ Ligation site | Partial sequence <sup>b</sup> (24 nt adjacent to ligation site) | Distance to next probe                |
|--|------------------|--------------------------|--------------------------------|---|---------------------------------------|
| <b>EBF1</b> gene, 5q33.3. Indicated locations are according to NM_001290360.3 and exon numbering according to MANE Select transcript. <i>EBF1</i> deletions are suggested to be an important factor in the relapse of ALL as 25% of relapsed ALL cases have deletions in the <i>EBF1</i> locus (Yang et al. 2008). Deletions of <i>EBF1</i> are suggested to be associated with a higher risk of relapse (Olsson et al. 2014).   |                  |                          |                                |   |                                       |
| 463  | 23047-L32673     | <b>EBF1</b> , ex 16      | 3043-3044                      | ATGGCCTCTCGA-GAGCTTGGTGGC                                       | 14.0 kb                               |
| 436  | 13868-L22053     | <b>EBF1</b> , ex 14      | 1888-1889                      | CCAATAAAGT-CCGTACCACGA  | 65.2 kb                               |
| 372  | 14059-L30509     | <b>EBF1</b> , ex 10      | 1416-1417                      | GTTGTGGAAGTC-ACACTGTCCTAC                                       | 322.2 kb                              |
| 226  | 12509-L13559     | <b>EBF1</b> , ex 1       | 278-279                        | ATTTGCTTTCCA-GCCCGCCTTGAT                                       | -                                     |
| <b>IKZF1</b> gene, 7p12.2. Indicated locations are according to NM_006060.6 and exon numbering according to MANE Select transcript. Exon numbering is different from the exon numbering used by Mullighan et al. (2009a), where exon 3-6 deletions correspond to exon 4-7 in the MANE exon numbering. Deletions or mutations of <i>IKZF1</i> are significantly associated with an increased risk of relapse and adverse events in ALL patients (Mullighan et al. 2009a; Martinelli et al. 2009; Iacobucci et al. 2009; Stanulla et al. 2020). In the absence of <i>ERG</i> deletions and in combination with <i>CDKN2A</i> , <i>CDKN2B</i> , <i>PAX5</i> and <i>PAR1</i> deletions, deletions of <i>IKZF1</i> ( <i>IKZF1</i> <sup>plus</sup> profile) are associated with a very poor prognosis (Stanulla et al. 2018). Additional <i>IKZF1</i> probes as well as probes for the <i>ERG</i> gene are present in the P202 <i>IKZF1-ERG</i> probemix.  |                  |                          |                                |   |                                       |
| 269 «  | 13877-L15918     | <b>IKZF1</b> , ex 1      | 187-188                        | TCTTGCCCCAA-AGCGCGACGCAC  | 14.2 kb                               |
| 208 «  | 14056-L15654     | <b>IKZF1</b> , ex 2      | 253-254                        | AGACATGTCCCA-AGTTTCAGGTGA                                       | 8.6 kb                                |
| 379  | 15427-L22113     | <b>IKZF1</b> , ex 3      | 355-356                        | GGGAGGACAGCA-AAGCTCCAAGAG                                       | 74.0 kb                               |
| 264  | 13873-L15917     | <b>IKZF1</b> , ex 4      | 477-478                        | TACGAATGCTTG-ATGCCTCGGGAG                                       | 5.9 kb                                |
| 142  | 21511-L30001     | <b>IKZF1</b> , ex 5      | 680-681                        | TGCGGGCCTCA-TTCACCCAGAAG  | 4.8 kb                                |
| 470  | 14061-L22112     | <b>IKZF1</b> , ex 6      | 812-813                        | TTTTCTGCAGTT-GGTAAACCTCAC                                       | 4.4 kb                                |
| 344  | 13869-L32515     | <b>IKZF1</b> , ex 7      | 994-995                        | CAAGATAGGATC-AGAGAGATCTCT                                       | 9.3 kb                                |
| 288  | 17109-L20256     | <b>IKZF1</b> , ex 8      | 2271-2272                      | GGTGTGCCGCCA-CCCAAGTGCCAA                                       | 25 Mb to <i>POR</i> (reference probe) |
| <b>CDKN2A/CDKN2B</b> genes, 9p21.3. Exon numbering of <i>CDKN2A</i> is based on the MANE Select transcript (NM_000077.5, p16 <sup>INK4a</sup> ) and the MANE Plus Clinical transcript (NM_058195.4, p14 <sup>ARF</sup> ). The <i>CDKN2A</i> exon numbering has been changed; the exon numbering (LRG) used in previous versions of this product description can be found in between brackets. Deletions of <i>CDKN2A</i> are present in 50-81% of T-ALL, in 15-35% of childhood B-ALL and in 30-45% of adult B-ALL cases (González-Gil et al. 2021). <i>CDKN2A/2B</i> deletions are more frequently found in high-risk patients of all ages, and the majority of studies on relapse ALL suggests that (homozygous) <i>CDKN2A/2B</i> deletions are more frequent at relapse than at diagnosis. In combination with <i>IKZF1</i> deletions and in the absence of <i>ERG</i> deletions ( <i>IKZF1</i> <sup>plus</sup> profile), deletions of the <i>CDKN2A</i> and/or <i>CDKN2B</i> genes are associated with a very poor prognosis (Stanulla et al. 2018). Additional <i>CDKN2A/2B</i> probes are present in the P419 <i>CDKN2A/2B-CDK4</i> and ME024 9p21 probemixes. The latter probemix detects both copy number and methylation changes of the <i>CDKN2A/2B</i> genes. |                  |                          |                                |   |                                       |

| Length (nt)   | SALSA MLPA probe | Gene / Exon <sup>a</sup> | Location (hg18)/ Ligation site                                       | Partial sequence <sup>b</sup> (24 nt adjacent to ligation site) | Distance to next probe             |
|---|------------------|--------------------------|--|---|------------------------------------|
| 157 ~   | 07452-L32509     | <b>JAK2</b>              | 9p24.1   | GAATCACTGACA-GAGAGCAAGTTT                                       | 17 Mb                              |
| 309   | 17814-L22631     | <b>CDKN2A</b> , ex 3 (4) | NM_000077.5 830-831; NM_058195.4; 904-905                            | TTGCGAGCCTCG-CAGCCTCCGGAA                                       | 7.1 kb                             |
| 252   | 10333- L30127    | <b>CDKN2A</b> , ex 1 (2) | NM_000077.5; 138 nt before exon 1; NM_058195.4; 3.8 kb before exon 2 | GCCTGGAAAGAT-ACCGCGGTCCCT                                       | 30.9 kb                            |
| 239   | 16059- L30167    | <b>CDKN2B</b> , ex 2     | 9p21.3   | GCCTGTCTGAGA-CTCACAGGAAGG                                       | 15 Mb to PAX5                      |
| <p><b>PAX5</b> gene, 9p13.2. Indicated locations are according to NM_016734.3 and exon numbering according to MANE Select transcript. PAX5 deletions in T-ALL are large, involving other genes and sometimes even extending into the CDKN2A/2B genes (Schwab and Harrison 2018). PAX5 deletions are present in 30% of B-progenitor ALL cases and in 50% of Ph+ ALL (Lejman et al. 2022; Li et al. 2021). Moreover, intragenic PAX5 amplifications were suggested to be present in a new subgroup of B-cell precursor ALL (~1%) and in 3% of B-other ALL cases, which is associated with poor outcome (Schwab et al. 2017). In combination with IKZF1 deletions and in the absence of ERG deletions (IKZF1<sup>plus</sup> profile), deletions of the PAX5 gene are associated with a very poor prognosis (Stanulla et al. 2018).</p> |                  |                          |  |   |                                    |
| 394   | 17842-L22633     | <b>PAX5</b> , ex 10      | 2018-2019  | CTCCTTCTTTAG-TATCTTTACGAG                                       | 42.0 kb                            |
| 274   | 17840-L22037     | <b>PAX5</b> , ex 8       | 2 nt after exon 8  | ATGGTGCCTGGT-GAGTTTGCCTG  | 41.4 kb                            |
| 337   | 17841-L22038     | <b>PAX5</b> , ex 7       | 1108-1109  | CTGACATCGGGA-GCAGTGTGCCAG                                       | 43.2 kb                            |
| 202   | 17839-L22036     | <b>PAX5</b> , ex 6       | 968-969  | GTTTGAGAGGCA-GCACTACTCAGA                                       | 36.1 kb                            |
| 282   | 13870-L15920     | <b>PAX5</b> , ex 5       | 756-757  | GTGAGCACGGAT-TCGGCCGGCTCG                                       | 18.0 kb                            |
| 172   | 14647-L15394     | <b>PAX5</b> , ex 2       | 378-379  | CTTGCTCATCAA-GGTGTCAGGCC  | 13.6 kb                            |
| 161   | 12501-L32513     | <b>PAX5</b> , ex 1       | 43 nt before exon 1  | CATCTTGATG-TTGGCGAGAACA   | 100 Mb to COL5A1 (reference probe) |
| <p><b>ETV6</b> gene, 12p13.2. Indicated locations are according to NM_001987.5 and exon numbering according to MANE Select transcript. ETV6 deletions are more frequent in B-ALL (51%) as compared to T-ALL (4%), and ETV6 is often involved in rearrangements and fusions detected in ALL patients (e.g. in ETV6-RUNX1 fusions) (Schwab et al. 2013). It is suggested that native ETV6 deletions in ETV6-RUNX1+ childhood ALL is associated with better prognosis (Ko et al. 2011). Although microdeletions often occur at the translocation breakpoints, this MLPA probemix will not detect all microdeletions in which ETV6 is involved.</p>   |                  |                          |  |   |                                    |
| 301 «   | 14058-L15656     | <b>ETV6</b> , ex 1       | 366-367  | AATGACCGGTC-TGGCTGGCCGTG  | 0.1 kb                             |
| 389 «   | 14060-L32514     | <b>ETV6</b> , ex 1       | 480-481  | TGCTCAGTGTAG-CATTAAGGTA   | 102.4 kb                           |
| 401   | 13875-L22014     | <b>ETV6</b> , ex 2       | 563-564  | TTCATGTTCCAG-TGCCTCGAGCGC                                       | 86.6 kb                            |
| 244   | 13874-L17160     | <b>ETV6</b> , ex 3       | 641-642  | TTTACTGGAGCA-GGGATGACGTAG                                       | 30.3 kb                            |
| 485   | 13871-L22009     | <b>ETV6</b> , ex 5       | 967-968  | AATGTGCACCAT-AACCCTCCCACC                                       | 22.3 kb                            |
| 196   | 17838-L22035     | <b>ETV6</b> , ex 8       | 2530-2531  | AGTCTTGGGGAT-TGTTGGCACCTA                                       | 79 Mb to BTG1                      |
| <p><b>BTG1</b> gene, 12q21.33. Indicated locations are according to NM_001731.3 and exon numbering according to MANE Select transcript. BTG1 deletions appeared to be more frequent in high risk ALL cases and are often combined with other deletions (Fang et al. 2018). BTG1 deletions are detected with high frequency in ALL patients with Down syndrome (Lundin et al. 2012). As many deletions of BTG1 extend into the centromeric area of BTG1 into a gene poor region (Waanders et al. 2012), two downstream probes were included for the BTG1-area.</p>   |                  |                          |  |   |                                    |
| 178 ~   | 18021-L22630     | <b>BTG1-area</b>         | Downstream of BTG1   | CACTAAAAATGT-GCATACTTCAAC                                       | 348.3 kb                           |
| 409 ~   | 18022-L22636     | <b>BTG1-area</b>         | Downstream of BTG1   | TGGAAAATGGGA-ATGTTCCAGGGT                                       | 55.7 kb                            |
| 330 #   | 12553-L22632     | <b>BTG1</b> , ex 2       | 1074-1075  | TTGCTAGGGAGG-GAAGTCCTAGGG                                       | 1.6 kb                             |
| 232   | 21378-L30126     | <b>BTG1</b> , ex 1       | 394-395  | GTTTCTCCGCAC-CAAGGGGCTCAC                                       | -                                  |
| <p><b>RB1</b> gene, 13q14.2. Indicated locations are according to NM_000321.3 and exon numbering according to MANE Select transcript. RB1 deletions have been reported to be more frequent in high risk ALL as compared to non-selected cases and often co-occur with other CNAs, in particular iAMP21 (Fang et al. 2018; Schwab et al. 2013; Steeghs et al. 2019). Deletions often involve one or more of the last 10 exons of this 27-exon gene (Schwab et al. 2013). Additional RB1 probes are available in the P047 RB1 probemix.</p>   |                  |                          |  |   |                                    |
| 220   | 01782-L01346     | <b>RB1</b> , ex 6        | 83 nt after exon 6   | ATTCCCAATTT-TTATTGAGTAAT  | 30.2 kb                            |

| Length (nt)  | SALSA MLPA probe | Gene / Exon <sup>a</sup> | Location (hg18)/ Ligation site | Partial sequence <sup>b</sup> (24 nt adjacent to ligation site) | Distance to next probe |
|--|------------------|--------------------------|--------------------------------|---|------------------------|
| 315  | 01789-L22025     | <b>RB1</b> , ex 14       | 268 nt before exon 14          | GCTTTTGTGTTG-TCTTGGCGGCCA                                       | 76.9 kb                |
| 358  | 01792-L22295     | <b>RB1</b> , ex 19       | 2041-2042                      | ATTCTACTGCAA-ATGCAGAGACAC                                       | 16.8 kb                |
| 427  | 01797-L16909     | <b>RB1</b> , ex 24       | 265 nt before exon 24          | GAAACTGCCTT-TGCCCTCCCTAA  | 3.8 kb                 |
| 445  | 01799-L01362     | <b>RB1</b> , ex 26       | 2852-2853                      | AGAGTCCAAATT-TCAGCAGAAACT                                       | –                      |
| <p><b>Xp22.33 / Yp11.32 (PAR1) region.</b> With exception of the ZFY probe, all these probes target the PAR1 region, which is present in two copies per cell, irrespective of gender. This region shows recurrent aberrations in ALL (Harvey et al. 2010). <i>CSF2RA/IL3RA</i> deletions are frequent in ALL (7% of ALL cases; up to 55% in Down syndrome-associated ALL). Focal deletions within the <i>P2RY8</i> and <i>CRLF2</i> genes, resulting in a fusion gene, are associated with a poor prognosis (Mullighan et al. 2009b). In combination with <i>IKZF1</i> deletions and in the absence of <i>ERG</i> deletions (<i>IKZF1</i><sup>plus</sup> profile), deletions in the PAR1 region (deletion of <i>CSF2RA</i> and <i>IL3RA</i> and retention of the <i>CRLF2</i> probe associated with <i>P2RY8-CRLF2</i> fusion) are associated with a very poor prognosis (Stanulla et al. 2018). Please note that loss in signal may also be due to the loss of the Y chromosome in male samples, which can be detected by comparing the PAR1 results to the ZFY probe at 120 nt and to the Y-specific control fragment at 105 nt. The <i>SHOX</i>-area probe is located in between the <i>SHOX</i> and <i>CRLF2</i> genes, to enable determination of the extent of a deletions within the PAR1 region. More probes for the PAR1 region are included in the P329 <i>CRLF2-CSF2RA-IL3RA</i> and P018 <i>SHOX</i> probemixes.</p> |                  |                          |                                |   |                        |
| 149  | 05648-L32510     | <b>SHOX-area</b>         | Xp22.33                        | TGGTGCTGAAAT-GAGGAAGCCCTG                                       | 510.7 kb               |
| 136  | 13889-L15427     | <b>CRLF2</b>             | Xp22.33                        | GGATCTCCTCTA-TGAGGTTTCAGTA                                      | 93.0 kb                |
| 166  | 13892-L16221     | <b>CSF2RA</b>            | Xp22.33                        | TTTCACTTACCA-GTAGGTTTTCCG                                       | 41.4 kb                |
| 351  | 13907-L22294     | <b>IL3RA</b>             | Xp22.33                        | GGAAGATATCAG-AAACATCCTAGG                                       | 129.4 kb               |
| 184  | 17837-L15740     | <b>P2RY8</b>             | Xp22.33                        | TTTACGCAAACA-TGTATTCCAGCA                                       | 1.3 Mb                 |
| 120 ~  | S0135-L27684     | <b>ZFY</b>               | Yp11.31                        | TCATAGAGGAGG-ATGTTTCAGTGCT                                      | –                      |

<sup>a</sup> See section Exon numbering on page 2 for more information.

<sup>b</sup> Only partial probe sequences are shown. Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com). Please notify us of any mistakes: [info@mrcholland.com](mailto:info@mrcholland.com).

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

– Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

# This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

**Table 3. Reference probes arranged according to chromosomal location.**

| Length (nt) | SALSA MLPA probe | Gene           | Chromosomal band (hg18) | Partial sequence (24 nt adjacent to ligation site) | Location (hg18) in kb |
|-------------|------------------|----------------|-------------------------|--|-----------------------|
| 214         | 13265-L32718     | <i>COL11A1</i> | 1p21                    | CCTCCTTAGGGA-TTTC AAGGCAAG                         | 01-103.204            |
| 504         | 09870-L19465     | <i>PEX13</i>   | 2p15                    | TGAGGATGACCA-TGTAGTTGCCAG                          | 02-061.126            |
| 258         | 04534-L22019     | <i>SCN1A</i>   | 2q24                    | ACTGTTCTCCAT-ATTGGTTAAAAG                          | 02-166.606            |
| 364         | 14675-L16327     | <i>BRK1</i>    | 3p25                    | GCAACACTAAAC-GAGAAATTGACA                          | 03-010.142            |
| 494         | 15203-L16978     | <i>GBE1</i>    | 3p12                    | GACCTAGAGGGA-CTCATGATCTTT                          | 03-081.775            |
| 478         | 14846-L22111     | <i>CPOX</i>    | 3q11                    | CAGAATTGAAAG-TATCTTGATGTC                          | 03-099.783            |
| 454         | 18691-L02476     | <i>SEMA5A</i>  | 5p15                    | GTCCATCACTGT-GTAGCTACCGTT                          | 05-009.094            |
| 124         | 15370-L13762     | <i>POR</i>     | 7q11                    | GATGGGAAGTGA-GTGCCACCCTG                           | 07-075.448            |
| 418         | 10063-L30260     | <i>VPS13B</i>  | 8q22                    | TCTTTATGGGAA-ACTTCTGAAACT                          | 08-100.252            |
| 295 ^ +     | 10435-L22110     | <i>COL5A1</i>  | 9q34                    | AGGGCCTTCCAA-GCCGGCTTCTCC                          | 09-136.850            |
| 192         | 06941-L32511     | <i>BEST1</i>   | 11q12                   | GCACCAGGACCT-GCCTCGGATGGA                          | 11-061.484            |
| 324         | 03918-L20270     | <i>FBN1</i>    | 15q21                   | CCTACAGATGTG-AATGCTTCCCTG                          | 15-046.585            |
| 130         | 13867-L15385     | <i>ABAT</i>    | 16p13                   | ACTTTGTGGAGA-AGCTCCGGCAGT                          | 16-008.765            |

<sup>^</sup> In comparison to focal chromosome 9p deletions, including the *JAK2*, *CDKN2A/2B* and *PAX5* genes, the *COL5A1* probe at 9q34 is also affected in case of complete chromosome 9 deletions.

± SNP rs556167410 could influence the 295 nt probe signal.

Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com).

## Related SALSA MLPA probemixes

|   |   |
|---|---|
| <b>P018 SHOX</b>                                      | Contains probes for the <i>SHOX</i> gene and Xp22 regions.  |
| <b>P047 RB1</b>                                       | Contains probes for 26 out of 27 exons of the <i>RB1</i> gene.  |
| <b>P202 IKZF1-ERG *</b>                               | Contains two probes for all exons and the regulatory regions of <i>IKZF1</i> transcript variant 1 (NM_006060.6), one probe for each exon of the <i>ERG</i> gene, and probes for <i>CDKN2A/2B</i> genes and the 14q32.33 region.   |
| <b>P327 iAMP21-ERG</b>                                | Contains probes for <i>RUNX1</i> , <i>ERG</i> genes and iAMP21 detection in ALL.  |
| <b>P329 CRLF2-CSF2RA-IL3RA</b>                        | Contains probes for <i>CRLF2</i> , <i>CSF2RA</i> , <i>IL3RA</i> and <i>SHOX</i> genes, involved in B-ALL.   |
| <b>P377 Hematologic malignancies</b>                  | Contains probes for the most common copy number alterations in ALL, AML, CLL, CML, MDS and various lymphomas.   |
| <b>P383 T-ALL</b>                                     | Contains probes for <i>STIL-TAL1</i> , <i>LEF1</i> , <i>CASP8AP2</i> , <i>MYB</i> , <i>EZH2</i> , <i>MLLT3</i> , <i>MTAP</i> , <i>CDKN2A/2B</i> , <i>NUP214-ABL1</i> , <i>PTEN</i> , <i>LMO1/2</i> , <i>NF1</i> , <i>SUZ12</i> , <i>PTPN2</i> and <i>PHF6</i> genes, involved in T-ALL. |
| <b>P419 CDKN2A/2B-CDK4</b>                            | Contains probes for <i>CDKN2A</i> , <i>CDKN2B</i> and <i>CDK4</i> genes.  |
| <b>ME024 9p21 CDKN2A/2B region</b>                    | Contains probes for the 9p21 region, including <i>CDKN2A</i> and <i>CDKN2B</i> genes for detection of both copy number and methylation status.  |
| <b>D007 Acute Lymphoblastic Leukemia <sup>^</sup></b> | digitalMLPA probemix that contains probes for 52 target genes and three chromosomal regions associated with ALL, including 16 probes for <i>IKZF1</i> .   |

\* Please note that probemix P202 can only be used to obtain additional information, but cannot be used to confirm results of single probe aberrations detected in the *IKZF1* gene.

<sup>^</sup> For probes that have different ligation sites than the P335 probes this probemix can be used for confirmation. This concerns the majority but not all of the *IKZF1* probes in D007. For more information, please contact [info@mrcholland.com](mailto:info@mrcholland.com).

## References

- Al Zaabi EA et al. (2010). Multiplex ligation-dependent probe amplification versus multiprobe fluorescence in situ hybridization to detect genomic aberrations in chronic lymphocytic leukemia: a tertiary center experience. *J Mol Diag*. 12:197-203.
- Bernt KM and Hunger SP (2014). Current concepts in pediatric Philadelphia chromosome-positive acute lymphoblastic leukemia. *Front Oncol*. 4:54.
- Coll-Mulet L et al. (2008). Multiplex ligation-dependent probe amplification for detection of genomic alterations in chronic lymphocytic leukaemia. *Br J Haematol*. 142:793-801.
- Ezzat S et al. (2003). Ikaros isoforms in human pituitary tumors: distinct localization, histone acetylation, and activation of the 5' fibroblast growth factor receptor-4 promoter. *Am J Pathol*. 163:1177-84.
- Fang Q et al. (2018). Prognostic significance of copy number alterations detected by multi-link probe amplification of multiple genes in adult acute lymphoblastic leukemia. *Oncol Lett*. 15:5359-67.
- González-Gil C et al. (2021). The Yin and Yang-like clinical implications of the *CDKN2A/ARF/CDKN2B* gene cluster in acute lymphoblastic leukemia. *Genes (Basel)*. 12.
- Harvey RC et al. (2010). Rearrangement of *CRLF2* is associated with mutation of *JAK* kinases, alteration of *IKZF1*, Hispanic/Latino ethnicity, and a poor outcome in pediatric B-progenitor acute lymphoblastic leukemia. *Blood*, 115:5312-21.
- Hömig-Hölzel C and Savola S (2012). Multiplex ligation-dependent probe amplification (MLPA) in tumor diagnostics and prognostics. *Diagn Mol Pathol*. 21:189-206.
- Iacobucci I et al. (2009). Identification and molecular characterization of recurrent genomic deletions on 7p12 in the *IKZF1* gene in a large cohort of BCR-ABL1 –positive acute lymphoblastic leukemia patients: on

- behalf of Gruppo Italiano Malattie Ematologiche dell'Adulto Acute Leukemia Working Party (GIMEMAAL WP). *Blood*. 114:2159-67.
- Kastner P et al. (2013). Function of Ikaros as a tumor suppressor in B cell acute lymphoblastic leukemia. *Am J Blood Res*. 3:1-13.
  - Ko DH et al. (2011). Native ETV6 deletions accompanied by ETV6-RUNX1 rearrangements are associated with a favourable prognosis in childhood acute lymphoblastic leukaemia: a candidate for prognostic marker. *Br J Haematol*. 155:530-3.
  - Kuiper RP et al. (2010). IKZF1 deletions predict relapse in uniformly treated pediatric precursor B-ALL. *Leukemia*. 24:1258-64.
  - Lejman M et al. (2022). Genetic biomarkers and their clinical implications in B-cell acute lymphoblastic leukemia in children. *Int J Mol Sci*. 23:2755.
  - Li J et al. (2021a). Emerging molecular subtypes and therapeutic targets in B-cell precursor acute lymphoblastic leukemia. *Front Med*. 15:347-71.
  - Lundin C et al. (2012). High frequency of BTG1 deletions in acute lymphoblastic leukemia in children with down syndrome. *Genes Chromosomes Cancer*. 51:196-206.
  - Malard F et al. (2020). Acute lymphoblastic leukaemia. *Lancet*. 395:1146-62.
  - Martinelli G et al. (2009). IKZF1 (Ikaros) deletions in BCR-ABL1-positive acute lymphoblastic leukemia are associated with short disease-free survival and high rate of cumulative incidence of relapse: a GIMEMA AL WP report. *J Clin Oncol*. 27:5202-7.
  - Mullighan CG et al. (2008). BCR-ABL1 lymphoblastic leukaemia is characterized by the deletion of Ikaros. *Nature*. 453:110-4.
  - Mullighan CG et al. (2009a). Deletion of IKZF1 and prognosis in acute lymphoblastic leukemia. *N Engl J Med*. 360: 470-80.
  - Mullighan CG et al. (2009b). Rearrangement of CRLF2 in B-progenitor- and Down syndrome-associated acute lymphoblastic leukemia. *Nat Genet*. 41:1243-6.
  - Olsson L et al. (2014). Deletions of IKZF1 and SPRED1 are associated with poor prognosis in a population-based series of pediatric B-cell precursor acute lymphoblastic leukemia diagnosed between 1992 and 2011. *Leukemia*. 28:302-10.
  - Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res*. 30:e57.
  - Schwab C et al. (2013). Genes commonly deleted in childhood B-cell precursor acute lymphoblastic leukemia: association with cytogenetics and clinical features. *Haematologica*. 98:1081-8.
  - Schwab C et al. (2017). Intragenic amplification of PAX5: a novel subgroup in B-cell precursor acute lymphoblastic leukemia? *Blood Adv*. 1:1473-7.
  - Schwab C et al. (2018). Advances in B-cell precursor acute lymphoblastic leukemia genomics. *Hemasphere*. 2:e53.
  - Stanulla M et al. (2018). IKZF1plus defines a new minimal residual disease-dependent very-poor prognostic profile in pediatric B-cell precursor acute lymphoblastic leukemia. *J Clin Oncol*. 36:1240-9.
  - Stanulla M et al. (2020). IKZF1 deletions in pediatric acute lymphoblastic leukemia: still a poor prognostic marker? *Blood*. 135:252-60.
  - Steeghs EMP et al. (2019). Copy number alterations in B-cell development genes, drug resistance, and clinical outcome in pediatric B-cell precursor acute lymphoblastic leukemia. *Sci Rep*. 9:4634.
  - Van der Sligte NE et al. (2015). Effect of IKZF1 deletions on signal transduction pathways in Philadelphia chromosome negative pediatric B-cell precursor acute lymphoblastic leukemia (BCP-ALL). *Exp Hematol Oncol*. 4:23.
  - Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem*. 421:799-801.
  - Waanders E et al. (2012). The origin and nature of tightly clustered BTG1 deletions in precursor B-cell acute lymphoblastic leukemia support a model of multiclonal evolution. *PLoS Genet*. 8:21002533.
  - Yang JJ et al. (2008). Genome-wide copy number profiling reveals molecular evolution from diagnosis to relapse in childhood acute lymphoblastic leukemia. *Blood*. 112:4178-83.



## Selected publications using SALSA MLPA Probemix P335 ALL-IKZF1

- Agarwal M et al. (2018). Cyclin dependent kinase inhibitor 2A/B gene deletions are markers of poor prognosis in Indian children with acute lymphoblastic leukemia. *Pediatr Blood Cancer*. 65:e27001.
- Ayon-Perez MF et al. (2019). IKZF1 Gene Deletion in Pediatric Patients Diagnosed with Acute Lymphoblastic Leukemia in Mexico. *Cytogenet Genome Res*. 158:10-6.
- Bhandari P et al. (2017). Molecular profiling of gene copy number abnormalities in key regulatory genes in high-risk B-lineage acute lymphoblastic leukemia: frequency and their association with clinicopathological findings in Indian patients. *Med Oncol*. 34:92.
- Boer JM et al. (2016). Prognostic value of rare IKZF1 deletion in childhood B-cell precursor acute lymphoblastic leukemia: an international collaborative study. *Leukemia*. 30:32-8.
- Chiaretti S et al. (2021). Philadelphia-like acute lymphoblastic leukemia is associated with minimal residual disease persistence and poor outcome. First report of the minimal residual disease-oriented GIMEMA LAL1913. *Haematologica*. 106:1559-68.
- Chung-Filho AA et al. (2020). MTHFR rs1801133 polymorphism is associated with increased risk of B-cell precursor lymphoblastic leukaemia with recurrent genetic aberrations of fetal origin. *Cancer Epidemiol*. 65:101693.
- Clappier E et al. (2014). An intragenic ERG deletion is a marker of an oncogenic subtype of B-cell precursor acute lymphoblastic leukemia with a favorable outcome despite frequent IKZF1 deletions. *Leukemia*. 28:70-7.
- den Boer M et al. (2021). Outcomes of paediatric patients with B-cell acute lymphocytic leukaemia with ABL-class fusion in the pre-tyrosine-kinase inhibitor era: a multicentre, retrospective, cohort study. *Lancet Haematol*. 8:e55-e66.
- Dorge P et al. (2013). IKZF1 deletion is an independent predictor of outcome in pediatric acute lymphoblastic leukemia treated according to the ALL-BFM 2000 protocol. *Haematologica*. 98:428-32.
- Harrison C et al. (2014). An international study of intrachromosomal amplification of chromosome 21 (iAMP21): cytogenetic characterization and outcome. *Leukemia*. 28:1015-21.
- Konialis C et al. (2014). Routine application of a novel MLPA-based first-line screening test uncovers clinically relevant copy number aberrations in haematological malignancies undetectable by conventional cytogenetics. *Hematology*. 19:217-24.
- Krentz S et al. (2013). Prognostic value of genetic alterations in children with first bone marrow relapse of childhood B-cell precursor acute lymphoblastic leukemia. *Leukemia*. 27:295-304.
- Lopes BA et al. (2016). COBL is a novel hotspot for IKZF1 deletions in childhood acute lymphoblastic leukemia. *Oncotarget*. 7:53064-73.
- Marks DI et al. (2013). The clinical characteristics, therapy and outcome of 85 adults with acute lymphoblastic leukemia and t (4; 11)(q21; q23)/MLL-AFF1 prospectively treated in the UKALLXII/ECOG2993 trial. *Haematologica*. 98:945-52.
- Moorman AV (2016). New and emerging prognostic and predictive genetic biomarkers in B-cell precursor acute lymphoblastic leukemia. *Haematologica*. 101:407-16.
- Moorman AV et al. (2014). A novel integrated cytogenetic and genomic classification refines risk stratification in pediatric acute lymphoblastic leukemia. *Blood*. 124:1434-44.
- Moorman AV et al. (2010). Prognostic effect of chromosomal abnormalities in childhood B-cell precursor acute lymphoblastic leukaemia: results from the UK Medical Research Council ALL97/99 randomised trial. *Lancet Oncol*. 11:429-38.
- Olsson L et al. (2015). The clinical impact of IKZF1 deletions in paediatric B-cell precursor acute lymphoblastic leukaemia is independent of minimal residual disease stratification in Nordic Society for Paediatric Haematology and Oncology treatment protocols used between 1992 and 2013. *Br J Haematol*. 170:847-58.
- Palmi C et al. (2013). What is the relevance of Ikaros gene deletions as prognostic marker in pediatric Philadelphia negative B-cell precursor acute lymphoblastic leukemia? *Haematologica*. 98:1226-31.
- Ribera J et al. (2015). Prognostic significance of copy number alterations in adolescent and adult patients with precursor B acute lymphoblastic leukemia enrolled in PETHEMA protocols. *Cancer*. 121:3809-17.



- Ryan S et al. (2016). The role of the RAS pathway in iAMP21-ALL. *Leukemia*. 30:1824-31.
- Schwab CJ et al. (2013). Genes commonly deleted in childhood B-cell precursor acute lymphoblastic leukemia: association with cytogenetics and clinical features. *Haematologica*. 98:1081-8.
- Schwab CJ et al. (2010). Evaluation of multiplex ligation-dependent probe amplification as a method for the detection of copy number abnormalities in B-cell precursor acute lymphoblastic leukemia. *Genes Chromosomes Cancer*. 49:1104-13.
- Singh M et al. (2018). High frequency of intermediate and poor risk copy number abnormalities in pediatric cohort of B-ALL correlate with high MRD post induction. *Leuk Res*. 66:79-84.
- Steeghs EMP et al. (2019). Copy number alterations in B-cell development genes, drug resistance, and clinical outcome in pediatric B-cell precursor acute lymphoblastic leukemia. *Sci Rep*. 9:4634.
- Zaliova M et al. (2016). Characterization of leukemias with ETV6-ABL1 fusion. *Haematologica*. 101:1082-93.

| P335 product history |   |
|----------------------|---|
| Version              | Modification  |
| C2                   | Lengths of several probes are adjusted, but no changes in the sequences detected.   |
| C1                   | One target probe and two reference probes replaced, length of several probes changed.   |
| B2                   | Length of one probe changed.  |
| B1                   | Several target probes replaced, one new probe included for PAX5 exon 7 and two new flanking probes included for 9p and Yp chromosome arms. All reference probes replaced. |
| A4                   | The 88 and 96 nt DNA denaturation control fragments replaced.   |
| A3                   | Length of two probes changed.   |
| A2                   | Two reference probes replaced.  |
| A1                   | First release.  |

| Implemented changes in the product description   |
|--|
| <p>Version C2-02 – 03 March 2023 (04P)</p> <ul style="list-style-type: none"> <li>- Section Clinical background and background information in Table 2 updated for multiple genes.</li> <li>- In sections Gene structure, Exon numbering and in Table 2, exon numbering is now according to MANE. Exon numbering for <i>CDKN2A</i> is now MANE transcripts.</li> <li>- Reference to the P335 specific note added to the section References samples.</li> <li>- In section Interpretation of results, the term “duplication” was replaced with “gain” as this is the term commonly used in ALL diagnostics.</li> <li>- P335 specific note rephrased to better explain the need for male reference samples and result interpretation.</li> <li>- In section <i>IKZF1</i> mutation databases, link to LOVD database updated.</li> <li>- Warning about SNP added for the 295 nt reference probe in Table 1 and 3.</li> <li>- digitalMLPA probemix D007 added to the list of related probemixes.</li> <li>- List of references updated.</li> <li>- Various minor textual and layout changes.</li> </ul> <p>Version C2-01 – 27 October 2021 (04P)</p> <ul style="list-style-type: none"> <li>- Product description rewritten and adapted to a new template.</li> <li>- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).</li> <li>- Intended purpose adapted to a new template</li> <li>- First point in “limitations of the procedure” adjusted to match text in intended use: “in the majority of patients, defects in the <i>IKZF1</i> gene are deletions, but point mutations can occur which will not be detected by MLPA”.</li> <li>- Term “amplification” changed to “heterozygous triplication/homozygous duplication” to be consistent with Copy number status table in section Interpretation of results.</li> <li>- Term “dosage quotient” changed to “final ratio”</li> <li>- Sentence added to section Reference samples: “ For example, when testing DNA samples derived from bone marrow, use reference samples from bone marrow in the same experiment for optimal data normalisation.”</li> </ul> |

- All probemix and positive sample data presented are from internal tests with the P335-C2 version.
- Explanation added to Table 2 for location of the *SHOX* area probe.
- Footnote added to table "Related SALSA MLPLA probemixes" regarding the use of the P202 probemix.
- (More recent) references added in clinical background, Table 2, and sections References and Selected Publications
- Costa Rica added to countries where P335 is registered as in-vitro diagnostic device.
- UK has been added to the list of countries in Europe that accept the CE mark.
- Various minor textual or layout changes.

## Version C1-03 – 04 November 2020 (02P)

- Positive sample information added on page 3.
- Ligation sites of the probes targeting the *PAX5*, *BTG1*, *ETV6*, *CDKN2A*, *RB1* and *EBF1* genes updated according to new version of the NM\_ reference sequences.
- Probemix P018 SHOX added to related probemixes on page 11.
- Selected publications using P335 list shortened to highlight the most informative publications.
- Additional publication added in Table 2 to clarify rearrangements in the PAR1 region.
- Product description adapted to a new template.
- Intended use was updated to a new template.
- Israel added as countries with IVD status.
- Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.
- Various minor textual or layout changes.


## Version C1-02 – 17 January 2019 (04)

- Regulatory status section updated.
- In Table 2 and throughout the document, the NM sequence and the ligation sites for the *IKZF1* gene were updated according to NM\_006060.6.
- Information and the name of the related mix P202-IKZF1 (IKAROS) has been changed to P202-IKZF1-ERG as this probemix has been revised (from version C1 onwards).

## Version C1-01 – 21 March 2018 (04)

- Product description completely rewritten and adapted to a new template
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).
- Small changes of probe lengths in Tables 1 and 2 in order to better reflect the true lengths of the amplification products.
- Added information about the 120 nt target-specific Y chromosome probe to improve data interpretation.
- Added information concerning the minimum percentage of tumour cells needed for reliable data analysis.
- Warning added to Table 2 for probe relying on its specificity on a single nucleotide difference between target and related gene or pseudogene.

**More information:** [www.mrcholland.com](http://www.mrcholland.com); [www.mrcholland.eu](http://www.mrcholland.eu)

|   |   |
|---|---|
|  | MRC Holland bv; Willem Schoutenstraat 1<br>1057 DL, Amsterdam, The Netherlands  |
| E-mail  | <a href="mailto:info@mrcholland.com">info@mrcholland.com</a> (information & technical questions)<br><a href="mailto:order@mrcholland.com">order@mrcholland.com</a> (orders) |
| Phone   | +31 888 657 200   |

|   |  |
|---|--|
|  | EUROPE* <br>MOROCCO<br>ISRAEL<br>COSTA RICA |
|  | ALL OTHER COUNTRIES  |

\*comprising EU (candidate) member states and members of the European Free Trade Association (EFTA), and the UK. The product is for RUO in all other European countries.