

# Product Description

## SALSA® MLPA® Probemix P038 CLL-2

To be used with the MLPA General Protocol.

### Version B2

As compared to version B1, two reference probes have been replaced and one probe has a change in length. For complete product history see page 13.

### Catalogue numbers:

- **P038-025R:** SALSA MLPA Probemix P038 CLL-2, 25 reactions.
- **P038-050R:** SALSA MLPA Probemix P038 CLL-2, 50 reactions.
- **P038-100R:** SALSA MLPA Probemix P038 CLL-2, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see [www.mrcholland.com](http://www.mrcholland.com)).

### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at [www.mrcholland.com](http://www.mrcholland.com).

### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: [www.mrcholland.com](http://www.mrcholland.com). It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

### General information

The SALSA MLPA Probemix P038 CLL-2 is a **research use only (RUO)** assay for the detection of deletions or duplications in the genes/chromosomal regions implicated in B-cell chronic lymphocytic leukemia (B-CLL) such as: 10q23 (*PTEN*), 11q (*ATM*), chromosome 12, 13q14 (*RB1*, *DLEU1/2*), 14q, 17p13 (*TP53*) and chromosome 19. This probemix can also be used to detect the presence of the *NOTCH1* c.7541-7542delCT (p.P2514Rfs\*4), *SF3B1* c.2098A>G (p.K700E) and *MYD88* c.794T>C (p.L265P, also known as p.L252P) point mutations.

B-CLL is the most common hematologic neoplasm in Western countries and results in the progressive accumulation of morphologically mature but functionally incompetent CD5(+) CD23(+) B lymphocytes in the bone marrow, blood, spleen and lymph nodes of the affected person. Chromosomal translocations are rare events in B-CLL. Copy number changes of certain chromosomal regions are, however, frequent. Some of these have been found to be highly prognostic markers of this disease.

More genomic regions affected by copy number alterations (CNAs) in CLL are targeted by the P037 CLL-1 probemix containing probes for 2p (*MYCN*, *ALK*, *REL*), 6q (*TNFAIP3*), 8p (*TNFRSF10A/B*), 8q (*EIF3H*, *MYC*), 9p21 (*CDKN2A/B*), 11q (*ATM*), chromosome 12, 13q14 (*RB1*, *MIR15A*, *DLEU2/7*) and *TP53*. Moreover, the P040 CLL probemix contains a selection of target genes and regions from P037 and P038 allowing copy number determination of 11q, chromosome 12, 13q14 and 17p13. Other related probemixes can be found on page 10.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK470433/>

**This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.**

### Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>  
For NM\_mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Matched Annotation from NCBI and EMBL-EBI (MANE): <http://www.ncbi.nlm.nih.gov/refseq/MANE/>

Tark – Transcript Archive: <http://tark.ensembl.org/>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

### Exon numbering

The *SF3B1*, *MYD88*, *NOTCH1*, *PTEN*, *ATM*, *RB* and *DLEU1* exon numbering used in this P038-B2 CLL-2 product description is the exon numbering from MANE project (release version 1.0) based on MANE Select transcripts, as indicated in Table 2. The *TP53* exon numbering is derived from the LRG\_321 sequence; the exon numbering derived from MANE Select for *TP53* can be found in between brackets in Table 2. As changes to the databases can occur after release of this product description, the NM\_ sequence and exon numbering may not be up-to-date. Please note that exon numbering for the same gene might be different from literature and in other MRC Holland product descriptions, where other resources used for exon numbering are indicated.

### Probemix content

The SALSA MLPA Probemix P038-B2 CLL-2 contains 51 MLPA probes with amplification products between 130 and 500 nucleotides (nt). This includes 35 probes for 10q23 (*PTEN*), 11q (*ATM*), chromosome 12, 13q14 (*RB1*, *DLEU1/2*), 14q, 17p13 (*TP53*) and chromosome 19. Furthermore, this probemix also contains three probes specific for the *NOTCH1* c.7541-7542delCT (p.P2514Rfs\*4), *SF3B1* c.2098A>G (p.K700E) and *MYD88* c.794T>C (p.L265P, also known as p.L252P) mutations which will only generate a signal when the mutation is present. In addition, 13 reference probes are included that target relatively copy number stable regions in CLL. The identity of the genes detected by the reference probes is available in Table 3. Complete probe sequences are available online ([www.mrcholland.com](http://www.mrcholland.com)).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at [www.mrcholland.com](http://www.mrcholland.com).

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

### MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola (2012).

### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

### Required specimens

Extracted DNA, which includes DNA derived from paraffin-embedded tissues, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

### Reference samples

A sufficient number ( $\geq 3$ ) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different healthy individuals without a history of cancer. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Samples from the Coriell Institute and Leibniz Institute DSMZ listed below have been tested with this P038-B2 probemix at MRC-Holland and can be used as positive control samples to detect various CNAs as described in the table below. The quality of cell lines can change; therefore samples should be validated before use.

Sample name	Source	Chromosomal position of CNA (hg18)*	Altered target genes in P038-B2	Expected copy number alteration / mutation
NA20125	Coriell Institute	10q23.31	<i>PTEN</i>	Heterozygous duplication
		14q.32.33	<i>IGHD</i> <sup>†</sup>	Heterozygous deletion
NA09596	Coriell Institute	11q22.3	<i>ATM, RDX</i>	Heterozygous deletion
		14q32.33	<i>IGHD</i> <sup>†</sup>	Heterozygous deletion
NA15099	Coriell Institute	11q22.3-q23.2	<i>ATM, RDX, PPP2R1B, CADM1</i>	Heterozygous duplication
NA00959	Coriell Institute	11q22.3-q23.2	<i>ATM, RDX, PPP2R1B, CADM1</i>	Heterozygous duplication
NA07981	Coriell Institute	12p13.31	<i>CD27</i>	Homozygous duplication
NA05832	Coriell Institute	13q14.2-q14.3	<i>RB1, DLEU2, KCNRG, DLEU1, ATP7B</i>	Heterozygous duplication
NA13721	Coriell Institute	13q14.2-q14.3	<i>RB1, DLEU2, KCNRG, DLEU1, ATP7B</i>	Heterozygous deletion
NA14164	Coriell Institute	13q14.2-q14.3	<i>RB1, DLEU2, KCNRG, DLEU1, ATP7B</i>	Heterozygous deletion
HNT-34	DSMZ	2q33.1	<i>SF3B1</i> c.2098A>G (p.K700E)	Point mutation

\* Indicated chromosomal bands accommodate genes targeted by MLPA probes, however, the whole extent of CNA present in this cell line cannot be determined by this P038-B2 CLL-2 probemix.

<sup>†</sup> *IGHD* is covered by one probe which targets a DNA sequence that undergoes V(D)J recombination in B-cells. Since Coriell Institute samples are cultured as cell lines, it can be expected that a specific recombination will be overrepresented when the sample is derived from a single B cell clone and thus CN-variation will be detected.

### SALSA Binning DNA SD009

The SD009 Binning DNA provided with this probemix can be used for binning of all probes including the three mutation-specific probes (*SF3B1* probe 17802-SP0549-L21900 c.2098A>G=p.K700E, *MYD88* probe 17803-L22642 c.794T>C=p.L265P (also known as c.755T>C=p.L252P) and *NOTCH1* probe 17801-SP0548-L22640 c.7541-7542delCT=p.P2514Rfs\*4). SD009 Binning DNA is a mixture of genomic DNA from healthy individuals and plasmid DNA that contains the target sequence detected by the above mentioned probes. Inclusion of one reaction with 5  $\mu$ l SD009 Binning DNA in initial MLPA experiments is essential as it can be used to aid in data binning of the peak pattern using Coffalyser.Net software. Furthermore, Binning DNA should be included in the experiment whenever changes have been applied to the set-up of the capillary electrophoresis device (e.g. when capillaries have been renewed). Binning DNA should never be used as a reference sample in the MLPA

data analysis, neither should it be used in quantification of mutation signal(s). It is strongly advised that all samples tested are extracted with the same method and derived from the same source of tissue. For further details, please consult the SD009 Binning DNA product description, available online: [www.mrcholland.com](http://www.mrcholland.com).

**This product is for research use only (RUO).**

### Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at [www.mrcholland.com](http://www.mrcholland.com). Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

### Interpretation of results

The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$ . When this criterion is fulfilled, the following cut-off values for the final ratio (FR) of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status indication	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	FR = 0
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication/gain	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication/gain	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

**Please note that these above mentioned final ratios are only valid for germline testing. Final ratios are affected both by percentage of tumour cells and by possible subclonality.**

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in subclonal cases.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.

- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

#### Limitations of the procedure

- In most populations, the major cause of genetic defects in cancer are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P038 CLL-2.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.
- MLPA analysis on tumour samples provides information on the *average* situation in the cells from which the DNA sample was purified. Gains or losses of genomic regions or genes may not be detected if the percentage of tumour cells is low. In addition, subclonality of the aberration affects the final ratio of the corresponding probe. Furthermore, there is always a possibility that one or more reference probes *do* show a copy number alteration in a patient sample, especially in solid tumours with more chaotic karyotypes.

#### Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

#### COSMIC and LOVD mutation databases

We strongly encourage users to deposit positive results in the COSMIC (<https://cancer.sanger.ac.uk/cosmic>) and LOVD (<https://databases.lovd.nl>) mutation databases. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report false positive results due to SNVs and unusual results (e.g., a deletion of *ATM* exons 46 and 63, but not exon 55) to MRC Holland: [info@mrcholland.com](mailto:info@mrcholland.com).

**Table 1. SALSA MLPA Probemix P038-B2-CLL-2**

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18)								
		reference	Mutation-specific	10q	11q	12p/q	13q	14q	17p	19p/q
64-105	Control fragments – see table in probemix content section for more information									
130	Reference probe 00797-L00463	5q31								
136	<b>STAT6 probe</b> 17492-L21300					12q13.3				
142	<b>CD27 probe</b> 16210-L18463					12p13.31				
148	<b>AKT1 probe</b> 17493-L21301							14q32.33		
154	Reference probe 06070-L06371	22q12								
160	<b>RB1 probe</b> 00845-L00378					13q14.2				
167	<b>TP53 probe</b> 01588-L06028								17p13.1	
172	Reference probe 07331-L06968	1q25								
178	<b>DLEU2 probe</b> 17494-L21302					13q14.3				
184	<b>ATM probe</b> 04044-L03849				11q22.3					
192	Reference probe 08181-L22377	10p13								
198	<b>TP53 probe</b> 01996-L22378								17p13.1	
203 + j	<b>IGHD probe</b> 16524-L21379							14q32.33		
211	Reference probe 04533-L03922	2q24								
216 § Δ Ж	<b>NOTCH1 probe</b> 17801-SP0548-L22640		c.7541-7542delCT (p.P2514Rfs*4)							
221	<b>KCNRG probe</b> 04017-L03414					13q14.3				
226	<b>ATM probe</b> 08442-L08312				11q22.3					
234 § Ж	<b>SF3B1 probe</b> 17802-SP0549-L21900		c.2098A>G (p.K700E)							
240	<b>PAH probe</b> 02334-L23006					12q23.2				
245 ~	<b>ATP7B probe</b> 03240-L23007						13q14.3			
252 §	<b>MYD88 probe</b> 17803-L22642		c.794T>C (p.L265P)							
258	Reference probe 11232-L21790	7p15								
265 ¥	<b>CDKN2D probe</b> 23374-L33075									19p13.2
274	<b>TP53 probe</b> 17419-L21141								17p13.1	
281 ~ #	<b>RDX probe</b> 01633-L01171				11q22.3					
288	<b>ATM probe</b> 08445-L21782				11q22.3					
294	Reference probe 01332-L16059	7q11								
301	<b>LDLR probe</b> 02320-L21325									19p13.2
310	<b>IGF1 probe</b> 09545-L22156					12q23.2				
319	<b>TP53 probe</b> 17421-L23008								17p13.1	
328	<b>AKT2 probe</b> 08238-L08116									19q13.2
339 ~	<b>PPP2R1B probe</b> 17496-L21657				11q23.1					
346	<b>TP53 probe</b> 17422-L21144								17p13.1	
353 *	Reference probe 05273-L25208	2p22								
361 ~	<b>CADM1 probe</b> 03816-L15763					11q23.2				
368 #	<b>PTEN probe</b> 03638-L08334				10q23.31					
373	<b>ATM probe</b> 08420-L08326				11q22.3					
382	<b>DLEU1 probe</b> 01589-L17521						13q14.3			
391 «	<b>CCNE1 probe</b> 15145-L22643									19q12
402	<b>DLEU1 probe</b> 00801-L21380						13q14.3			
409 *	Reference probe 19953-L33032	4p16								
418	<b>HMG A2 probe</b> 15074-L16832					12q14.3				
427	<b>MIR498 probe</b> 14249-L15099									19q13.41
436	Reference probe 13363-L14793	15q24								
445	<b>MTA1 probe</b> 16357-L15668							14q32.33		
454	Reference probe 13254-L21630	1p21								
463	<b>ATM probe</b> 02674-L02141				11q22.3					
472 #	<b>PTEN probe</b> 13696-L22095				10q23.31					
480 ±	<b>RB1 probe</b> 01800-L22096						13q14.2			
490	Reference probe 14883-L22098	14q11								
500	Reference probe 15203-L20113	3p12								

\* New in version B2.

¥ Changed in version B2. Minor alteration, no change in sequence detected.

+ This probe targets a DNA sequence that undergoes V(D)J recombination in B-cells, and therefore it can be expected that a specific recombination will be overrepresented when the DNA sample is derived from a single B cell clone and thus CN-variation will be detected.

] Gene name changed to *IGHD* (was *KIAA0125*).

§ Mutation-specific probe. This probe will only generate a signal when the corresponding mutation is present. It has been tested on artificial DNA but not on positive human samples!

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

± SNP rs566016487 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

# This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

**Table 2. P038-B2 target probes arranged according to chromosomal location**

Length (nt)	SALSA MLPA probe	Gene/Exon <sup>a</sup>	Location / Ligation site	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Location (hg18) in kb	Distance to next probe
<b>SF3B1 c.2098A&gt;G = p.K700E mutation</b>						
Recurrent mutations of the splicing factor gene <i>SF3B1</i> occur mainly together with 11q deletions (Wang et al. 2011) and are associated with rapid disease progression and poor overall survival (Quesada et al. 2011, Oscier et al. 2013). <i>SF3B1</i> mutations are recurrent in fludarabine-refractory cases (Rossi et al. 2011). Recurrent somatic mutations of <i>SF3B1</i> have been found in other haematological malignancies including MDS (~28%) and AML (~5%) (Malcovati et al. 2011, Je et al. 2013, Damm et al. 2012). K700E amino acid substitution in the protein encoded by <i>SF3B1</i> is the most frequent mutation found in CLL.						
The ligation site and exon numbering for the SF3B1 (at 2q33.1) probe is according to MANE Select transcript NM_012433.4.						
234 § Ж	17802-SP0549-L21900	<b>SF3B1</b> , exon 15 c.2098A>G (p.K700E)	2127-2126 and 22 nt before exon 15, reverse	GGTCCGA <del>ACTTC</del> -42 nt spanning oligo-TAACTATGCCCC	02-197,975	-
<b>MYD88 c.794T&gt;C = p.L265P (also known as c.755T&gt;C = p.L252P) mutation</b>						
<i>MYD88</i> mutations (mainly c.794T>C = p.L265P substitution, also known as p.L252P) are found in 2-4% Caucasian and 8% Asian CLL patients predominantly in cases with a mutated immunoglobulin gene. Patients with the <i>MYD88</i> mutation are diagnosed at a younger age and there are controversial reports about the prognostic parameters for <i>MYD88</i> -mutated cases (Martinez-Trillos et al. 2014, Qin et al. 2017, Improgo et al. 2019). Moreover, the <i>MYD88</i> c.794T>C = p.L265P mutation has also been found in 29% of ABC subtype of DLBCL (not in other subtypes) (Ngo et al. 2011) and it is also a commonly recurring mutation in patients with Waldenström's macroglobulinemia and can be used in the differential diagnosis of this disease (Treon et al. 2012, Varettoni et al. 2013).						
The ligation site and exon numbering for the MYD88 (at 3p22.2) probe is according to MANE Select transcript NM_002468.5.						
252 §	17803-L22642	<b>MYD88</b> , exon 5 c.794T>C (p.L265P)	795-796	TCAGAAGCGACC-GATCCCCATCAA	03-038,158	43,6 Mb to ref <i>GBE1</i>
<b>NOTCH1 c.7541-7542delCT = p.P2514Rfs*4 (also known as p.2514*fs) mutation</b>						
A 2 bp frameshift deletion of <i>NOTCH1</i> (c.7541-7542delCT, also known as c.7544-7545delCT, p.2514*fs) is found in ~5% of CLL patients and is correlated with unfavourable prognosis (Villamor et al. 2013, Del Giudice et al. 2012, Sportoletti et al. 2010, Rossi et al. 2012 <sup>1</sup> ). <i>NOTCH1</i> mutations are more prevalent in progressive and/or chemorefractory CLL (Fabbri						

Length (nt)	SALSA MLPA probe	Gene/Exon <sup>a</sup>	Location / Ligation site	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Location (hg18) in kb	Distance to next probe
et al. 2011) and might predict transformation of CLL into a clonally related Richter syndrome (Rossi et al. 2012 <sup>2</sup> ). <i>NOTCH1</i> mutations are associated with poor overall survival in mantle cell lymphoma as well (Kridel et al. 2012). The ligation site and exon numbering for the NOTCH1 probe is according to MANE Select transcript NM_017617.5.						
216 § Δ Ж	17801-SP0548-L22640	<b>NOTCH1</b> , exon 34 c.7541-7542delCT (p.P2514Rfs*4)	7804-7805 and 7849-7850	TCACCCCGTCCC-45 nt spanning oligo-TCCGACTGGTCC	09-138,510	-
<b>PTEN loss at 10q23.31</b>						
The loss of the tumour suppressor <i>PTEN</i> is common in various kinds of tumours including leukemia. Low expression level of <i>PTEN</i> predicts poor prognosis in chronic lymphocytic leukemia (Zou et al. 2013). Ligation sites and exon numbering for PTEN probes are according to MANE Select transcript NM_000314.8.						
368 #	03638-L08334	<b>PTEN</b> , exon 5	1250-1251	GGTGTAAATGATA-TGTGCATATTTA	10-089,683	32,4 kb
472 #	13696-L22095	<b>PTEN</b> , exon 9	2171-2170, reverse	AGAGAATTGTTC-CTATAACTGGTA	10-089,715	-
<b>11q deletion</b>						
11q deletion, which results in loss of the <i>ATM</i> gene at 11q22.3, is found in 15-20% of CLL cases. Deletion of 11q22-q23 as well as <i>ATM</i> mutations are associated with aggressive disease and short median survival (Döhner et al. 1997, Neilson et al. 1997, Guarini et al. 2012). Larger 11q deletions also occur, including tumour suppressor genes <i>PPP2R1B</i> , <i>CADM1</i> and <i>RDX</i> (Gunn et al. 2009, Döhner et al. 1997). Ligation sites and exon numbering for ATM probes are according to MANE Select transcript NM_000051.4.						
184	04044-L03849	<b>ATM</b> , exon 1	65-66	GAGGAGTCGGGA-TCTGCGCTGCAG	11-107,599	30,8 kb
373	08420-L08326	<b>ATM</b> , exon 13	2174-2175	AGAAAAGCACCA-GTCCAGTATTGG	11-107,630	66,0 kb
463	02674-L02141	<b>ATM</b> , exon 44	6520-6521	AAGGAACCAGTT-ACCATGAATCAT	11-107,696	15,0 kb
226	08442-L08312	<b>ATM</b> , exon 55	8186-8187	AGAATATGGAAA-TCTGGTACTAT	11-107,711	30,4 kb
288	08445-L21782	<b>ATM</b> , exon 63	9192-9193	ATGAGACTACAA-GAGAAACTGAAA	11-107,741	1,9 Mb
281 ~ #	01633-L01171	<b>RDX</b>	11q22.3	GAAGCAGAACGA-CTTGAAAAGGAG	11-109,613	1,5 Mb
339 ~	17496-L21657	<b>PPP2R1B</b>	11q23.1	TGAGACAGATCT-CCCAGGAGCATA	11-111,137	3,7 Mb
361 ~	03816-L15763	<b>CADM1</b>	11q23.2	CAGCGCATGTCA-TTAGCATCTCAT	11-114,880	-
<b>Trisomy 12</b>						
Trisomy 12 is the third most common cytogenetic abnormality in CLL detected in 10-20% of patients; it confers intermediate or favourable treatment response and overall survival (Hallek et al. 2010). Atypical lymphocyte morphology is observed in some trisomy 12 cases (Matutes et al. 1996). A large portion of trisomy 12 cases harbouring <i>NOTCH1</i> mutations exhibit aggressive CLL (Balatti et al. 2012).						
142	16210-L18463	<b>CD27</b>	12p13.31	CCATCACTGCCA-ATGCTGAGTGTG	12-006,430	4,9 Mb
136	17492-L21300	<b>STAT6</b>	12q13.3	CCGACGCCTTCT-GCTGCAACTTGG	12-055,788	8,7 Mb
418	15074-L16832	<b>HMGA2</b>	12q14.3	GACCCAGGGGAA-GACCCAAAGGCA	12-064,508	36,8 Mb
310	09545-L22156	<b>IGF1</b>	12q23.2	GAGGCTGGAGAT-GTATTGCGCACC	12-101,337	457,9 kb
240	02334-L23006	<b>PAH</b>	12q23.2	GGTCCCAAGAA-CCATTCAAGAGC	12-101,795	-
<b>13q14 deletion</b>						
Interstitial deletion at 13q14 is the most common (~50%) chromosomal aberration in CLL. The <i>DLEU/miR15A/16-1</i> cluster, as well as the <i>RB1</i> gene, are important tumour suppressor candidates within the 13q14 deletion region (Klein et al. 2010, Palamarchuk et al. 2010). Deletion of 13q14 represents a CLL group with the best prognosis and, when it is the sole abnormality, also with the highest overall survival. In P037 CLL-1 probemix more probes on 13q14 are included to define the deletion size. Ligation sites and exon numbering for RB1 (13q14.2) probes are according to MANE Select transcript NM_000321.3. Ligation sites and exon numbering for DLEU1 (13q14.3) probes are according to NR_109973.1.						
160	00845-L00378	<b>RB1</b> , exon 17	1689-1690	CTTGATTCTGGA-ACAGATTTGTCT	13-047,853	99,1 kb
480 ±	01800-L22096	<b>RB1</b> , exon 27	3270-3271	GAGTCCTGATAA-CCCAGGCCTGTC	13-047,953	1,5 Mb
178	17494-L21302	<b>DLEU2</b>	13q14.3	CCAATCTCAAGC-CTGTACATTGTT	13-049,455	32,8 kb
221	04017-L03414	<b>KCNRG</b>	13q14.3	CTCTAGTTTGAA-GTGAGGGAAGAA	13-049,488	89,2 kb
402	00801-L21380	<b>DLEU1</b> , exon 3	474-475	GAAGAACAGAAC-CTTCAGGAATTG	13-049,577	205,5 kb
382	01589-L17521	<b>DLEU1</b> , downstream	203 kb after exon 3	CCTTTTAATAGG-ATCTCTCTGGA	13-049,782	1,7 Mb
245 ~	03240-L23007	<b>ATP7B</b>	13q14.3	GGTTGCCTTGAT-GGCAGGAAAGGC	13-051,441	-
<b>14q deletion</b>						



Length (nt)	SALSA MLPA probe	Gene/Exon <sup>a</sup>	Location / Ligation site	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Location (hg18) in kb	Distance to next probe
14q24.1-q32.33 loss is detected in ~8% of CLL patients (Kay et al. 2010, Pospisilova et al. 2007, Wlodarska et al. 2007). 14q deletions are frequently associated with trisomy 12, <i>NOTCH1</i> mutations and unmutated <i>IGHV</i> genes (Cosson et al. 2014). The co-occurrence of 14q deletion and trisomy 12 is significantly associated with shorter overall survival and time to treatment in comparison to non 14p deletion cases (Reindl et al. 2010).						
148	17493-L21301	<b>AKT1</b>	14q32.33	AGAATGATGGCA-CCTTCATTGGCT	14-104,318	665,3 kb
445	16357-L15668	<b>MTA1</b>	14q32.33	ACAGACGGCCAA-TGGGAACGTGGA	14-104,983	400,1 kb
203 + j	16524-L21379	<b>IGHD</b>	14q32.33	TCCGTGACTGTC-ACCTGGTACATG	14-105,383	-
<b>17p deletion</b>						
17p deletions are detected in 5-10% of newly diagnosed CLL resulting in a loss of <i>TP53</i> tumour suppressor gene at 17p13.1. Del(17p) and also <i>TP53</i> mutations are associated with a more aggressive clinical course, worse prognosis and short overall survival, and belong to ultra-high risk CLL (Mougalian and O'Brien 2011). <i>TP53</i> deletion/mutations predict no response to treatment with purine analogues and are thus important for therapy selection (Stilgenbauer et al. 2010, Schetelig et al. 2008, Dreger et al. 2010).						
Ligation sites for <i>TP53</i> (17p13.1) probes are according to MANE Select transcript NM_000546.6, and the exon numbering is according to LRG_321. The MANE Select transcript NM_000546.6 exon numbering is indicated between brackets.						
346	17422-L21144	<b>TP53</b> , exon 10	1188-1189	TGAGGCCCTTGGAACTCAAGGATGC	17-007,515	3,6 kb
274	17419-L21141	<b>TP53</b> , exon 6 (7)	831-832	CTCTGACTGTAC-CACCATCCACTA	17-007,518	0,7 kb
319	17421-L23008	<b>TP53</b> , exon 5 (6)	735-736	TATCCGAGTGGAAAGGAAATTTGCG	17-007,519	1,7 kb
198	01996-L22378	<b>TP53</b> , exon 2a (2)	117-118	CTCTTGCAGCAG-CCAGACTGCCTT	17-007,521	10,8 kb
168	01588-L06028	<b>TP53</b> , exon 1	58-59	TCCGGGGACACT-TTGC GTTCGGGC	17-007,531	-
<b>Trisomy 19</b>						
Trisomy 19 has been detected in 1.6-5% of CLL cases and is concurrent with trisomy 12 (Gunnarsson et al. 2011, Groenen et al. 2011). In addition to trisomy 12, trisomy 19 has been associated with mutated <i>IGHV</i> genes (Sellmann et al. 2007).						
265	23374-L33075	<b>CDKN2D</b>	19p13.2	AGCACGGGGCTG-ATGTCAACGTGC	19-010,539	543,5 kb
301	02320-L21325	<b>LDLR</b>	19p13.2	AATGACCTTAAG-ATCGGCTACGAG	19-011,082	23,9 Mb
391 «	15145-L22643	<b>CCNE1</b>	19q12	TTGTCTGAACAA-AATAGGGCTTCT	19-035,006	10,4 Mb
328	08238-L08116	<b>AKT2</b>	19q13.2	CGACTATCTCAA-ACTCCTTGCCAA	19-045,440	13,4 Mb
427	14249-L15099	<b>MIR498</b>	19q13.41	GCTGTGATTTC-AGCCAGGGGGCG	19-058,869	-

<sup>a</sup> See section Exon numbering on page 2 for more information.

<sup>b</sup> Only partial probe sequences are shown. Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com). Please notify us of any mistakes: [info@mrcholland.com](mailto:info@mrcholland.com).

§ Mutation-specific probe. This probe will only generate a signal when the corresponding mutation is present. It has been tested on artificial DNA but not on positive human samples!

+ This probe targets a DNA sequence that undergoes V(D)J recombination in B-cells, and therefore it can be expected that a specific recombination will be overrepresented when the DNA sample is derived from a single B cell clone and thus CN-variation will be detected.

j Gene name changed to *IGHD* (was *KIAA0125*).

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

± SNP rs566016487 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

# This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

**Table 3. P038-B2 reference probes arranged according to chromosomal location.**

Length (nt)	SALSA MLPA probe	Gene	Location (hg18)	Location (hg18) in kb
454	13254-L21630	<i>COL11A1</i>	1p21	01-103,251
172	07331-L06968	<i>SERPINC1</i>	1q25	01-172,145
353	05273-L25208	<i>SPAST</i>	2p22	02-032,215
211	04533-L03922	<i>SCN1A</i>	2q24	02-166,607
500	15203-L20113	<i>GBE1</i>	3p12	03-081,775
409	19953-L33032	<i>EVC</i>	4p16	04-005,857
130	00797-L00463	<i>IL4</i>	5q31	05-132,038
258	11232-L21790	<i>GARS</i>	7p15	07-030,609
294	01332-L16059	<i>STX1A</i>	7q11	07-072,755
192	08181-L22377	<i>OPTN</i>	10p13	10-013,208
490	14883-L22098	<i>RPGRIP1</i>	14q11	14-020,873
436	13363-L14793	<i>SEMA7A</i>	15q24	15-072,490
154	06070-L06371	<i>CACNG2</i>	22q12	22-035,290

Complete probe sequences are available at [www.mlpa.com](http://www.mlpa.com).

### Related SALSA MLPA probemixes

<b>P037 CLL-1</b>	Contains probes for 2p, 6q, 8p/q, 9p21, <i>ATM</i> , chromosome 12, 13q14 and <i>TP53</i>
<b>P040 CLL</b>	Contains probes for 11q, 12p/q, 13q14 and <i>TP53</i>
<b>P041 and P042 ATM</b>	Contain probes for all <i>ATM</i> exons on 11q
<b>P056 TP53</b>	Contains probes for all <i>TP53</i> exons
<b>P047 RB1</b>	Contains more probes for 13q14.2
<b>P252 NB mix 2</b>	Contains more probes for the 2p region
<b>P323 CDK4-HMGA2-MDM2</b>	Contains more probes for chromosome 12
<b>P377 Hematologic Malignancies</b>	Contains more probes for 6q, 8q, 9p, 11q, 12p/q, 13q14 and <i>TP53</i>

### References

- Balatti V et al. (2012). NOTCH1 mutations in CLL associated with trisomy 12. *Blood*. 119:329-31.
- Cosson A et al. (2014). 14q deletions are associated with trisomy 12, NOTCH1 mutations and unmutated IGHV genes in chronic lymphocytic leukemia and small lymphocytic lymphoma. *Genes Chromosomes Cancer*. 53:657-66.
- Damm F et al. (2012). SF3B1 mutations in myelodysplastic syndromes: clinical associations and prognostic implications. *Leukemia*. 26:1137-40.
- Del Giudice I et al. (2012). NOTCH1 mutations in +12 chronic lymphocytic leukemia (CLL) confer an unfavorable prognosis, induce a distinctive transcriptional profiling and refine the intermediate prognosis of +12 CLL. *Haematologica*. 97:437-41.
- Döhner H et al. (1997). 11q deletions identify a new subset of B-cell chronic lymphocytic leukemia characterized by extensive nodal involvement and inferior prognosis. *Blood*. 89:2516-22.
- Dreger P et al. (2010). Allogeneic stem cell transplantation provides durable disease control in poor-risk chronic lymphocytic leukemia: long-term clinical and MRD results of the German CLL Study Group CLL3X trial. *Blood*. 116:2438-47.
- Fabbri G et al. (2011). Analysis of the chronic lymphocytic leukemia coding genome: role of NOTCH1 mutational activation. *J Exp Med*. 208:1389-401.
- Groenen PJTA et al. (2011). High prevalence of adverse prognostic genetic aberrations and unmutated IGHV genes in small lymphocytic lymphoma as compared to chronic lymphocytic leukemia. *Journal of Hematopathology*. 4:189-97.
- Guarini A et al. (2012). *ATM* gene alterations in chronic lymphocytic leukemia patients induce a distinct gene expression profile and predict disease progression. *Haematologica*. 97:47-55.
- Gunn SR et al. (2009). Atypical 11q deletions identified by array CGH may be missed by FISH panels for prognostic markers in chronic lymphocytic leukemia. *Leukemia*. 23:1011-7.

- Gunnarsson R et al. (2011). Array-based genomic screening at diagnosis and during follow-up in chronic lymphocytic leukemia. *Haematologica*. 96:1161-9.
- Hallek M et al. (2010). Addition of rituximab to fludarabine and cyclophosphamide in patients with chronic lymphocytic leukaemia: a randomised, open-label, phase 3 trial. *Lancet*. 376:1164-74.
- Hömig-Hölzel C and Savola S (2012). Multiplex ligation-dependent probe amplification (MLPA) in tumor diagnostics and prognostics. *Diagn Mol Pathol*. 21:189-206.
- Improgo MR et al. (2019). MYD88 L265P mutations identify a prognostic gene expression signature and a pathway for targeted inhibition in CLL. *Br J Haematol*. 184:925-36.
- Je EM et al. (2013). Mutational analysis of splicing machinery genes SF3B1, U2AF1 and SRSF2 in myelodysplasia and other common tumors. *Int J Cancer*. 133:260-5.
- Kay NE et al. (2010). Progressive but previously untreated CLL patients with greater array CGH complexity exhibit a less durable response to chemoimmunotherapy. *Cancer Genet Cytogenet*. 203:161-8.
- Klein U et al. (2010). The DLEU2/miR-15a/16-1 cluster controls B cell proliferation and its deletion leads to chronic lymphocytic leukemia. *Cancer Cell*. 17:28-40.
- Kridel R et al. (2012). Whole transcriptome sequencing reveals recurrent NOTCH1 mutations in mantle cell lymphoma. *Blood*. 119:1963-71.
- Malcovati L et al. (2011). Clinical significance of SF3B1 mutations in myelodysplastic syndromes and myelodysplastic/myeloproliferative neoplasms. *Blood*. 118:6239-46.
- Martinez-Trillos A et al. (2014). Mutations in TLR/MYD88 pathway identify a subset of young chronic lymphocytic leukemia patients with favorable outcome. *Blood*. 123:3790-6.
- Matutes E et al. (1996). Trisomy 12 defines a group of CLL with atypical morphology: correlation between cytogenetic, clinical and laboratory features in 544 patients. *Br J Haematol*. 92:382-8.
- Mougalian SS, O'Brien S (2011). Adverse prognostic features in chronic lymphocytic leukemia. *Oncology (Williston Park)*. 25:692-6, 699.
- Neilson JR et al. (1997). Deletions at 11q identify a subset of patients with typical CLL who show consistent disease progression and reduced survival. *Leukemia*. 11:1929-32.
- Ngo VN et al. (2011). Oncogenically active MYD88 mutations in human lymphoma. *Nature*. 470:115-9.
- Oscier DG et al. (2013). The clinical significance of NOTCH1 and SF3B1 mutations in the UK LRF CLL4 trial. *Blood*. 121:468-75.
- Palamarchuk A et al. (2010). 13q14 deletions in CLL involve cooperating tumor suppressors. *Blood*. 115:3916-22.
- Pospisilova H et al. (2007). Interstitial del(14)(q) involving IGH: a novel recurrent aberration in B-NHL. *Leukemia*. 21:2079-83.
- Qin SC et al. (2017). MYD88 mutations predict unfavorable prognosis in Chronic Lymphocytic Leukemia patients with mutated IGHV gene. *Blood Cancer J*. 7:651.
- Quesada V et al. (2011). Exome sequencing identifies recurrent mutations of the splicing factor SF3B1 gene in chronic lymphocytic leukemia. *Nat Genet*. 44:47-52.
- Reindl L et al. (2010). Biological and clinical characterization of recurrent 14q deletions in CLL and other mature B-cell neoplasms. *Br J Haematol*. 151:25-36.
- Rossi D et al. (2011). Mutations of the SF3B1 splicing factor in chronic lymphocytic leukemia: association with progression and fludarabine-refractoriness. *Blood*. 118:6904-8.
- Rossi D et al. (2012)<sup>1</sup>. Mutations of NOTCH1 are an independent predictor of survival in chronic lymphocytic leukemia. *Blood*. 119:521-9.
- Rossi D et al. (2012)<sup>2</sup>. Different impact of NOTCH1 and SF3B1 mutations on the risk of chronic lymphocytic leukemia transformation to Richter syndrome. *Br J Haematol*. 158:426-9.
- Schetelig J et al. (2008). Allogeneic hematopoietic stem-cell transplantation for chronic lymphocytic leukemia with 17p deletion: a retrospective European Group for Blood and Marrow Transplantation analysis. *J Clin Oncol*. 26:5094-100.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res*. 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat*. 28:205.
- Sellmann L et al. (2007). Trisomy 19 is associated with trisomy 12 and mutated IGHV genes in B-chronic lymphocytic leukaemia. *Br J Haematol*. 138:217-20.
- Sportoletti P et al. (2010). NOTCH1 PEST domain mutation is an adverse prognostic factor in B-CLL. *Br J Haematol*. 151:404-6.

- Stilgenbauer S et al. (2010). Understanding and managing ultra high-risk chronic lymphocytic leukemia. *Hematology Am Soc Hematol Educ Program*. 2010:481-8.
- Treon SP et al. (2012). MYD88 L265P somatic mutation in Waldenstrom's macroglobulinemia. *N Engl J Med*. 367:826-33.
- Varettoni M et al. (2013). Prevalence and clinical significance of the MYD88 (L265P) somatic mutation in Waldenstrom's macroglobulinemia and related lymphoid neoplasms. *Blood*. 121:2522-8.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem*. 421:799-801.
- Villamor N et al. (2013). NOTCH1 mutations identify a genetic subgroup of chronic lymphocytic leukemia patients with high risk of transformation and poor outcome. *Leukemia*. 27:1100-6.
- Wang L et al. (2011). SF3B1 and other novel cancer genes in chronic lymphocytic leukemia. *N Engl J Med*. 365:2497-506.
- Wlodarska I et al. (2007). Telomeric IGH losses detectable by fluorescence in situ hybridization in chronic lymphocytic leukemia reflect somatic VH recombination events. *J Mol Diagn*. 9:47-54.
- Zou ZJ et al. (2013). Low expression level of phosphatase and tensin homolog deleted on chromosome ten predicts poor prognosis in chronic lymphocytic leukemia. *Leuk Lymphoma*. 54:1159-64.

### Selected publications using SALSA MLPA Probemix P038 CLL-2

- Abdool A et al. (2010). Detection, analysis and clinical validation of chromosomal aberrations by multiplex ligation-dependent probe amplification in chronic leukemia. *PLoS One*. 5:e15407.
- Al Zaabi EA et al. (2010). Multiplex ligation-dependent probe amplification versus multiprobe fluorescence in situ hybridization to detect genomic aberrations in chronic lymphocytic leukemia: a tertiary center experience. *J Mol Diagn*. 12:197-203.
- Balla B et al. (2023). Copy Number Variations and Gene Mutations Identified by Multiplex Ligation-Dependent Probe Amplification in Romanian Chronic Lymphocytic Leukemia Patients. *Journal of Personalized Medicine*. 13:1239.
- Coll-Mulet L et al. (2008). Multiplex ligation-dependent probe amplification for detection of genomic alterations in chronic lymphocytic leukaemia. *Br J Haematol*. 142:793-801.
- Delfau-Larue M-H et al. (2015). High-dose cytarabine does not overcome the adverse prognostic value of CDKN2A and TP53 deletions in mantle cell lymphoma. *Blood*. 126:604-11
- Eid OM et al. (2021). Evaluation of MLPA as a comprehensive molecular cytogenetic tool to detect cytogenetic markers of chronic lymphocytic leukemia in Egyptian patients. *J Genet Eng Biotechnol*. 19:98.
- Fabris S et al. (2011). Multiplex ligation-dependent probe amplification and fluorescence in situ hybridization to detect chromosomal abnormalities in chronic lymphocytic leukemia: a comparative study. *Genes Chromosomes Cancer*. 50:726-34.
- Gángó A et al. (2022). Morphologic and molecular analysis of Richter syndrome in chronic lymphocytic leukaemia patients treated with ibrutinib or venetoclax. *Pathology*. 54:95-103.
- Groenen PJTA et al. (2011). High prevalence of adverse prognostic genetic aberrations and unmutated IGHV genes in small lymphocytic lymphoma as compared to chronic lymphocytic leukemia. *Journal of Hematopathology*. 4:189-97.
- Simonin M et al. (2023). Prognostic value and oncogenic landscape of TP53 alterations in adult and pediatric T-ALL. *Blood*. 141:1353-1358.
- Srinivasan VK et al. (2020). Genomic alterations in chronic lymphocytic leukemia and their correlation with clinico-hematological parameters and disease progression. *Blood Res*. 55:131-138.
- Stevens-Kroef M et al. (2009). Identification of chromosomal abnormalities relevant to prognosis in chronic lymphocytic leukemia using multiplex ligation-dependent probe amplification. *Cancer Genet Cytogenet*. 195:97-104.
- Stevens-Kroef MJ et al. (2014). Identification of prognostic relevant chromosomal abnormalities in chronic lymphocytic leukemia using microarray-based genomic profiling. *Mol Cytogenet*. 7:3.
- Véronèse L et al. (2013). Contribution of MLPA to routine diagnostic testing of recurrent genomic aberrations in chronic lymphocytic leukemia. *Cancer Genet*. 206:19-25.

P038 product history	
Version	Modification
B2	Two reference probes have been replaced and one probe has a minor change in length but not in the sequence targeted.
B1	19 target probes have been replaced and 11 new probes added. Moreover, 12 reference probes have been replaced. In addition, the 88 and 96 nt control fragments have been replaced (QDX2).
A2	Extra control fragments at 88-96-100-105 nt have been added.
A1	First Release

Implemented changes in the product description
<p>Version B2-01 – 16 January 2024 (04P)</p> <ul style="list-style-type: none"> <li>- Product description adapted to a new product version (version number changed, changes in Tables 1-3).</li> <li>- Probe KIAA0125 (16524-L21379) was renamed to IGHD throughout the document.</li> <li>- Footnotes added for IGHD probe 16524-L21379 in the positive sample table on page 3, Tables 1 and 2.</li> <li>- Samples L-363, U-266, NA0330 and NA08123 removed from Positive control DNA samples table.</li> <li>- "Heterozygous duplication" corrected to "Heterozygous deletion" for NA13721 in Positive control DNA samples table.</li> <li>- Footnote "« Probe located in or near a GC-rich region" removed for 01332-L16059, 15074-L16832 and 16357-L15668 probes in Tables 1-3.</li> <li>- 17496-L21657 probe length changed from 338 to 339 in Table 1 and 2.</li> <li>- 01589-L17521 probe length changed from 381 to 382 in Table 1 and 2.</li> <li>- Various minor textual changes.</li> </ul> <p>Version B1-02 – 12 April 2023 (04P)</p> <ul style="list-style-type: none"> <li>- Product description rewritten and adapted to a new template.</li> <li>- Transcript numbers of the SF3B1, MYD88, ATM, RB1 and TP53 genes have been updated according to the new version of MANE Select.</li> <li>- Ligation sites of the probes targeting SF3B1, MYD88, ATM, RB1 and TP53 have been updated according to their respective NM_ sequence according to the new version of the MANE Select.</li> <li>- The annotation for the NOTCH1 mutation-specific probe was changed from NOTCH1 p.P2514*fs to NOTCH1 p.P251Rfs*4 and for the MYD88 mutation-specific probe alternative names c.755T&gt;C and p.L252P were added next to existing annotations c.794T&gt;C and p.L265P.</li> <li>- Footnotes were added to Tables 1 and 2 to indicate flanking probes targeting 11q.</li> <li>- New references added in the Selected publications section.</li> <li>- The related SALSA Probemixes section was updated (now includes P377 Hematologic malignancies).</li> <li>- Various minor textual or layout changes.</li> </ul> <p>Version B1-01 – 06 February 2020 (02P)</p> <ul style="list-style-type: none"> <li>- Product description rewritten and adapted to a new template.</li> <li>- Product description separated from P037 probemix description.</li> <li>- Various minor textual or layout changes.</li> <li>- Ligation sites of the probes targeting the SF3B1, MYD88 genes updated according to new version of the NM_ reference sequence in Table 2a.</li> <li>- Ligation sites and exon numbering added for DLEU1 probes in Table 2a.</li> <li>- NM_ reference sequence version updated for NOTCH1 and PTEN probes in Table 2a.</li> <li>- Warning added to Table 1 and Table 2a for probe specificity relying on a single nucleotide difference between target gene and related gene or pseudogene.</li> <li>- New references added in Table 2a and on pages 8-10.</li> <li>- For uniformity, the chromosomal locations and bands in this document are now all based on hg18 (NCBI36).</li> </ul>

More information: <a href="http://www.mrcholland.com">www.mrcholland.com</a> ; <a href="http://www.mrcholland.eu">www.mrcholland.eu</a>	
	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	<a href="mailto:info@mrcholland.com">info@mrcholland.com</a> (information & technical questions) <a href="mailto:order@mrcholland.com">order@mrcholland.com</a> (orders)
Phone	+31 888 657 200