

Product Description

SALSA® MLPA® Probemix P265-B4 PROC

To be used with the MLPA General Protocol.

Version B4

For complete product history see page 6.

Catalogue numbers:

- **P265-025R:** SALSA MLPA Probemix P265 PROC, 25 reactions.
- **P265-050R:** SALSA MLPA Probemix P265 PROC, 50 reactions.
- **P265-100R:** SALSA MLPA Probemix P265 PROC, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P265 PROC is a **research use only (RUO)** assay for the detection of deletions or duplications in the *PROC* gene, which is associated with protein C deficiency.

Protein C deficiency is usually inherited as an autosomal dominant trait that is associated with an increased risk of venous thrombosis. A very rare severe homozygous (or compound heterozygous) form of protein C deficiency was first described as a cause of neonatal purpura fulminans, resulting in rapid death caused by disseminated intravascular coagulation (DIC).

Protein C (PROC) is a vitamin K-dependent serine protease zymogen. Purified human activated protein C selectively destroys factors Va and VIII in human plasma and thus has an important anticoagulant role. Regulatory DNA elements responsible for *PROC* gene expression have been identified in the upstream promoter region and first (untranslated) exon of the *PROC* gene.

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>
For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>
Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering

The *PROC* exon numbering used in this P265-B4 PROC product description is the exon numbering from the LRG_599 sequence. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P265-B4 PROC contains 21 MLPA probes with amplification products between 135 and 310 nucleotides (nt). This includes ten probes for the *PROC* gene, with one probe for all nine exons and two probes for exon 1. In addition, eleven reference probes are included which detect eleven different autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola (2012).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of protein C deficiency. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	$FR = 0$
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun

the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *PROC* gene are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P265 PROC.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

PROC mutation database

<https://databases.lovd.nl/shared/genes/PROC>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *PROC* exons 6 and 8 but not exon 7) to MRC Holland: info@mrcholland.com.

Table 1. SALSA MLPA Probemix P265-B4 PROC

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a	
		Reference	PROC
64-105	Control fragments – see table in probemix content section for more information		
135 Δ	Reference probe 16316-L21434	3q	
142	Reference probe 13375-L14832	6q	
157	PROC probe 15027-L16778		Exon 4
166	Reference probe 09890-L10303	16p	
175	PROC probe 08211-L08144		Exon 5
181	PROC probe 09196-L10264		Exon 2
188	PROC probe 08215-L12576		Exon 9
196	Reference probe 19754-L26537	9q	
205	PROC probe 08207-L16917		Exon 1
211	Reference probe 02518-L01949	17q	
220	Reference probe 07405-L07052	12q	
229	PROC probe 08208-L08141		Exon 1
238	Reference probe 07030-L06641	14q	
247	PROC probe 08210-L08143		Exon 3
256	PROC probe 08212-L08145		Exon 6
265	Reference probe 17063-L20141	7q	
273	PROC probe 08213-L08146		Exon 7
283	PROC probe 08214-L08147		Exon 8
292	Reference probe 20858-L28876	1q	
301	Reference probe 14479-L16199	4q	
310	Reference probe 16509-L18970	5q	

^a See section Exon numbering on page 1 for more information.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2. PROC probes arranged according to chromosomal location

Length (nt)	SALSA MLPA probe	PROC exon ^a	Ligation site NM_000312.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	88-90 (Exon 2)		
205	08207-L16917	Exon 1	7 nt before exon 1	ATATTTGTGGTT-ATGGATTAAGTC	0.3 kb
229	08208-L08141	Exon 1	176 nt after exon 1	CCACAGTCTCAG-GTCCCTTTGCCA	1.3 kb
181	09196-L10264	Exon 2	102-103	TGGCAGCTCACA-AGCCTCTGCTG	1.4 kb
247	08210-L08143	Exon 3	280-281	AGGAGATCTGTG-CTTCGAGGAGG	1.5 kb
157	15027-L16778	Exon 4	340-341	CCTTCTGGTCCA-AGCACGTCGGTG	0.1 kb
175	08211-L08144	Exon 5	364-365	GTGACCAAGTGCT-TGGTCTTGCCCT	0.3 kb
256	08212-L08145	Exon 6	542-543	GCATTACTGCCT-AGAGGAGGTGGG	2.8 kb
273	08213-L08146	Exon 7	699-700	ACAGAAGACCAA-GAAGACCAAGTA	1.1 kb
283	08214-L08147	Exon 8	2 nt after exon 8	GTCAGGCTTGGT-ATGGGCTGGAGC	1.4 kb
188	08215-L12576	Exon 9	1198-1199	TGGTCCCGCACA-ATGAGTGCAGCG	
		<i>stop codon</i>	1471-1473 (Exon 9)		

^a See section Exon numbering on page 1 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Related SALSA MLPA probemixes

- P112 PROS1 Contains probes for the *PROS1* gene, involved in protein S deficiency associated with venous thrombosis.
- P227 SERPINC1 Contains probes for the *SERPINC1* gene, involved in antithrombin (III) deficiency.

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P265 PROC

- Caspers M et al. (2012). Deficiencies of antithrombin, protein C and protein S practical experience in genetic analysis of a large patient cohort. *Thromb Haemost* 108:247-257.
- Fidalgo T et al. (2015). Familial thrombotic risk based on the genetic background of Protein C Deficiency in a Portuguese Study. *Eur J Haematol* 95(4):294-307.
- Limperger V et al. (2014). Clinical and laboratory characteristics of children with venous thromboembolism and protein C-deficiency: an observational Israeli-German cohort study. *Br J Haematol* 163(3):385-393.
- Martos L et al. (2019). Identification of 58 mutations (26 Novel) in 94 of 109 symptomatic Spanish probands with protein C deficiency. *Thromb Haemost*, 119(09), 1409-1418.
- Reda S et al. (2021). PC Deficiency Testing: Thrombin-Thrombomodulin as PC Activator and Aptamer-Based Enzyme Capturing Increase Diagnostic Accuracy. *Front Cardiovasc Med*, 8.
- Togashi T et al. (2020). Molecular genetic analysis of inherited protein C deficiency caused by the novel large deletion across two exons of PROC. *Thromb Res*, 188, 115-118.
- Vrtel P et al. (2022). Detection of Unknown and Rare Pathogenic Variants in Antithrombin, Protein C and Protein S Deficiency Using High-Throughput Targeted Sequencing. *Diagnostics*, 12(5), 1060.
- Winther-Larsen A et al. (2020). Protein C deficiency; PROC gene variants in a Danish population. *Thromb Res*, 185, 153-159.
- Wypasek E et al. (2017). Genetic characterization of antithrombin, protein C, and protein S deficiencies in Polish patients. *Pol Arch Intern Med*, 127(7-8), 512-523.

P265 product history	
Version	Modification
B4	Two reference probes have been replaced and one reference probe removed.
B3	One reference probe has been replaced and four reference probes have been added.
B2	Two reference probes have been replaced and the control fragments have been adjusted (QDX2).
B1	One probe detecting the <i>PROC</i> gene and one reference probe have been added. Four reference probes have been replaced.
A1	First release.

Implemented changes in the product description


Version B4-02 – 28 September 2022 (04P)

- Product description rewritten and adapted to a new template.
- Ligation sites of the probes targeting the *PROC* gene updated according to new version of the NM_ reference sequence.

Version B4-01 – 18 April 2019 (01P)

- Product description restructured and adapted to a new template.
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).

More information: www.mrcholland.com; www.mrcholland.eu

	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mrcholland.com (information & technical questions) order@mrcholland.com (orders)
Phone	+31 888 657 200