

# Product Description SALSA<sup>®</sup> MLPA<sup>®</sup> Probemix P147-C1 1p36

To be used with the MLPA General Protocol.

#### Version C1

As compared to version B2, seven target probes were replaced, nine reference probes were replaced and one reference probe was added. Two target probes were changed in length but not in the sequence detected. For complete product history see page 7.

#### Catalogue numbers:

- P147-025R: SALSA MLPA Probemix P147 1p36, 25 reactions.
- **P147-050R:** SALSA MLPA Probemix P147 1p36, 50 reactions.
- P147-100R: SALSA MLPA Probemix P147 1p36, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

#### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

#### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

#### General information

The SALSA MLPA Probemix P147 1p36 is a **research use only (RUO)** assay for the detection of deletions or duplications in the 1p36 subtelomeric region, which are associated with the 1p36 deletion syndrome, also referred to as "monosomy 1p36". The P147 1p36 can be used to confirm and further characterise abnormalities detected by P036 Subtelomeres Mix 1 and/or P070 Subtelomeres Mix 2B.

1p36 deletion syndrome is a chromosome disorder where the end of the short arm of one of the two chromosomes is lost. The breakpoints for this cytogenetic syndrome are variable and range from bands 1p36.13 to 1p36.33 (Slavotinek et al. 1999). It is considered to be one of the most common chromosome terminal deletion syndromes and is a frequent cause of intellectual disability, although there is phenotypic variability based on the location and the extent of the deletions (Jordan et al. 2015). The incidence has been estimated to be 1 in 5,000 to 1 in 10,000 live born children (Gajecka et al. 2007). To date, more affected females than males have been reported.

# This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

#### Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene For NM\_ mRNA reference sequences: http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide Locus Reference Genomic (LRG) database: http://www.lrg-sequence.org/ Matched Annotation from NCBI and EMBL-EBI (MANE): https://www.ncbi.nlm.nih.gov/refseq/MANE

#### **Probemix content**

The SALSA MLPA Probemix P147-C1 1p36 contains 48 MLPA probes with amplification products between 124 and 500 nucleotides (nt). This includes 37 probes for the 1p36 subtelomeric region. In addition, 11 reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

## MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

#### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

#### **Required specimens**

Extracted DNA, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

#### **Reference samples**

A sufficient number (≥3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of 1p36 deletion syndrome. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

#### **Positive control DNA samples**

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

#### Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

#### Interpretation of results

The standard deviation of each individual reference probe over all the reference samples should be  $\leq 0.10$  and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:



Copy number status	Final ratio (FR)
Normal	0.80 < FR < 1.20
Homozygous deletion	FR = 0
Heterozygous deletion	0.40 < FR < 0.65
Heterozygous duplication	1.30 < FR < 1.65
Heterozygous triplication/homozygous duplication	1.75 < FR < 2.15
Ambiguous copy number	All other values

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- <u>Arranging probes</u> according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- <u>False positive results</u>: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- <u>Normal copy number variation</u> in healthy individuals is described in the database of genomic variants: <u>http://dgv.tcag.ca/dgv/app/home</u>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- <u>Copy number changes detected by reference probes</u> or flanking probes are unlikely to have any relation to the condition tested for.
- <u>False results can be obtained if one or more peaks are off-scale</u>. For example, a duplication of one or more
  exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale
  peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net
  software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun
  the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of
  sample by diluting PCR products.
- Interpretation of abnormal copy number findings in subtelomeric regions is complicated. Subtelomeric copy number changes can also occur in unaffected individuals and the effect of a deletion or duplication will depend on the genes involved. A considerable number of abnormalities detected by a single probe may not be the cause of any phenotypic effect but can be due to a rare polymorphism or a copy number change which is also present in one of the parents. For some chromosome arms, even large subtelomeric deletions or duplications (>1 Mb) can be inherited without a clear phenotypic effect. For all abnormalities detected, we strongly recommend testing parents to determine whether the copy number aberration in the patient is de novo.

# P147 specific note:

Please note that deletions or duplications of part of 1p36 can be present in healthy persons! In case of positive results, it is therefore strongly recommended to also test the parents. As an example, a duplication of approximately 1 Mb was detected in a DNA sample from a healthy individual with the P036 and P070 Subtelomeres probemixes. This duplication was further characterized with this P147 probemix. The duplication included all probes from *TNFRSF4* up to *GABRD* (Kathleen Claes, Ghent, personal communication).

# Limitations of the procedure

- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false
  positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe
  signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe
  oligonucleotide to the sample DNA.

# Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

## Database of Genomic Variants

We strongly encourage users to deposit positive results in the Database of Genomic Variants (http://dgv.tcag.ca/dgv/app/home). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results to MRC Holland: info@mrcholland.com.



# Table 1. SALSA MLPA Probemix P147-C1 1p36

Longth (nt)	SALSA MI DA probo	Chromosomal position (hg18)		Location (hg19) in kh
Length (III)	SALSA MLPA probe	Reference	Target region	
64-105	Control fragments – see table in prob	pemix content section	on for more informa	tion
124 *	Reference probe 15370-L13762	7q11		
130 «	TNFRSF4 probe 02269-L01761		1p36.33	01-001,137
136 *	CHD5 probe 09115-L09175		1p36.31	01-006,088
142 *	CCNL2 probe 23079-L32576		1p36.33	01-001,323
148 «	GABRD probe 13354-L14784		1p36.33	01-001,949
154 *	Reference probe 13560-L15872	19p13	-	
160 «	GABRD probe 04690-L04068		1p36.33	01-001,946
166	TNFRSF1B probe 00553-L00122		1p36.22	01-012,171
178	GNB1 probe 02890-L02511		1p36.33	01-001,747
184	CASP9 probe 02880-L02347		1p36.21	01-015,704
191 *	Reference probe 06057-L05512	4p16		
196 *	EPHA8 probe 23080-L32577	i	1p36.12	01-022.794
202	SCNN1D probe 04692-L04070		1p36.33	01-001.207
208 *	<b>RER1 probe</b> 23081-L32578		1p36.32	01-002,318
214	AJAP1 probe 04704-L04082		1p36.32	01-004.672
220 *	Reference probe 16368-L18761	12a13		
229 «	<b>AGRN probe</b> 04687-1 04065		1p36.33	01-000.948
238	ICMT probe 04698-1 04076		1p36.31	01-006,208
247 «	PANK4 probe 01122-1 00680		1p36.32	01-002 443
255 * «	<b>TNERSE18 probe</b> 23082-1 32579		1p36.33	01-001,130
263	<b>GNB1 probe</b> 13358-I 15767		1p36.33	01-001 766
270 *	Reference probe 17398-L32045	3n21		01 001,700
276	MTOR probe 04679-1 15769	0021	1n36.22	01-011 091
283	<b>NPHP4 probe</b> 04700-1 04078		1p36.31	01-005 969
292 *	ACTRT2 probe 23083-132580		1p36.32	01-002 929
304 ¥	<b>KIF1B probe</b> 23090-1 04451		1p36.22	01-010 358
313 *	Reference probe 06580-1 24038	2a24		
319	<b>DEEB probe</b> 04696-1 04074	-9	1n36.32	01-003 790
327 «	<b>TP73 probe</b> 01682-1 01262		1n36.32	01-003 558
337	PLOD1 probe 04685-1 04063		1p36.22	01-011.943
346	<b>TNERSE14 probe</b> 04693-1 04071		1p36.32	01-002.480
355 «	<b>PRDM2 probe</b> 04702-1 04080		1p36.21	01-013.904
364	<b>CAMTA1 probe</b> 04695-1 04073		1p36.23	01-007.728
373 *	Reference probe 20867-1 28885	21a21		0.007,720
382 «	PRDM2 probe 03423-1 04457		1p36.21	01-013.899
391	<b>SLC45A1 probe</b> 04697-L04075		1p36.23	01-008.327
399 «	<b>TNFRSF4 probe</b> 09198-L15771		1p36.33	01-001.138
410	<b>TNFRSF9 probe</b> 02185-L15770		1p36.23	01-007.923
418 *	Reference probe 09793-L25209	15a15		
427	<b>MTHFR probe</b> 04684-L04062		1p36.22	01-011.773
436 «	<b>SKI probe</b> 02891-L02359		1p36.33	01-002.227
445 *	Reference probe 16286-L18578	13a14		,
454	<b>PARK7 probe</b> 02189-L02365		1p36.23	01-007.968
463	<b>KIF1B probe</b> 04680-1 04059		1p36.22	01-010.215
469 *	Reference probe 20128-1 27639	18n11		
481 ¥	<b>PRDM16 probe</b> 04703-1 25622		1n36.32	01-003 151
490 * «	ISG15 probe 23084-1 32581		1n36.33	01-000 939
500 ×	Reference probe 17001-1 22047	20a11	1,000.00	01000,000
000		20411		

\* New in version C1.

¥ Changed in version C1. Minor alteration, no change in sequence detected.

# Table 2. P147-C1 probes arranged according to chromosomal location

Length	SALSA MLPA	Cons /Even	Partial sequence <sup>a</sup>	Distance to	Distance to
(nt)	probe	Gene/Exon	(24 nt adjacent to ligation site)	p-telomere (hg18)	next probe
490 «	23084-L32581	ISG15	CAAGGGCCGCAG-CAGCACCTACGA	939 kb	8 kb
229 «	04687-L04065	AGRN	CCCCTCATCTGT-GACAACCAGGTG	947 kb	182 kb
255 «	23082-L32579	TNFRSF18	CACAGTCGATAC-ACTGGAAGCCAA	1.129 kb	8 kb
130 «	02269-L01761	<i>TNFRSF4</i> , exon 5 NM_003327.4	GCCGGCCAGCAA-TAGCTCGGACGC	1.137 kb	1 kb
399 «	09198-L15771	<i>TNFRSF4</i> , exon 3 NM_003327.4	AACTCCAGGCTT-GTAGCTGTCCAG	1.138 kb	69 kb
202	04692-L04070	SCNN1D	AGTGACGAAGCT-GTGATTCACACA	1.207 kb	116 kb
142	23079-L32576	CCNL2	CAAGACGCATAC-GGGACGTCATCA	1.323 kb	423 kb
178	02890-L02511	<i>GNB1,</i> exon 3 NM_002074.5	CTAAGATCGGAA-GATGAGTGAGCT	1.746 kb	19 kb
263	13358-L15767	<i>GNB1</i> , intron 1 NM_002074.5	AGTGGTGCACTT-ATGTGTTTCCCA	1.765 kb	181 kb
160 «	04690-L04068	GABRD, exon 2 NM_000815.5	CGGCGACTACGT-GGGCTCCAACCT	1.946 kb	3 kb
148 «	13354-L14784	GABRD, exon 6 NM_000815.5	TCATCGGAGGAC-ATCGTCTACTAC	1.949 kb	278 kb
436 «	02891-L02359	SKI	AACGAGAAGAAG-ATGAAAGAGGCC	2.227 kb	91 kb
208	23081-L32578	RER1	TACATGATTCGA-GTTTACCTGCTG	2.318 kb	125 kb
247 «	01122-L00680	PANK4	CTATTCAACGGT-ACAGCACAAAGT	2.443 kb	36 kb
346	04693-L04071	TNFRSF14	CAATACCCTCAT-TCACGGGGAGGA	2.479 kb	450 kb
292	23083-L32580	ACTRT2	CGGAGGTCCCAA-ACTCCTTGAAGT	2.929 kb	221 kb
481	04703-L25622	PRDM16	ACGGACGTGGAA-GTGTCGCCCCAG	3.150 kb	408 kb
327 «	01682-L01262	TP73	GAGACCCGGGTG-TCAGGAAAGATG	3.558 kb	232 kb
319	04696-L04074	DFFB	TGCACATTGTCT-GCCATAAGAAAA	3.790 kb	0,9 <b>M</b> b
214	04704-L04082	AJAP1	TGATAGCCATGT-TTCAGCTCGCCG	4.671 kb	1,3 <b>M</b> b
283	04700-L04078	NPHP4	GGATGAACGACT-GGCACAGGATCT	5.968 kb	119 kb
136	09115-L09175	CHD5	TGCAGCACTGAT-GTCTCTTTACCG	6.087 kb	120 kb
238	04698-L04076	ICMT	CAGCTATGCCCT-GACAGTGTGGCG	6.207 kb	1,5 <b>M</b> b
364	04695-L04073	CAMTA1	AATGAGCTGGCT-GGCCAGTTATCT	7.727 kb	195 kb
410	02185-L15770	TNFRSF9	GGTCCTCAACTT-TGAGAGGACAAG	7.922 kb	45 kb
454	02189-L02365	PARK7	AGAGCAGCGAAC-TGCGACGATCAC	7.967 kb	359 kb
391	04697-L04075	SLC45A1	CAACGGGGTGAT-GTACTTCTCCAG	8.326 kb	1,9 <b>M</b> b
463	04680-L04059	<i>KIF1B,</i> exon 2 NM_001365951.3	CTCAGTGAAGGT-GGCTGTCCGGGT	10.214 kb	143 kb
304	23090-L04451	<i>KIF1B,</i> exon 48 (46) NM_001365951.3	CGTGGGGTCCTT-TTGCAGGCCCTC	10.357 kb	733 kb
276	04679-L15769	MTOR	TCCAACGCAAGT-TGAGCTGCTCAT	11.090 kb	683 kb
427	04684-L04062	MTHFR	TGACTTCCCACT-GGACAACTGCCT	11.773 kb	170 kb
337	04685-L04063	PLOD1	GCATGGCAGCGA-GTACCAGTCTGT	11.943 kb	228 kb
166	00553-L00122	TNFRSF1B	GGCTCAGAGAAT-ACTATGACCAGA	12.171 kb	1,7 <b>M</b> b
382 «	03423-L04457	<i>PRDM2</i> , upstream NM_012231.5	GCCATTGGGCGA-CGGCGCAGGGTC	13.899 kb	5 kb
355 «	04702-L04080	<i>PRDM2</i> , exon 1 NM_012231.5	TTGACCTTCCCT-CCACTCTTACAG	13.904 kb	1,8 <b>M</b> b
184	02880-L02347	CASP9	GGTCGAGAAGAT-TGTGAACATCTT	15.703 kb	7,1 <b>M</b> b
196	23080-L32577	EPHA8	GCCTGACGCTCA-TCACGGGCCTGG	22.794 kb	

<sup>a</sup> Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

The exon numbering used in previous versions of this product description can be found in between brackets.



SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

# Related products

For related products, see the product page on our website.

# References

- Gajecka M et al. (2007). Monosomy 1p36 deletion syndrome. Am J Med Genet Part C. 145C:346-356.
- Jordan VK et al. (2015). 1p36 deletion syndrome: an update. Appl Clin Genet. 8:189-200.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Slavotinek A et al. (1999). Monosomy 1p36. J Med Genet. 36:657-663.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

# Selected publications using SALSA MLPA Probemix P147 1p36

- D'Angelo CS et al. (2006) Prader-Willi-like phenotype: investigation of 1p36 deletion in 41 patients with delayed psychomotor development, hypotonia, obesity and/or hyperphagia, learning disabilities and behavioral problems. *Eur J Med Genet*. 49(6):451-60.
- D'Angelo CS et al. (2009). Extending the phenotype of monosomy 1p36 syndrome and mapping of a critical region for obesity and hyperphagia. *Am J Med Genet Part A*. 152A:102-110.
- Hirschfeldova K et al. (2011). Cryptic chromosomal rearrangements in children with idiopathic mental retardation in the Czech population. *Genetic Testing and Molecular Biomarkers*. 15:607-611.
- Xu F et al. (2014) The first patient with a pure 1p36 microtriplication associated with severe clinical phenotypes. *Mol Cytogenet*. 7(1):64.

P147 product history		
Version	Modification	
C1	Seven target probes and nine reference probes have been replaced and one reference probe has been added. Two target probes changed in length, but not in sequence detected.	
B2	Four reference probes were replaced and the control fragments were adjusted (QDX2).	
B1	Four probes on 1p36 were removed and five new probes were included. In addition, four extra control fragments at 88, 96, 100 and 105 nt were added.	
A1	First release.	

## Implemented changes in the product description

Version C1-02- 28 April 2025 (04P)

- Gene structure and transcript variants section: link to MANE website added.
- NM sequence of the KIF1B probes updated.
- Exon numbering of the 304 nt probe changed.
- Removed Related SALSA MLPA products section.
- Small layout changes.

Version C1-01- 09 August 2021 (04P)

- Product description rewritten and adapted to a new template.

- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).



More information: www.mrcholland.com; www.mrcholland.eu		
	MRC Holland bv; Willem Schoutenstraat 1	
	1057 DL, Amsterdam, The Netherlands	
E-mail	info@mrcholland.com (information & technical questions)	
	order@mrcholland.com (orders)	
Phone	+31 888 657 200	