

## Instructions for Use SALSA<sup>®</sup> MLPA<sup>®</sup> Probemix P190 CHEK2

L See also the MLPA General Protocol, the product descriptions of the SALSA<sup>®</sup> MLPA<sup>®</sup> Reagent Kit and SALSA<sup>®</sup> Binning DNA SD078, and the Coffalyser.Net Reference Manual.

Visit the SALSA® MLPA® Probemix P190 CHEK2 product page on our website to find Certificates of Analysis and a list of related products.

| Broduct Namo         | SALSA <sup>®</sup> MLPA <sup>®</sup> Probemix |
|----------------------|---|
| FIGUELINAME          | P190 CHEK2                                    |
| Version              | D1  |
| Catalogue<br>numbers | P190-025R (25 reactions)                      |
|                      | P190-050R (50 reactions)                      |
|                      | P190-100R (100 reactions)                     |
| Basic UDI-DI         | 872021148P19064                               |
|                      | Synthetic oligonucleotides,                   |
| Ingredients          | oligonucleotides purified from bacteria,      |
|                      | Tris-HCl, EDTA                                |

| Additional Test Components                       | Catalogue<br>Numbers |
|--|----------------------|
|  | EK1-FAM              |
|  | EK1-CY5              |
| SALSA <sup>®</sup> MLPA <sup>®</sup> Reagent Kit | EK5-FAM              |
|  | EK5-CY5              |
|  | EK20-FAM             |
| SALSA <sup>®</sup> Binning DNA SD078             | SD078                |

#### Storage and Shelf Life

| Recommended conditions | -25°C | ** |
|------------------------|-------|----|
|------------------------|-------|----|

A shelf life of until the expiry date is guaranteed, also after opening when stored in the original packaging under recommended conditions. For the exact expiry date, see the label on the vial. This product should not be exposed to more than 25 freeze-thaw cycles. Do not use the product if the packaging is damaged or opened. Leave chemicals in original containers. Waste material must be disposed of in accordance with the national and local regulations. 
 Regulatory Status

 IVD
 EUROPE
 2797

 ISRAEL
 ALL OTHER COUNTRIES

| Label Symbols    |  |                  |                           |
|------------------|--|------------------|---------------------------|
| IVD              | In Vitro Diagnostic  | RUO              | Research Use Only         |
| More in<br>www.m | formation:<br>archolland.com   |                  |                           |
|                  | MRC Holland BV; Willer<br>1057 DL, Amsterdam, t  | n Scho<br>he Net | outenstraat 1<br>herlands |
| E-mail           | info@mrcholland.com (information & technical<br>questions);<br>order@mrcholland.com (orders) |                  |                           |
| Phone            | +31 888 657 200  |                  | •                         |

Any serious incident that has occurred in relation to this product should be reported to MRC Holland and the competent authority of the Member State in which the user and/or the patient is located.

### **Changes in this Product Version**

Compared to C1 version

16 target probes have been added and several target probes have been adjusted, replaced or removed. Almost all reference probes have been replaced and two were added. Several probes have a small change in length but no change in sequence detected.

### 1. Intended Purpose

The SALSA MLPA Probemix P190 CHEK2 is an in vitro diagnostic  $(IVD)^1$  or research use only (RUO) semi-quantitative manual assay<sup>2</sup> for the detection of deletions or duplications in the human *CHEK2*, and the presence of the *CHEK2* c.1100delC mutation in genomic DNA isolated from human peripheral whole blood specimens. P190 CHEK2 is intended to confirm a potential cause for and clinical diagnosis of a predisposition to breast cancer and other *CHEK2*-related cancer types. In addition, deletions and duplications in the human *ATM* and *TP53* genes can be detected with this probemix for differential diagnosis of *ATM*- or *TP53*-related predisposition to cancer, respectively, in patients originally suspected of a *CHEK2*-related cancer predisposition. This product can also be used to determine increased cancer susceptibility in at-risk family members.

Copy number variations (CNVs) detected with P190 CHEK2 should be confirmed with a different technique. In particular, CNVs detected by only a single probe always require confirmation by another method. Most defects in the *CHEK2*, *ATM* and *TP53* genes are point mutations, of which only the *CHEK2* c.1100delC mutation will be detected by MLPA. It is therefore recommended to use this assay in combination with sequence analysis.

Please note that this probemix covers all exons of *CHEK2* but not of *ATM* and *TP53*. For the latter two genes, the SALSA MLPA Probemixes P041/P042 ATM and SALSA MLPA Probemix P056 TP53 provide a better coverage and may detect aberrations that are not detected by this P190 CHEK2 probemix.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes, pre-implantation or prenatal testing, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations, e.g. from DNA extracted from formalin-fixed paraffin embedded (FFPE) or fresh tumour materials.

<sup>1</sup> Please note that this probemix is for IVD use in the countries specified on page 1 of this product description. In all other countries, this is a RUO product.

<sup>2</sup> To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

### 2. Sample Requirements

| Specimen             | 50-250 ng purified human genomic DNA, dissolved in 5 $\mu$ l TE <sub>0.1</sub> buffer, pH 8.0-8.5   |
|----------------------|---|
| Collection<br>Method | standard methods  |
| Extraction<br>Method | <ul> <li>Methods tested by MRC Holland:</li> <li>QIAGEN Autopure LS (automated) and<br/>QIAamp DNA mini/midi/maxi kit (manual)</li> <li>Promega Wizard Genomic DNA Purification<br/>Kit (manual)</li> <li>salting out (manual)</li> </ul> |

| Sample Types                                   |  |   |  |
|--|--|---|--|
| Test Sample                                    | <ul> <li>provided by user</li> </ul>   |   |  |
| Reference<br>Samples<br>(Required)             | <ul> <li>provided by user</li> <li>extraction method, tissue type, DNA concentration and treatment as similar as possible in all test and reference samples.</li> <li>have a normal copy number and ≤0.10 standard deviation for all probes except for the mutation-specific probes.</li> <li>at least three* independent reference samples required in each experiment for proper data normalisation. Derived from unrelated individuals from families without a history of cancer predisposition.</li> </ul> |   |  |
| No-DNA<br>Control<br>(Preferably)              | <ul> <li>provided by user</li> <li>TE<sub>0.1</sub> buffer instead of DNA</li> <li>to check for DNA contamination</li> </ul>   |   |  |
| Binning<br>Sample<br>(Initial<br>Experiment)   | <ul> <li>SALSA Binning DNA SD078, provided by<br/>MRC Holland</li> <li>required in initial experiment to determine<br/>suitable bin set</li> <li>should never be used as a reference sample</li> </ul>   |   |  |
|  | <ul> <li>provided by user,</li> </ul>  | or  |  |
| Positive<br>Control<br>Samples<br>(Preferably) | Available from<br>third parties  | See the table of positive<br>samples on the<br>probemix product page<br>on our website. |  |
| Validation<br>Samples<br>(Required)            | <ul> <li>In the validation experiments of this<br/>probemix, the peaks of the mutation-specific<br/>probes are expected to be absent in the<br/>majority of samples from healthy<br/>individuals.</li> </ul>   |   |  |

\*When testing >21 samples, include one extra reference for each 7 test samples.





### 3. Test Procedure

See the <u>MLPA General Protocol</u>.

### 4. Quality Control, Data Analysis, and Troubleshooting

| Quality Control Fragments in the Probemix |  |
|---|--|
| Length (nt)                               | Function                               |
| 64-70-76-82                               | DNA quantity control fragments         |
| 88-96                                     | DNA denaturation control fragments     |
| 92  | Benchmark fragment                     |
| 100                                       | Chromosome X presence control fragment |
| 105                                       | Chromosome Y presence control fragment |

<u>Coffalyser.Net</u> should be used for data analysis in combination with the appropriate product and lot-specific Coffalyser sheet. See the <u>Coffalyser.Net Reference Manual</u> for details on data analysis and quality control.

For troubleshooting help, see the additional resources offered on our support portal.

### 5. Interpretation of Results

# Determining Typical Values in Normal and Affected Populations

The typical final ratio (FR) values stated in the copy number tables were determined in a validation study with samples containing abnormal copy numbers. The standard deviation of each individual probe over all the reference samples was  $\leq 0.10$ .

Expected Results of Reference Probes

| Final Ratio<br>(FR) | Copy<br>Number | Description |
|---------------------|----------------|-------------|
| 0.80 - 1.20         | 2              | Normal      |

### Typical Results of Probes Targeting Two Copies

| Final Ratio<br>(FR) | Copy<br>Number | Description   |
|---------------------|----------------|---|
| 0                   | 0              | Homozygous deletion                                       |
| 0.40 - 0.65         | 1              | Heterozygous deletion                                     |
| 0.80 - 1.20         | 2              | Normal  |
| 1.30 - 1.65         | 3              | Heterozygous duplication                                  |
| 1.75 – 2.15         | 4              | Homozygous duplication<br>or<br>Heterozygous triplication |
| All other<br>values | -              | Ambiguous   |

The tables illustrate the relationship between final probe ratio and corresponding copy number. Test results are expected to center around these values. Ambiguous values can indicate a technical problem, but may also reflect a biological cause such as mosaicism or a SNV influencing a single probe. It is important to use Coffalyser.Net to determine the significance of values found.

#### Possible Results of Mutation-Specific Probes

| Signal Strength     | Mutation Status                      |
|---------------------|--------------------------------------|
| ≥10% median         |                                      |
| peak neight         | Mutation CHEK2 C. I ToudelC detected |
| probes              | (expected only in positive samples)  |
| <10% median         | Mutation CHEK2 c 1100delC not        |
| peak height         | detected (expected in most samples   |
| reference<br>probes | from healthy individuals)            |



### 6. Performance Characteristics

| Study  | Description  |   |                          |   |
|--|--|---|--------------------------|---|
| Expected values for copy<br>numbers in normal and<br>affected populations                | To determine the expected values in normal and affected populations a study was conducted on over 1500 MLPA reactions using samples with and without abnormal copy numbers. When the standard deviation of each individual probe over all the reference samples is ≤0.10, the ranges stated in the copy number table above can be used. Cut-off values for copy number determination were verified with SALSA MLPA Probemix P190 CHEK2 in 42 samples from healthy individuals with normal copy number and seven samples with known CNVs. The expected FRs for the corresponding copy number were found in all samples tested.  |   |                          |   |
| Expected values for point<br>mutation detection in<br>normal and affected<br>populations | The mutation-specific probe will only generate a signal when the <i>CHEK2</i> c.1100delC mutation is present.<br>Please note that background signals of the mutation-specific probe can be expected above the threshold<br>in some cases. Users should always compare the relative peak height of the mutation-specific probe in<br>mutation-positive samples to the relative peak height in reference samples. A clear signal (at least 10%<br>of the median peak height of all reference probes in that sample) indicates that the mutation is present.<br>It is not possible to determine the copy number of mutation-specific probes.<br>The expected values for the mutation-specific probe were verified with SALSA MLPA Probemix P190<br>CHEK2 using one mutation positive sample, six samples positive for other aberrations detected by P190<br>CHEK2, and 42 samples from healthy individuals without the <i>CHEK2</i> c.1100delC mutation, and the<br>expected signal percentage was found in all tested samples. |   |                          |   |
| Limit of detection   | A study using representative probemixes was conducted to evaluate the minimum and maximum amount of DNA acceptable as the assay input. Results support the use of 50-250 ng of human DNA as the recommend input amount. The use of insufficient or too much sample DNA can affect performance. These lower and higher limits of detection were verified using SALSA MLPA Probemix P190 CHEK2 on five samples with known CNVs/mutations and on one sample without any mutation. Expected results were obtained in all samples using both the lower and upper input amount of DNA.   |   |                          |   |
| Interfering substances   | SNVs or other pol  | ymorphisms (e.g. indel<br>FDTA and hemoglobin | s) in the DNA targ       | et sequence and impurities in the DNA sample  |
|  | A study using SALSA MLPA Probemix P190 CHEK2 was performed to assess the potential for interference of endogenous and exogenous substances on genomic DNA on six samples with known CNVs/mutations status. For most probes, expected FRs/signals (FRs within the expected cut-off category) were obtained even in the presence of potential interferents at concentrations shown in the table below.   |   |                          |   |
|  | Interferent  | Source  | Testing<br>Concentration | Results*  |
|  | EDTA   | Exogenous –<br>specimen collection<br>tubes   | 1.5 mM                   | <u>Copy number</u> : Expected FR for 730/738<br>measurements<br><u>Mutation</u> : Expected signal for 18/18<br>measurements                 |
|  | NaCl   | Exogenous – DNA<br>extraction                 | 40 mM                    | <u>Copy number</u> : Expected FR for 731/738<br>measurements<br><u>Mutation</u> : Expected signal for 18/18<br>measurements                 |
|  | Fe³⁺ (FeCl₃)   | Exogenous – DNA<br>extraction                 | 1 µM                     | <u>Copy number</u> : Expected FR for 736/738<br>measurements<br><u>Mutation</u> : Expected signal for 18/18<br>measurements                 |
|  | Heparin  | Exogenous –<br>specimen collection<br>tubes   | 0.02 U/mL                | <u>Copy number</u> : Expected FR for 735/738<br>measurements<br><u>Mutation</u> : Expected signal for 18/18<br>measurements                 |
|  | Hemoglobin   | Endogenous –<br>blood sample                  | 0.02 µg/µl               | <u>Copy number</u> : Expected FR for 661/738<br>measurements<br><u>Mutation</u> : Expected signal for 18/18<br>measurements                 |
|  | <ul> <li>* Results are summarised for all probes across all four samples tested in triplicate.</li> <li>NaCl, FeCl<sub>3</sub> and heparin did not interfere with copy number determination, while an effect on the FRs was observed for a low number of probes with EDTA. Hemoglobin had the largest effect on the FRs ambiguous as well as false results were obtained. None of the interferents had an effect on the determination of mutation status.</li> </ul>   |   |                          | mples tested in triplicate.   |
|  |  |   |                          | determination, while an effect on the FRs was<br>noglobin had the largest effect on the FRs;<br>ne of the interferents had an effect on the |
|  |  |   |                          |   |



| Cross-reactivity   | Cross-reactivity is the potential for probes to bind to homologous regions (e.g. pseudogenes) or other cross-reactive sequences. Seven probes for <i>CHEK2</i> target sequences for which homologous regions were identified. Quality tests were carried out to determine whether probes are specific to their target sequence and all probes met the quality criteria for specificity.   |
|--------------------|---|
| Accuracy           | Results of accuracy are derived from trueness and precision studies. For trueness, eight previously genotyped samples were tested using SALSA MLPA Probemix P190 CHEK2 and found to have the expected results. Assay precision was tested by repeatedly testing samples with known copy number/mutations over multiple days, and by multiple operators. Results showed a correct call in 100% between triplicates for both copy number determination and determination of mutation status. Reproducibility between operators and days is 98.6% for copy number determination and 100% for determination of the mutation status. Overall, precision is >99%. |
| Clinical validity* | <i>CHEK2</i> : Deletions/duplications in the <i>CHEK2</i> gene are detected in 0-3.5% of individuals with a predisposition to <i>CHEK2</i> -related cancer.   |
|                    | The <i>CHEK2</i> c.1100delC mutation is the most common founder mutation of this gene and it is of Northern and Eastern European origin. Therefore, the frequency of <i>CHEK2</i> c.1100delC mutation varies across populations.  |
|                    | <i>ATM</i> : Deletions/duplications in the <i>ATM</i> gene are detected in <1% of individuals with a predisposition to <i>ATM</i> -related cancer.  |
|                    | <i>TP53</i> : Deletions/duplications in <i>TP53</i> gene are detected in <1% of individuals with a predisposition to <i>TP53</i> -related cancer.   |
|                    | *Based on a 2008-2024 literature review   |

Summary of Safety and Performance (SSP) The SSP is available in the European database on medical devices (Eudamed), <u>https://ec.europa.eu/tools/eudamed</u>, or upon request.



#### Content - Probe Details Sorted by Chromosomal Position

| Chr.               | Target    | Exon                                  | Distance to       | Mutation   | Length (nt) | Probe number         | Warnings     |
|--------------------|-----------|---------------------------------------|-------------------|------------|-------------|----------------------|--------------|
| position           | A TA 4    | Even 1                                | 10.0 kb           |            | 200         | 21665120207          |              |
| 11q22.3            | ATM       | Exon F                                | 12.8 KD           |            | 328         | 21005-L30287         | +            |
| 11q22.3            | ATM       | EXON 9                                | 11.2 KD           |            | 430         | 21008-L30290         |              |
| 11q22.3            | ATM       | Exon 0                                | 2.1 KU<br>9.5 kb  |            | 292         | 21003-L30525         |              |
| 11q22.3            | ATM       | EXON 9                                | 0.0 KD            |            | 238         | 21000-L30282         |              |
| 11q22.3            | ATM       | Exon 13                               | 9.8 KD            |            | 483         | 210/2-L30294         |              |
| 11q22.3            | ATM       | Exon 20                               | 4.0 KD            |            | 280         | 21650120284          |              |
| 11q22.3            |           | Exon 24                               | 9.0 KD<br>9.5 kb  |            | 106         | 21657-120270         |              |
| 11q22.3            |           | Exon 20                               | 0.3 KD<br>7.7 kb  |            | 190         | 21637-230279         |              |
| 11q22.3            |           | Exon 22                               | 10.6 kb           |            | 210         | 21071-230293         |              |
| 11q22.3            |           | Exon 29                               | 0.5 kb            |            | 160         | 21652-120274         |              |
| 11q22.3            |           | Exon 42                               | 9.J KU<br>2.9 kb  |            | 100         | 21651-120272         |              |
| 11q22.3            |           | Exon 45                               | 3.0 KD<br>7.7 kb  |            | 130         | 21031-230273         |              |
| 11q22.3            |           | Exon 40                               | 7.7 KD<br>6.0 kb  |            | 262         | 21661-120291         |              |
| 11q22.3            |           | Exon 55                               | 0.0 KD            |            | 19/         | 21655-120277         |              |
| 11q22.3            |           | Exon 59                               | 0.0 kb            |            | 104         | 21650-120272         |              |
| 11q22.3            |           | Exon 61                               | 9.0 KD            |            | 202         | 21658-130272         |              |
| 11q22.3            |           | Exon 62                               | 0.3 KD            |            | 202         | 21636-L30280         |              |
| 11q22.3            |           | Exon 63                               | 0.3 KD            |            | 100         | 21656-1 30840        |              |
| 17q22.3            | TD52      | Exon 11                               | 30 kh             |            | 346         | 00345-131314         |              |
| 17p13.1<br>17p13.1 | TP53      | Exon 9 (8)                            | 2.9 KD            |            | 340<br>166  | 21652-130275         |              |
| 17p13.1            | TD53      | Exon $f(0)$                           | 2.7 KD<br>11.2 kb |            | 226         | 01007-1 31312        |              |
| 17p13.1            | TD53      | Exon 1                                | 11.2 KU           |            | 220         | 21581-1 25082        |              |
| 22a12.1            | CHEK2     | Exon 15                               | 10 kb             |            | 201         | 07281-SP0800-1 30/57 | <i>"</i> ¥ # |
| 22q12.1            | CHEK2     | Intron 14 (Evon 14)                   | 0.4 kb            |            | 230         | 06636-1 30524        | «#Ø          |
| 22q12.1            | CHEK2     | Evon 1/                               | 4.7 kb            |            | 382         | 21666-1 30288        | «#±          |
| 22q12.1<br>22q12.1 | CHEK2     | Exon 13                               | 4.7 KD            |            | 178         | 21654-1 30276        | <u>«#П</u>   |
| 22q12.1            | CHEK2     | Exon 13                               | 0.1 kb            |            | 355         | 19654-1 26320        | «#+          |
| 22q12.1            | CHEK2     | Exon 12                               | 0.6 kb            |            | 143         | 21419-1 29916        | «#+          |
| 22g12.1            | CHEK2     | Exon 11                               | 0.2 kb            |            | 172         | 21418-1 29915        | «#+          |
| 22g12.1            | CHEK2     | Exon 11                               | 1.1 kb            | c.1100delC | 313         | 22034-SP0468-L31261  | 8 Ж » П      |
| 22g12.1            | CHEK2     | Exon 10                               | 2.9 kb            |            | 418         | 06632-L06190         | <u></u> «П   |
| 22g12.1            | CHEK2     | Exon 9                                | 3.6 kb            |            | 400         | 02579-L02041         | «[»          |
| 22g12.1            | CHEK2     | Exon 8                                | 6.5 kb            |            | 427         | 06631-L06189         | «П           |
| 22g12.1            | CHEK2     | Exon 7                                | 2.0 kb            |            | 154         | 06630-L07119         |              |
| 22g12.1            | CHEK2     | Exon 6                                | 7.5 kb            |            | 208         | 06629-L30453         |              |
| 22g12.1            | CHEK2     | Exon 5                                | 5.6 kb            |            | 373         | 06628-L06186         |              |
| 22g12.1            | CHEK2     | Exon 4                                | 0.2 kb            |            | 274         | 06627-L06185         |              |
| 22q12.1            | CHEK2     | Exon 3                                | 9.2 kb            |            | 306         | 06626-L30841         |              |
| 22q12.1            | CHEK2     | Exon 2                                | 7.3 kb            |            | 391         | 06625-L06183         |              |
| 22q12.1            | CHEK2     | Exon 1                                | 0.2 kb            |            | 337         | 06624-L24131         | »            |
| 22q12.1            | CHEK2     | Upstream (Exon 1)                     | 1.9 kb            |            | 220         | 06623-L31306         | Ø            |
| 22q12.1            | HSCB      | · · · · · · · · · · · · · · · · · · · |                   |            | 268         | 06800-L30458         | ~            |
| 2p                 | Reference |                                       |                   |            | 250         | 17871-L22467         |              |
| 2q                 | Reference |                                       |                   |            | 148         | 14199-L23450         |              |
| 3p                 | Reference |                                       |                   |            | 298         | 15388-L17790         |              |
| 4p                 | Reference |                                       |                   |            | 490         | 20096-L27538         |              |
| 5p                 | Reference |                                       |                   |            | 454         | 18691-L02476         |              |
| 5q                 | Reference |                                       |                   |            | 124         | 18709-L21056         |              |
| 6р                 | Reference |                                       |                   |            | 214         | 10730-L30523         |              |
| 9q                 | Reference |                                       |                   |            | 409         | 08725-L08736         |              |
| 11p                | Reference |                                       |                   |            | 364         | 18676-L24030         |              |
| 20q                | Reference |                                       |                   |            | 500         | 17001-L22947         |              |

Probe lengths may vary slightly depending on capillary electrophoresis instrument settings. Please see the most up-to-date Coffalyser sheet for exact probe lengths obtained at MRC Holland.

The *CHEK2*, *ATM*, and *TP53* exon numbers are derived from MANE project and are based on MANE Select transcript. For more information, see the probe sequences document available on the product page at <u>www.mrcholland.com</u>. Annotations of several probes with a target at the edge of or slightly outside the coding region are altered. The exon numbering from the previous version of this product description is disclosed between brackets.

Chromosomal bands are based on: hg18.



### 7. Precautions and Warnings

#### Probe warnings

- S This probe will only generate a signal when the CHEK2 c.1100delC mutation is present. It has been tested on artificial DNA and on positive human samples. However, the probe can give an extra signal due to a simultaneous activity of the ligase and polymerase enzymes.
- This probe is a flanking probe, included to help determine the extent of a deletion/duplication. Copy number alterations of flanking probes are unlikely to be related to the condition tested.
- This probe is located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation.
- X This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.
- # The specificity of this probe relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.
- » Detects the same sequence as the CHEK2 probes in SALSA® MLPA® Probemix P045 BRCA2/CHEK2.
- Π Detects the same sequence as the CHEK2 probes in SALSA<sup>®</sup> MLPA<sup>®</sup> Probemix P056 TP53.
- Probe insensitive to depurination. A high signal of the 400 nt probe can be due to depurination of the sample DNA, e.g. due to insufficient buffer concentration in the DNA sample or a prolonged denaturation time. Reduced signals of other probes caused by sample depurination lead to seemingly high signals of the 400 nt probe.
- + The ligation site of this probe is >20 nt away from the nearest exon. For more information, download the probe sequence sheet from the probemix-specific page on www.mrcholland.com.
- Ø This probe targets a sequence outside of the known coding region. Copy number alterations of only this probe are of unknown clinical significance.

Probemix-specific precautions

- This product is not known to contain any harmful agents. Based on the concentrations present, none of the ingredients are hazardous as defined by the Hazard Communication Standard. A Safety Data Sheet (SDS) is not required for this product: none of the ingredients contain dangerous substances at concentrations requiring distribution of an SDS (as per Regulation (EC) No 1272/2008 [EU-GHS/CLP] and 1907/2006 [REACH] and amendments).
- 2. Sample or technical artefacts may appear as a (mosaic) copy number change of the whole/partial gene. Whole/partial gene deletions or duplications should therefore be confirmed by analysis of an independent DNA sample, to exclude false positive results.
- 3. Small changes (e.g. SNVs, small indels) in the sequence targeted by a probe can cause false positive results, even when >20 nt from the probe ligation site. Sequence changes can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA. Deviations detected by this product should be confirmed, and single-probe deviations always require confirmation. Sequencing of the target region is recommended. Please contact MRC Holland for more information: info@mrcholland.com.

- 4. CHEK2 c.1100delC mutation-specific probe: We have received reports of experiments in which a peak for the CHEK2 c.1100delC probe appeared in *all* samples. This was due to simultaneous ligase and polymerase activity caused by either incomplete heat inactivation of Ligase-65 or contamination of ligase mastermix with polymerase mastermix or vice versa. For more information on this issue, please contact info@mrcholland.com. Results obtained with this CHEK2 c.1100delC probe should therefore be treated with caution.
- 5. Deletions of the <u>last exons of ATM</u> (exon 62-63) are encountered relatively frequently (own validation observations, Micol et al. 2011, Nakamura et al. 2012, Podralska et al. 2014, Susswein et al. 2016, Tung et al. 2015). *Duplication* of these exons might not result in inactivation of that gene copy and should therefore be interpreted with caution. In the ClinVar database (https://www.ncbi.nlm.nih.gov/clinvar/), these variants are classified as variant of uncertain significance (VUS)/likely benign.
- 6. Copy number alterations of reference probes are unlikely to be related to the condition tested.

### Technique-specific precautions

See the <u>MLPA General Protocol</u>.

### 8. Limitations

Probemix-specific limitations

1. The *CHEK2* c.1100delC mutation-specific probe is only intended to determine the presence (or absence) of the mutation.

Technique-specific limitations

See the MLPA General Protocol.

### 9. References Cited in this IFU

- 1. Micol R et al. (2011). Morbidity and mortality from Ataxia-Telangiectasia are associated with ATM genotype. J Allergy Clin Immunol. 128:382-9. E381.
- Nakamura K et al. (2012). Functional characterization and targeted correction of ATM mutations identified in Japanese patients with Ataxia-Telangiectasia. Hum Mutat. 33:198-208.
- Podralska MJ et al. (2014). Ten new ATM alterations in Polish patients with Ataxia-Telangiectasia. Mol Genet Genomic Med. 2:504-11.
- Susswein LR et al. (2016). Pathogenic and likely pathogenic variant prevalence among the first 10,000 patients referred for next-generation cancer panel testing. Genet Med. 18:823-32.
- 5. Tung N et al. (2015). Frequency of mutations in individuals with breast cancer referred for BRCA1 and BRCA2 testing using next-generation sequencing with a 25-gene panel. Cancer. 121:25-33.



#### Implemented changes in the product description

Version D1-07 – 21 March 2025 (03S)

- Updated to a new template.
- Intended purpose was updated, specifying the testing population and clinical application of *ATM* and *TP53* copy number determination, and specifying assay is manual.
- Reference to SALSA Binning DNA SD078 removed from the intended purpose footnote.
- Exon numbering updated for *TP53* probes.
- Description of probe targets at the edge of or slightly outside the coding region has been adjusted. No change in actual target sites.
- SNV rs564605612 can affect the probe signal. However, the warning for this SNV present in previous product description versions has been replaced by a general warning for small sequence changes in the section Precautions and warnings.
   Probemix is now IVDR certified.

MRC Holland, SALSA, MLPA, digitalMLPA, Coffalyser.Net, Coffalyser digitalMLPA, and their logos are trademarks or registered trademarks of MRC Holland BV. All other brands and names herein are the property of their respective owners.