

# Product Description

## SALSA® MLPA® Probemix P488-B1 RS1

To be used with the MLPA General Protocol.

### Version B1

For complete product history see page 6.

### Catalogue numbers

- **P488-025R:** SALSA® MLPA® Probemix P488 RS1, 25 reactions
- **P488-050R:** SALSA® MLPA® Probemix P488 RS1, 50 reactions
- **P488-100R:** SALSA® MLPA® Probemix P488 RS1, 100 reactions

SALSA® MLPA® Probemix P488 RS1 (hereafter: P488 RS1) is to be used in combination with:

1. SALSA® MLPA® Reagent Kit (Cat. No: EK1-FAM, EK1-CY5, EK5-FAM, EK5-CY5, EK20-FAM),
2. Data analysis software Coffalyser.Net™ (Cat. No: n.a.)

### Volumes and ingredients

Volumes			Ingredients
P488-025R	P488-050R	P488-100R	
40 µl	80 µl	160 µl	Synthetic oligonucleotides, oligonucleotides purified from bacteria, Tris-HCl, EDTA

The MLPA probemix is not known to contain any harmful agents. Based on the concentrations present, none of the ingredients are hazardous as defined by the Hazard Communication Standard. A Safety Data Sheet (SDS) is not required for this product: none of the ingredients contain dangerous substances at concentrations requiring distribution of an SDS (as per Regulation (EC) No 1272/2008 [EU-GHS/CLP] and 1907/2006 [REACH] and amendments).

### Storage and handling

Recommended storage conditions		
--------------------------------	---	---

A shelf life of until the expiry date is guaranteed, when stored in the original packaging under recommended conditions. For the exact expiry date, see the label on the vial. This product should not be exposed to more than 25 freeze-thaw cycles. Do not use the product if the packaging is damaged or opened. Leave chemicals in original containers. Waste material must be disposed of in accordance with the national and local regulations.

### Certificate of Analysis

Information regarding quality tests and a sample electropherogram from the current sales lot is available at [www.mrcholland.com](http://www.mrcholland.com).

### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: [www.mrcholland.com](http://www.mrcholland.com). It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

### General information

SALSA® MLPA® Probemix P488 RS1 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *RS1* gene, which is associated with X-linked juvenile retinoschisis.

X-linked juvenile retinoschisis is a retinal dystrophy that leads to schisis (splitting) of the neural retina leading to reduced visual acuity in affected men. The split in the retina occurs predominantly within the inner retinal layers and is very different from retinal detachment, which is a split between the neural retina and the retinal

pigment epithelium. This condition is caused by defects in the *RS1* gene that lead to a dysfunctional retinoschisin protein which causes impaired adhesion of retinal cells. The prevalence of X-linked juvenile retinoschisis is estimated between 1 in 15,000 and 1 in 30,000 (Sikkink et al. 2007).

The *RS1* gene (6 exons) spans ~32.4 kb of genomic DNA and is located on Xp22.13, 18.6 Mb from the p-telomere.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK1222/>.

**This product is not CE/FDA registered for use in diagnostic procedures. The SALSA® MLPA® technique is covered by US patent 6,955,901 and corresponding patents outside the US. The purchase of this product includes a license to use only this amount of product solely for the purchaser's own use.**

#### Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <https://www.ncbi.nlm.nih.gov/gene>

For NM\_ mRNA reference sequences: <https://www.ncbi.nlm.nih.gov/nucleotide>

Matched Annotation from NCBI and EMBL-EBI (MANE): <https://www.ncbi.nlm.nih.gov/refseq/MANE>

#### Exon numbering

The *RS1* exon numbering used in this P488-B1 *RS1* product description is the exon numbering derived from MANE project based on MANE Select transcript NM\_000330.4. From description version B1-02 onwards, we have adopted the NCBI exon numbering that is present in the MANE Select transcript for this gene. The exon numbering used in previous versions of this product description can be found in between brackets in Table 1 and 2. As changes to the databases can occur after release of this product description, the NM\_ sequence and exon numbering may not be up-to-date.

#### Probemix content

P488-B1 *RS1* contains 19 MLPA probes with amplification products between 157 and 298 nucleotides (nt). This includes six probes for the *RS1* gene (one probe for each exon of the gene), and four probes (of which three downstream) targeting the *CDKL5* gene. In addition, nine reference probes are included that detect other locations on the X-chromosome. The probes targeting *CDKL5* could be informative on the size of the deletion or duplication. Complete probe sequences and the identity of the genes detected by the reference probes are available online ([www.mrcholland.com](http://www.mrcholland.com)).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at [www.mrcholland.com](http://www.mrcholland.com).

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

#### MLPA technique

The principles of MLPA (Schouten et al. 2002) are described in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

#### MLPA technique validation

Internal validation using 16 different DNA samples from healthy individuals of the same sex is required, in particular when using MLPA for the first time, or when changing the sample type or the sample handling

procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

### Required specimens

Extracted DNA, free from impurities known to affect MLPA reactions. MRC Holland has tested and can recommend the following extraction methods:

- QIAGEN Autopure LS (automated) and QIAamp DNA mini/midi/maxi kit (manual)
- Promega Wizard Genomic DNA Purification Kit (manual)
- Salting out (manual)

All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

### Reference samples

A sufficient number ( $\geq 3$ ) of different reference samples from unrelated individuals should be included in each MLPA experiment for data normalisation. Reference samples should be derived from individuals who are from families without a history of hereditary eye or vision abnormalities. It is recommended to use samples of the same sex to facilitate interpretation. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore deviations to the indicated copy number alteration (CNA) findings might occur.

### Data analysis

Coffalyser.Net should be used for data analysis in combination with the appropriate lot-specific Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net is freely downloadable at [www.mrcholland.com](http://www.mrcholland.com). Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

### Interpretation of results

The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$  and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results:

Copy Number status: Male samples	Final ratio
Normal	$0.80 < FR < 1.20$
Deletion	$FR = 0$
Duplication	$1.65 < FR < 2.25$
Ambiguous copy number	All other values

Copy Number status: Female samples	Final ratio
Normal	$0.80 < FR < 1.20$
Homozygous deletion	$FR = 0$
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of Coffalyser.Net (Calculations, cut-offs and interpretation remain unchanged.) Please note that Coffalyser.Net also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. single nucleotide variants (SNVs), point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination in the DNA sample) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <https://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

### Limitations of the procedure

- In most populations, the major cause of genetic defects in the *RS1* gene are small (point) mutations, none of which will be detected by using P488 RS1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

### Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a

heterozygous mutation or polymorphism in sequence data indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

### RS1 mutation database

<https://databases.lovd.nl/shared/genes/RS1>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <https://varnomen.hgvs.org>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *RS1* exons 4 and 6 but not exon 5) to MRC Holland: [info@mrcholland.com](mailto:info@mrcholland.com).

**Table 1. P488-B1 RS1**

Length (nt)	MLPA probe	Chromosomal position (hg18) <sup>a</sup>		
		Reference	<i>RS1</i>	<i>CDKL5</i>
64-105	Control fragments – see table in probemix content section for more information			
157	Reference probe 20449-L28192	Xp11		
166	Reference probe 04423-L05579	Xp22		
172	<b>RS1 probe</b> 21245-L29081		<b>Exon 1</b>	
178	Reference probe 05153-L04557	Xq22		
184 Ø	<b>CDKL5 probe</b> 21246-L29116			NM_001323289.2: <b>Downstream</b> (NM_003159.2: <b>Exon 20</b> )
190	Reference probe 07653-L07359	Xp11		
195	<b>RS1 probe</b> 21247-L29115		<b>Exon 5</b>	
211 Ø	<b>CDKL5 probe</b> 21255-L31466			<b>Exon 18</b>
220	<b>RS1 probe</b> 21248-L29086		<b>Exon 2</b>	
226	Reference probe 13519-L14318	Xq21		
232 Ø	<b>CDKL5 probe</b> 21257-L20734			NM_001323289.2: <b>Downstream</b> (NM_003159.2: <b>Exon 21</b> )
238	Reference probe 08248-L07945	Xq26		
246	<b>RS1 probe</b> 21249-L29085		<b>Exon 6</b>	
256	<b>RS1 probe</b> 21250-L29117		<b>Exon 3</b>	
267 Ø	<b>CDKL5 probe</b> 21251-L29084			NM_001323289.2: <b>Downstream</b> (NM_003159.2: <b>Exon 19</b> )
274	Reference probe 19689-L26419	Xq23		
283	<b>RS1 probe</b> 21254-L29087		<b>Exon 4</b>	
291	Reference probe 01370-L01287	Xp21		
298	Reference probe 05614-L05891	Xq28		

<sup>a</sup> See section Exon numbering on page 2 for more information.

Ø Probe detects the *CDKL5* gene, which partially overlaps with the *RS1* gene, according to alternative *CDKL5* transcript NM\_003159.2. *CDKL5* exons 19, 20 and 21 are situated in *RS1* intronic sequences. According to the MANE Select transcript NM\_001323289.2, these probes are downstream of the *CDKL5* gene.

The probe lengths in the table above may vary slightly depending on the capillary electrophoresis machine settings. Please see the most up-to-date Coffalyser sheet for exact probe lengths obtained at MRC Holland.

SNVs located in the target sequence of a probe can influence probe hybridisation and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

**Table 2. RS1/CDKL5 probes arranged according to chromosomal location**

Length (nt)	MLPA probe	RS1 exon <sup>a</sup>	Ligation site <sup>b</sup> NM_000330.4	Partial sequence <sup>c</sup> (24 nt adjacent to ligation site)	Distance to next probe
		start codon	41-43 (Exon 1)		
172	21245-L29081	Exon 1	2 nt before exon 1 reverse	TTACTGAACTGG-AGAGCTGGCCCA	14.5 kb
220	21248-L29086	Exon 2	105-106	CACATTGGGATT-ATCGTCTACCGA	0.9 kb
256	21250-L29117	Exon 3	173-174	ATTGCCAAGGAG-GACCCAATGCTC	3.1 kb
232 Ø	21257-L20734	CDKL5 gene; Downstream (Exon 21)	NM_001323289.2; 13.5 kb after exon 18 (NM_003159.23402-3403)	CTCATGGAAGAA-CCAATTAACACC	3.1 kb
184 Ø	21246-L29116	CDKL5 gene; Downstream (Exon 20)	NM_001323289.2; 10.4 kb after exon 18 reverse (NM_003159.2; 3148-3147 reverse)	ACCTGGAGTACT-GGATAGGGAGCA	3.1 kb
283	21254-L29087	Exon 4	61 nt before exon 4	TTGTTATTTTCAG-TCACCTGGTGCT	1.4 kb
267 Ø	21251-L29084	CDKL5 gene; Downstream (Exon 19)	NM_001323289.5; 9 kb after exon 18 reverse (NM_003159.2; 72 nt before exon 19 reverse)	AAACATATTCTA-TTGGCCACTTCT	1.3 kb
195	21247-L29115	Exon 5	44 nt before exon 5	CTTTTTCAGAC-ATGCACAGGGGG	3.1 kb
		stop codon	713-715 (Exon 6)		
246	21249-L29085	Exon 6	1095-1096	ACCACATTCAAG-ACAGAAACATAC	13.1 kb
211 Ø	21255-L31466	CDKL5 gene; Exon 18	NM_001323289.5; 2861-2862 (NM_003159.2; 2864-2865)	CCAAAAATTCCT-TCTCAGAAATTC	

<sup>a</sup> See section Exon numbering on page 2 for more information.

<sup>b</sup> Ligation sites are relative to the start of the NM\_ sequence, and not relative to the coding sequence.

<sup>c</sup> Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com). Please notify us of any mistakes: [info@mrcholland.com](mailto:info@mrcholland.com).

Ø Probe detects the *CDKL5* gene, which partially overlaps with the *RS1* gene, according to alternative *CDKL5* transcript NM\_003159.2. *CDKL5* exons 19, 20 and 21 are situated in *RS1* intronic sequences. According to the MANE Select transcript NM\_001323289.2, these probes are downstream of the *CDKL5* gene.

SNVs located in the target sequence of a probe can influence probe hybridisation and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

### Related products

For related products, see the [product page](#) on our website.

### References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Sikkink SK et al. (2007). X-linked retinoschisis: an update. *J Med Genet.* 44:225-232.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

P488 product history	
Version	Modification
B1	First release.

**Implemented changes in the product description**


Version B1-02 – 02 December 2025 (05P)

- Product description adapted to a new template.
- Gene structure and transcript variants section: link to MANE website added.
- Figure 1 removed.
- Table 1 and 2: updated the naming, exon numbering, ligation sites and NM for the CDKL5 probes according to MANE.
- Exon numbering on page 2: updated to MANE for RS1.

Version B1-01 – 31 January 2020 (02P)

- Not applicable, new document.

**More information: [www.mrcholland.com](http://www.mrcholland.com); [www.mrcholland.eu](http://www.mrcholland.eu)**

	MRC Holland BV; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	<a href="mailto:info@mrcholland.com">info@mrcholland.com</a> (information & technical questions) <a href="mailto:order@mrcholland.com">order@mrcholland.com</a> (orders)
Phone	+31 888 657 200

MRC Holland, SALSA, MLPA, digitalMLPA, Coffalyser.Net, Coffalyser digitalMLPA, and their logos are trademarks or registered trademarks of MRC Holland BV. All other brands and names herein are the property of their respective owners.