

Instructions for Use


SALSA® MLPA® Probemix P137 SCN1A



See also the MLPA General Protocol, the product description of the SALSA® MLPA® Reagent Kit and the Coffalyser.Net Reference Manual.

Visit the SALSA® MLPA® Probemix P137 SCN1A product page on our website to find Certificates of Analysis and a list of related products.

Product Name	SALSA® MLPA® Probemix P137 SCN1A
Version	C1
Catalogue numbers	P137-025R (25 reactions) P137-050R (50 reactions) P137-100R (100 reactions)
Basic UDI-DI	872021148P1375Y
Ingredients	Synthetic oligonucleotides, oligonucleotides purified from bacteria, Tris-HCl, EDTA

Regulatory Status	
IVD	EUROPE  2797
RUO	ALL OTHER COUNTRIES


Additional Test Components	Catalogue Numbers
SALSA® MLPA® Reagent Kit	EK1-FAM EK1-CY5 EK5-FAM EK5-CY5 EK20-FAM

Label Symbols			
IVD	In Vitro Diagnostic	RUO	Research Use Only

Storage and Shelf Life

Recommended conditions		
------------------------	--	--

A shelf life of until the expiry date is guaranteed, also after opening when stored in the original packaging under recommended conditions. For the exact expiry date, see the label on the vial. This product should not be exposed to more than 25 freeze-thaw cycles. Do not use the product if the packaging is damaged or opened. Leave chemicals in original containers. Waste material must be disposed of in accordance with the national and local regulations.

More Information:	
www.mrcholland.com	
	MRC Holland BV; Willem Schoutenstraat 1 1057 DL, Amsterdam, the Netherlands
E-mail	info@mrcholland.com (information & technical questions); order@mrcholland.com (orders)
Phone	+31 888 657 200

Any serious incident that has occurred in relation to this product should be reported to MRC Holland and the competent authority of the Member State or country in which the user and/or the patient is located.

Changes in this Product Version

As compared to version B3, probes targeting SCN1A exon 6 and 18, and two reference probes have been replaced in version C1.

1. Intended Purpose

The SALSA MLPA Probemix P137 SCN1A is an in vitro diagnostic (IVD)¹ or research use only (RUO) semi-quantitative manual assay² for the detection of deletions or duplications in the human *SCN1A* gene in genomic DNA isolated from human peripheral whole blood specimens. P137 SCN1A is intended to confirm a potential cause for and clinical diagnosis of Dravet syndrome (DS) and other *SCN1A*-related seizure disorders, and for molecular genetic testing of at-risk family members. In addition, results of this assay can aid in patient management decisions.

Copy number variations (CNVs) detected with P137 SCN1A should be confirmed with a different technique. In particular, CNVs detected by only a single probe always require confirmation by another method. Most defects in the *SCN1A* gene are point mutations, none of which will be detected by MLPA. It is therefore recommended to use this assay in combination with sequence analysis.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, parental evaluation, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes, pre-implantation or prenatal testing, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations.

¹ Please note that this probemix is for IVD use in the countries specified on page 1 of this product description. In all other countries, this is a RUO product.

² To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

2. Sample Requirements

Specimen	50-250 ng purified human genomic DNA, dissolved in 5 µl TE _{0.1} buffer, pH 8.0-8.5
Collection Method	Standard methods
Extraction Method	Methods tested by MRC Holland: <ul style="list-style-type: none"> • QIAGEN Autopure LS (automated) and QIAamp DNA mini/midi/maxi kit (manual) • Promega Wizard Genomic DNA Purification Kit (manual) • Salting out (manual)

Sample Types		
Test Sample	<ul style="list-style-type: none"> • Provided by user 	
Reference Samples (Required)	<ul style="list-style-type: none"> • Provided by user • Extraction method, tissue type, DNA concentration and treatment as similar as possible in all test and reference samples. • Have a normal copy number and ≤0.10 standard deviation for all probes. • At least three* independent reference samples required in each experiment for proper data normalisation. Derived from unrelated individuals from families without a history of DS or other <i>SCN1A</i>-related seizure disorders. 	
No-DNA Control (Preferably)	<ul style="list-style-type: none"> • Provided by user • TE_{0.1} buffer instead of DNA • To check for DNA contamination 	
Positive Control Samples (Preferably)	<ul style="list-style-type: none"> • Provided by user, or 	
	Available from third parties	See the table of positive samples on the probemix product page on our website.

*When testing >21 samples, include one extra reference for each 7 test samples.

3. Test Procedure

See the [MLPA General Protocol](#).

4. Quality Control, Data Analysis, and Troubleshooting

Quality Control Fragments in the Probemix	
Length (nt)	Function
64-70-76-82	DNA quantity control fragments
88-96	DNA denaturation control fragments
92	Benchmark fragment
100	Chromosome X presence control fragment
105	Chromosome Y presence control fragment

[Coffalyser.Net](#) should be used for data analysis in combination with the appropriate product and lot-specific Coffalyser sheet. See the [Coffalyser.Net Reference Manual](#) for details on data analysis and quality control.

For troubleshooting help, see the additional resources offered on our [support portal](#).

5. Interpretation of Results

Determining Typical Values in Normal and Affected Populations

The typical final ratio (FR) values stated in the copy number tables were determined in a validation study with samples containing abnormal copy numbers. The standard deviation of each individual probe over all the reference samples was ≤ 0.10 .

Expected Results of Reference Probes

Final Ratio (FR)	Copy Number	Description
0.80 – 1.20	2	Normal

Typical Results of Probes Targeting Two Copies (*SCN1A* gene)

Final Ratio (FR)	Copy Number	Description
0	0	Homozygous deletion
0.40 – 0.65	1	Heterozygous deletion
0.80 – 1.20	2	Normal
1.30 – 1.65	3	Heterozygous duplication
1.75 – 2.15	4	Homozygous duplication or Heterozygous triplication
All other values	-	Ambiguous

The tables illustrate the relationship between final probe ratio and corresponding copy number. Test results are expected to center around these values. Ambiguous values can indicate a technical problem, but may also reflect a biological cause such as mosaicism or a SNV influencing a single probe. It is important to use Coffalyser.Net to determine the significance of values found.

6. Performance Characteristics

Study	Description
Expected values for copy number in normal and affected populations	To determine the expected values in normal and affected populations a study was conducted on over 1500 MLPA reactions using samples with and without abnormal copy numbers. When the standard deviation of each individual probe over all the reference samples is ≤ 0.10 , the ranges stated in the copy number table in the product description can be used. Cut-off values for copy number determination were verified with SALSA MLPA Probemix P137 SCN1A in 45 samples from healthy individuals with normal copy number and 8 samples with known CNVs. The expected FRs for the corresponding copy number were found in all samples tested.
Limit of Detection	A study using representative probemixes was conducted to evaluate the minimum and maximum amount of DNA acceptable as the assay input. Results support the use of 50-250 ng of human DNA as the recommend input amount. The use of insufficient or too much sample DNA can affect performance. These lower and higher limits of detection were verified using SALSA MLPA Probemix P137 SCN1A on two samples with known CNVs and on one sample without any mutation and expected results were obtained using both the lower and upper input amount of DNA.
Interfering substances	SNVs or other polymorphisms (e.g. indels) in the DNA target sequence and impurities in the DNA sample (e.g. NaCl or KCl, EDTA and hemoglobin) can affect the MLPA reaction. A study using SALSA MLPA Probemix P137 SCN1A was performed to assess the potential for interference of endogenous and exogenous substances on genomic DNA on samples with known CNVs and one wild-type sample. For most probes, expected FRs (FRs within the expected cut-off category) were obtained even in the presence of potential interferents at concentrations shown in the table below.

Study	Description																											
	<table border="1" data-bbox="459 230 1436 600"> <thead> <tr> <th data-bbox="467 241 627 286">Interferent</th> <th data-bbox="635 241 834 286">Source</th> <th data-bbox="842 241 1010 286">Testing Concentration</th> <th data-bbox="1018 241 1436 286">Results*</th> </tr> </thead> <tbody> <tr> <td data-bbox="467 297 627 365">EDTA</td> <td data-bbox="635 297 834 365">Exogenous – specimen collection tubes</td> <td data-bbox="842 297 1010 365">1.5 mM</td> <td data-bbox="1018 297 1436 365">Copy number: Expected FR for 270/270 measurements</td> </tr> <tr> <td data-bbox="467 376 627 421">NaCl</td> <td data-bbox="635 376 834 421">Exogenous – DNA extraction</td> <td data-bbox="842 376 1010 421">40 mM</td> <td data-bbox="1018 376 1436 421">Copy number: Expected FR for 250/270 measurements</td> </tr> <tr> <td data-bbox="467 432 627 477">Fe3+ (FeCl3)</td> <td data-bbox="635 432 834 477">Exogenous – DNA extraction</td> <td data-bbox="842 432 1010 477">1 µM</td> <td data-bbox="1018 432 1436 477">Copy number: Expected FR for 270/270 measurements</td> </tr> <tr> <td data-bbox="467 488 627 555">Heparin</td> <td data-bbox="635 488 834 555">Exogenous – specimen collection tubes</td> <td data-bbox="842 488 1010 555">0.02 U/mL</td> <td data-bbox="1018 488 1436 555">Copy number: Expected FR for 270/270 measurements</td> </tr> <tr> <td data-bbox="467 566 627 600">Hemoglobin</td> <td data-bbox="635 566 834 600">Endogenous – blood sample</td> <td data-bbox="842 566 1010 600">0.02 µg/µl</td> <td data-bbox="1018 566 1436 600">Copy number: Expected FR for 270/270 measurements</td> </tr> </tbody> </table> <p data-bbox="459 600 1436 633">* Results are summarised for all 30 SCN1A probes across three samples tested in triplicate.</p> <p data-bbox="459 645 1436 779">NaCl had the largest effect on copy number determination, as FRs within an ambiguous or incorrect range were found in all samples. DNA extraction methods from blood remove NaCl and during testing of 45 samples extracted from blood, the expected final ratios were found. Therefore, it is only when NaCl is present in excess that deviating results for these probes can be found. Importantly, Coffalyser.Net issued warnings for the samples in which NaCl showed an effect, as well as lowered quality scores.</p> <p data-bbox="459 790 1436 869">To minimise variability across samples, all samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible.</p>				Interferent	Source	Testing Concentration	Results*	EDTA	Exogenous – specimen collection tubes	1.5 mM	Copy number: Expected FR for 270/270 measurements	NaCl	Exogenous – DNA extraction	40 mM	Copy number: Expected FR for 250/270 measurements	Fe3+ (FeCl3)	Exogenous – DNA extraction	1 µM	Copy number: Expected FR for 270/270 measurements	Heparin	Exogenous – specimen collection tubes	0.02 U/mL	Copy number: Expected FR for 270/270 measurements	Hemoglobin	Endogenous – blood sample	0.02 µg/µl	Copy number: Expected FR for 270/270 measurements
Interferent	Source	Testing Concentration	Results*																									
EDTA	Exogenous – specimen collection tubes	1.5 mM	Copy number: Expected FR for 270/270 measurements																									
NaCl	Exogenous – DNA extraction	40 mM	Copy number: Expected FR for 250/270 measurements																									
Fe3+ (FeCl3)	Exogenous – DNA extraction	1 µM	Copy number: Expected FR for 270/270 measurements																									
Heparin	Exogenous – specimen collection tubes	0.02 U/mL	Copy number: Expected FR for 270/270 measurements																									
Hemoglobin	Endogenous – blood sample	0.02 µg/µl	Copy number: Expected FR for 270/270 measurements																									
Cross-reactivity	Cross-reactivity is the potential for probes to bind to homologous regions (e.g. pseudogenes) or other cross-reactive sequences. Quality tests were carried out to determine whether probes are specific to their target sequence and all probes met the quality criteria for specificity.																											
Accuracy	Results of accuracy are derived from trueness and precision studies. For trueness, five previously genotyped samples were tested using SALSA MLPA Probemix P137 SCN1A and found to have the expected results. Assay precision was tested by repeatedly testing samples with known copy number status over multiple days, and by multiple operators. Results showed a correct call in 1288/1290 data points, leading to a precision of >99%.																											
Clinical validity*	Large deletions and duplications in the <i>SCN1A</i> gene explain 2 – 3% of the DS cases, and 0.3 – 2% of cases with other <i>SCN1A</i> -related seizure disorders. *(Based on a 2006-2025 literature review)																											

Summary of Safety and Performance (SSP)

The SSP is available in the European database on medical devices (Eudamed), <https://ec.europa.eu/tools/eudamed>, or upon request.

Content – Probe Details Sorted by Chromosomal Position

Chr. position	Target	Exon	Distance to next probe	Length (nt)	Probe number	Warnings
2q24.3	SCN1A	Exon 29 (26)	2.5 kb	445	04547-L18519	
2q24.3	SCN1A	Exon 28 (25)	1.8 kb	364	17380-L04745	
2q24.3	SCN1A	Exon 27 (24)	2.1 kb	319	04545-L18535	
2q24.3	SCN1A	Exon 26 (23)	1.6 kb	280	04544-L04900	
2q24.3	SCN1A	Exon 25 (22)	2.9 kb	243	04543-L18518	
2q24.3	SCN1A	Exon 24 (21)	7.2 kb	202	04542-L03931	
2q24.3	SCN1A	Exon 23 (20)	2.4 kb	172	04541-L03930	
2q24.3	SCN1A	Exon 22 (19)	1.7 kb	149	04540-L03929	#
2q24.3	SCN1A	Exon 21 (18)	1.9 kb	427	04539-L03928	
2q24.3	SCN1A	Exon 20 (17)	20.5 kb	400	04538-L03927	
2q24.3	SCN1A	Exon 19 (16)	1.9 kb	373	04537-L03926	
2q24.3	SCN1A	Exon 18 (15)	1.4 kb	328	22542-L31769	
2q24.3	SCN1A	Exon 17 (14)	1.9 kb	292	04535-L03924	
2q24.3	SCN1A	Exon 16 (13)	1.0 kb	250	04534-L30856	
2q24.3	SCN1A	Exon 15 (12)	1.5 kb	211	04533-L03922	
2q24.3	SCN1A	Exon 14 (11)	1.4 kb	178	04532-L03921	
2q24.3	SCN1A	Exon 13 (10)	1.5 kb	154	04531-L05030	
2q24.3	SCN1A	Exon 12 (9)	1.0 kb	337	07367-L07014	
2q24.3	SCN1A	Exon 11 (8)	1.2 kb	391	04530-L03919	
2q24.3	SCN1A	Exon 10 (7)	3.0 kb	355	04529-L03918	
2q24.3	SCN1A	Exon 9 (6)	1.0 kb	310	04528-L03917	#
2q24.3	SCN1A	Exon 8 (5)	1.7 kb	274	04527-L04899	
2q24.3	SCN1A	Exon 7 (4)	1.8 kb	229	04526-L03915	
2q24.3	SCN1A	Exon 6 (3)	2.2 kb	196	22541-L31721	#
2q24.3	SCN1A	Exon 5 (2)	14.9 kb	166	04524-L03913	
2q24.3	SCN1A	Exon 4 (1)	54.3 kb	160	04523-L18517	
2q24.3	SCN1A	Exon 1 (hB)	0.2 kb	265	15942-L18068	∅
2q24.3	SCN1A	Exon 1 (hB)	21.1 kb	409	15943-L18069	∅
2q24.3	SCN1A	Upstream (hA)	0.1 kb	183	15940-L18066	∅
2q24.3	SCN1A	Upstream (hA)		236	15941-L18067	∅
1q	Reference			436	14775-L16472	
3p	Reference			221	07223-L21127	
4q	Reference			256	10808-L27953	
6q	Reference			454	14954-L16687	
7q	Reference			142	09258-L11422	
8q	Reference			190	06743-L06347	
11q	Reference			301	08313-L08182	
12q	Reference			382	04278-L23577	
14q	Reference			418	11057-L11726	
20p	Reference			346	05982-L05407	

Probe lengths may vary slightly depending on capillary electrophoresis instrument settings. Please see the most up to date Coffalyser sheet for exact probe lengths obtained at MRC Holland.

The SCN1A exon numbers are derived from the MANE project and are based on the MANE Select transcript. For more information, see the probe sequences document available on the product page at www.mrcholland.com. The previous exon numbering using RefSeq transcript NM_06920.4 is disclosed between brackets.

Chromosomal bands are based on: hg18.

7. Precautions and Warnings

Probe warnings

- ∅ These probes target sequences outside of the known coding region. Copy number alterations of only these probes are of unknown clinical significance. The significance of deletions/duplications of these probes (previously referred to as exons hA and hB) is not clear as these exons are non-coding. However, heterozygous microdeletions in this region have been found in two DS patients, and a point mutation in the upstream region has been described in a patient with partial epilepsy with antecedent febrile seizures (Gao et al. 2017, Nakayama et al. 2010).
- # The specificity of these probes relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only these probes can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

Probemix-specific precautions

1. This product is not known to contain any harmful agents. Based on the concentrations present, none of the ingredients are hazardous as defined by the Hazard Communication Standard. **A Safety Data Sheet (SDS) is not required for this product:** none of the ingredients contain dangerous substances at concentrations requiring distribution of an SDS (as per Regulation (EC) No 1272/2008 [EU-GHS/CLP] and 1907/2006 [REACH] and amendments).
2. Sample or technical artefacts may appear as a (mosaic) copy number change of the whole/partial gene. Whole/partial gene deletions or duplications should therefore be confirmed by analysis of an independent DNA sample, to exclude false positive results.
3. Small changes (e.g. SNVs, small indels) in the sequence targeted by a probe can cause false positive results, even when >20 nt from the probe ligation site. Sequence changes can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a

probe oligonucleotide to the sample DNA. Deviations detected by this product should be confirmed, and single-probe deviations always require confirmation. Sequencing of the target region is recommended. Please contact MRC Holland for more information: info@mrcholland.com.

4. Copy number alterations of reference probes are unlikely to be related to the condition tested.
5. The *SCN1A* gene is located in a complicated 2q24 region, since several highly homologous genes are also present (*SCN2A*, *SCN3A*, *SCN7A* and *SCN9A*). In rare cases, apparent duplications might therefore be due to sequence changes in the other similar genes.
6. Mosaicism has been reported in individuals with DS (de Lange et al. 2018, Nakayama et al. 2018). Mosaic *SCN1A* cases obtained with the P137 *SCN1A* probemix must be confirmed by analysis of a second, independently collected DNA sample or by a different technique, in order to exclude a false positive mosaic result.
7. Large deletions are often detected with this probemix. The slope correction algorithm in Coffalyser.Net may confuse a large *SCN1A* deletion for sloping, which may lead to an incorrectly applied slope correction and ambiguous results for multiple probes. Please contact MRC Holland if you suspect that the slope correction has been incorrectly applied: info@mrcholland.com.

Technique-specific precautions

See the [MLPA General Protocol](#).

8. Limitations

Probemix-specific limitations

1. The clinical significance of the following findings is not yet clear/clearly established: deletions or duplications in exon 1 and/or the upstream region (Gao et al. 2017, Nakayama et al. 2010).
2. The significance of exons 1 to 3 deletions is not clear as these exons are non-coding and alternative transcript variants using other transcription start sites are known.
3. Probes for non-coding exons 2 and 3 are not included in this mix, therefore their association with *SCN1A*-related seizure disorders cannot be determined with this probemix.

Technique-specific limitations

See the [MLPA General Protocol](#).

9. References Cited in this IFU

1. de Lange IM et al. (2018). Mosaicism of de novo pathogenic *SCN1A* variants in epilepsy is a frequent phenomenon that correlates with variable phenotypes. *Epilepsia*. 59:690-703.
2. Gao QW et al. (2017). A Point Mutation in *SCN1A* 5' Genomic Region Decreases the Promoter Activity and Is Associated with Mild Epilepsy and Seizure Aggravation Induced by Antiepileptic Drug. *Mol Neurobiol*. 54:2428-2434.
3. Madia F et al. (2006). Cryptic chromosome deletions involving *SCN1A* in severe myoclonic epilepsy of infancy. *Neurology*. 67:1230-1235.
4. Marini C et al. (2009). *SCN1A* duplications and deletions detected in Dravet syndrome: implications for molecular diagnosis. *Epilepsia*. 50:1670-1678.
5. Nakayama T et al. (2010). Deletions of *SCN1A* 5' genomic region with promoter activity in Dravet syndrome. *Hum Mutat*. 31:820-829.
6. Nakayama T et al. (2018). Somatic mosaic deletions involving *SCN1A* cause Dravet syndrome. *Am J Med Genet A*. 176:657-662.
7. Wang JW et al. (2008). Microchromosomal deletions involving *SCN1A* and adjacent genes in severe myoclonic epilepsy in infancy. *Epilepsia*. 49:1528-1534.

Implemented changes in the product description

Version C1-04 – 13 January 2026 (03S)

- Intended purpose was updated. It was specified that the assay is manual, the intended population was extended to include molecular genetic testing of at-risk family members, and it was specified that the assay can aid in patient management decisions. In addition, parental evaluation was added. Familial hemiplegic migraine type 3 (FHM3) was removed from the intended purpose due to lack of clinical evidence.
- Exon numbering in "Changes in this Product Version" adjusted to adhere to MANE Select transcript exon numbering.
- Section 6. Performance Characteristics was updated with the addition of a summary table presenting performance data.
- Previous exon numbering was based on RefSeq transcript NM_06920.4. This is now mentioned in the legend below the content table.
- Precaution 7 was added to specify follow-up steps when signal sloping is present.
- "Upstream probe" was changed to "upstream region" in probemix-specific limitation 1.
- Probemix is now IVDR certified.

MRC Holland, SALSA, MLPA, digitalMLPA, Coffalyser.Net, Coffalyser digitalMLPA, and their logos are trademarks or registered trademarks of MRC Holland BV. All other brands and names herein are the property of their respective owners.