

# Product Description

## SALSA® MLPA® Probemix P225-E1 PTEN

To be used with the MLPA General Protocol.

### Version E1

For complete product history see page 11.

### Catalogue numbers:

- **P225-025R:** SALSA MLPA Probemix P225 PTEN, 25 reactions.
- **P225-050R:** SALSA MLPA Probemix P225 PTEN, 50 reactions.
- **P225-100R:** SALSA MLPA Probemix, P225 PTEN, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see [www.mrcholland.com](http://www.mrcholland.com)).

### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at [www.mrcholland.com](http://www.mrcholland.com).

### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: [www.mrcholland.com](http://www.mrcholland.com). It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

### Intended purpose

The SALSA MLPA Probemix P225 PTEN is an in vitro diagnostic (IVD)<sup>1</sup> or research use only (RUO) semi-quantitative assay<sup>2</sup> for the detection of deletions or duplications in *PTEN* in genomic DNA isolated from human peripheral whole blood specimens. P225 PTEN is intended to confirm a potential cause for and clinical diagnosis of PTEN Hamartoma Tumour Syndrome (PHTS) and for molecular genetic testing of at-risk family members. PHTS includes Cowden syndrome (CS), Bannayan-Riley-Ruvalcaba syndrome (BRRS), *PTEN*-related Proteus syndrome (PS), and Proteus-like syndrome (PLS). This probemix can also be used for the detection of deletions or duplications in the *PTEN* pseudogene (*PTENP1*) in a research setting.

Copy number variations (CNVs) detected with P225 PTEN should be confirmed with a different technique. In particular, CNVs detected by only a single probe always require confirmation by another method. Most defects in *PTEN* are point mutations, none of which will be detected by MLPA. It is therefore recommended to use this assay in combination with sequence analysis.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, parental evaluation, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes, pre-implantation or prenatal testing, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations. In a research setting this assay can be used on tumour tissue-derived DNA.

<sup>1</sup> Please note that this probemix is for In Vitro Diagnostic use (IVD) in the countries specified at the end of this product description. In all other countries, the product is for Research Use Only (RUO).

<sup>2</sup> To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

### Clinical background

Phosphatase and tensin homolog (*PTEN*) is a tumour suppressor gene that is mutated in a large number of cancers at high frequency. Defects in the *PTEN* gene are the main cause of PTEN Hamartoma Tumour

Syndrome (PHTS), which is a dominantly inherited cancer predisposition syndrome, characterized by multiple hamartomas in several areas of the body. PHTS includes Cowden syndrome (CS), Bannayan-Riley-Ruvalcaba syndrome (BRRS), *PTEN*-related Proteus syndrome (PS), and Proteus-like syndrome (GeneReviews: <https://www.ncbi.nlm.nih.gov/books/NBK1488/>).

Cowden syndrome (CS; OMIM #158350) is inherited in an autosomal dominant manner and comprises 85% of PHTS cases. The incidence of CS is estimated to be 1:200,000. Affected individuals usually have macrocephaly, trichilemmomas, and papillomatous papules, and present by their late 20s. The lifetime risk of developing breast cancer is 85%, with an average age of diagnosis between 38 and 46 years. The lifetime risk for thyroid cancer (usually follicular, rarely papillary, but never medullary thyroid cancer) is approximately 35%. The risk for endometrial cancer may approach 28%.

Bannayan-Riley-Ruvalcaba syndrome (BRSS) is inherited in an autosomal dominant manner. It is present at birth and is characterized by macrocephaly, intestinal hamartomatous polyposis, lipomas, pigmented macules of the glans penis, intellectual disability (50% of the cases) and development delay. The risk of developing cancer in BRRS patients with a *PTEN* pathogenic variant is similar to patients with CS.

*PTEN*-related Proteus syndrome (PS) is a highly variable, severe disorder characterized by progressive segmental or patchy overgrowth of diverse tissues of all germ layers, affecting the skeleton, skin, adipose tissue and central nervous systems. PS is a rare condition with an incidence of less than 1 in 1 million people worldwide and is associated with tumours, pulmonary complications, and deep vein thrombosis.

Proteus-like syndrome (PLS) describes individuals that do not meet the diagnostic criteria of Proteus syndrome, but share many of the characteristic signs and symptoms associated with this condition. Inheritance is autosomal dominant in those with a *PTEN* mutation.

*PTENP1* can regulate cellular levels of *PTEN* (via binding to mRNAs that target *PTEN*) and thereby suppresses cell growth. It has been shown that *PTENP1* is selectively lost in sporadic colon cancer (Poliseno et al. 2010). A specific *PTENP1* deletion (not *PTEN*) was also demonstrated in human melanoma (Poliseno et al. 2011). Furthermore, the importance of *PTENP1* as a tumour suppressor has been recently shown in head and neck squamous cell carcinoma (Liu et al. 2017). *PTENP1* copy number detection is of great importance in cancer research, however, the clinical validity of this gene is not yet fully established.

### Gene structure

The *PTEN* gene spans 108 kilobases (kb) on chromosome 10q23.31 and contains 9 exons. The *PTEN* LRG\_311 is available at [www.lrg-sequence.org](http://www.lrg-sequence.org) and is identical to GenBank NG\_007466.2.

The *PTENP1* gene spans 3.9 kilobases (kb) on chromosome 9p13.3 and contains 1 exon. The *PTEN* NCBI sequence is NR\_023917.1. No LRG sequence is available for this gene.

### Transcript variants

For *PTEN*, multiple transcript variants have been described. Transcript variant 1 encodes multiple isoforms due to the use of alternative translation initiation codons and is a reference standard in the RefSeqGene project. The *PTEN* transcript variant 1 (NM\_000314.8; 8515 nt; coding sequence 846-2057; <https://www.ncbi.nlm.nih.gov/gene/5728>) encodes three isoforms: PTEN-L (or PTENalpha), the longest isoform resulting from an upstream non-AUG (CUG) start codon; and two shorter isoforms resulting from downstream AUG start codons. PTEN, the most abundant isoform, is derived from the 5'-most AUG start codon. Transcript variant 2 (NM\_001304718.2; 8628 nt) has an additional exon in intron 5 in the 3' coding region.

For *PTENP1*, information for this long-noncoding RNA was obtained from the NR\_023917.1 sequence (3932 nt) (<https://www.ncbi.nlm.nih.gov/gene/11191>).

### Exon numbering

The *PTEN* exon numbering used in this P225-E1 *PTEN* product description is the exon numbering from the LRG\_311 sequence. The *PTENP1* exon numbering used is the exon numbering from the RefSeq NR\_023917.1.

The exon numbering of the NM\_ sequence or NR\_sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM\_ sequence or NR\_sequence and exon numbering may not be up-to-date.

### Probemix content

The SALSA MLPA Probemix P225-E1 PTEN contains 49 MLPA probes with amplification products between 130 and 496 nucleotides (nt). This includes 22 probes for the *PTEN* gene, at least one probe per exon, 10 flanking probes, and 2 probes for the pseudogene *PTENP1*. In addition, 15 reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online ([www.mrcholland.com](http://www.mrcholland.com)).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at [www.mrcholland.com](http://www.mrcholland.com).

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

**Probemix P225 is not suitable for methylation analysis:** From version E1 onwards, the P225 probemix is not suitable for methylation analysis since PTEN/KLLN methylation-specific probes have been removed. For methylation analysis of this region, MRC Holland provides the SALSA MS-MLPA probemix ME001 Tumour suppressor mix-1, which contains three PTEN/KLLN methylation-specific probes for the shared CpG-island for *KLLN* and *PTEN* genes.

### MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola (2012).

### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation of  $\leq 0.10$  for all probes over the experiment.

### Required specimens

Extracted DNA from human peripheral whole blood specimens, free from impurities known to affect MLPA reactions. In a research setting, this assay can be used on tumour tissue-derived DNA. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

### Reference samples

A sufficient number ( $\geq 3$ ) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a

PHTS history. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### Positive control DNA samples

MRC Holland cannot provide positive DNA samples. The inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID number NA20125 from the Coriell Institute has been tested with this P225-E1 probemix at MRC Holland and can be used as a positive control sample to detect a heterozygous duplication of the complete *PTEN* gene. Several flanking probes (from *BMPR1A* until *HTRA1*) are also affected. Sample ID number ACC-50 (OPM-2) from the Leibniz Institute DSMZ has been tested at MRC Holland and can be used as a positive control sample to detect a homozygous *PTEN* exon 3-7 deletion. The quality of cell lines can change; therefore samples should be validated before use.

### Performance characteristics

*PTEN* copy number changes explain 4-12% of CS and 11% of BRSS cases (GeneReviews <https://www.ncbi.nlm.nih.gov/books/NBK1488/>). It is unknown how many PLS and PS cases can be explained by copy number changes in *PTEN*, however, the association between *PTEN* mutations and PLS and PS is well established (Eng 2003, Loffeld et al. 2006, Orloff and Eng 2008, Smith et al. 2002, Zhou et al. 2001). Analytical performance for the detection of deletions/duplications in the *PTEN* gene is very high and can be considered >99% (based on a 2009-2020 literature review).

Analytical performance can be compromised by: SNVs or other polymorphisms in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.

### Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at [www.mrcholland.com](http://www.mrcholland.com). The use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

### Interpretation of results

The expected results of *PTEN* and *PTENP1* specific MLPA probes are allele copy numbers of 2 (normal), 1 (heterozygous deletion), 3 (heterozygous duplication), 0 (homozygous deletion) or 4 (homozygous duplication).

The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$ . When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	FR = 0
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders

as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

**Please note that these above mentioned final ratios are only valid for germline testing. Final ratios are affected both by percentage of tumour cells and by possible subclonality.**

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for the correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or *ITIH5*, *CELF2*, *ZNF25*, *PCDH15*, *ANXA7*, *LG11*, *SUFU*, *HTRA1* flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.
- Probemix P225 does not cover the additional exon of transcript variant 2 NM\_001304718.2 located in *PTEN* intron 5.

#### **P225 specific note:**

- In a research setting, when using DNA samples isolated from tumour tissues, reference probes are more prone to have deviating copy number results as compared to blood-derived germline samples. When regions targeted by reference probes are affected by copy number alterations, it can help to turn the slope correction off in Coffalyser.Net analysis to get the correct copy number interpretation on the target region.

#### **Limitations of the procedure**

- In most populations, the major cause of genetic defects in the *PTEN* gene are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P225 PTEN.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

- When used on tumour DNA (for research use only): MLPA analysis on tumour samples provides information on the *average* situation in the cells from which the DNA sample was purified. Gains or losses of genomic regions or genes may not be detected if the percentage of tumour cells is low. In addition, subclonality of the aberration affects the final ratio of the corresponding probe. Furthermore, there is always a possibility that one or more reference probes *do* show a copy number alteration in a patient sample.

### **Confirmation of results**

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained. SALSA MS-MLPA probemix ME001 contains three PTEN/KLLN methylation-specific probes for the shared CpG-island in the *KLLN* and *PTEN* genes.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

### **PTEN mutation databases**

Cosmic *PTEN* mutation database <https://cancer.sanger.ac.uk/cosmic/gene/analysis?ln=PTEN> and the LOVD *PTEN* mutation database <https://databases.lovd.nl/shared/genes/PTEN>. We strongly encourage users to deposit positive results in the Cosmic Database and in the LOVD Database. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *PTEN* exons 6 and 8 but not exon 7) to MRC Holland: [info@mrcholland.com](mailto:info@mrcholland.com).

**Table 1. SALSA MLPA Probemix P225-E1 PTEN**

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) <sup>a</sup>			
		Reference	PTEN	PTENP1	Flanking
64-105	Control fragments – see table in probemix content section for more information				
130 *	Reference probe 19616-L26704	4p13			
137 *	Reference probe 03797-L04594	21q22			
142 *	<b>PTEN probe</b> 21999-L30837		<b>Exon 4</b>		
148 *	Reference probe 14199-L23450	2q13			
155	<b>PTEN probe</b> 13690-L15159		<b>Exon 6</b>		
162 ¥ Δ	<b>PTEN probe</b> 07685-L31034		<b>Exon 8</b>		
168 * « ~	<b>SUFU probe</b> 21051-L31033				10q24.32
173 *	<b>PTEN probe</b> 22000-L30838		<b>Exon 2</b>		
178	<b>PTEN probe</b> 17314-L20922		<b>Exon 3</b>		
184 * Δ ~	<b>ZNF25 probe</b> 05760-L06666				<b>10p11.1</b>
190 *	Reference probe 10710-L31035	6p12			
196 ¥ ±	<b>PTEN probe</b> 06729-L31036		<b>Exon 2</b>		
201 ¥	<b>PTEN probe</b> 18254-L31140		<b>Exon 1</b>		
208 *	<b>PTEN probe</b> 22001-L30839		<b>Exon 4</b>		
214	<b>PTEN probe</b> 07686-L15591		<b>Exon 9</b>		
222 ¥	<b>PTEN probe</b> 17387-L30897		<b>Exon 1</b>		
229 *	Reference probe 00967-L31037	3q12			
239 ¥ ~	<b>PCDH15 probe</b> 08751-L22240				10q21.1
246 *	Reference probe 08715-L30393	9q21			
252 ~	<b>CELF2 probe</b> 17393-L22030				<b>10p14</b>
259 ~	<b>BMPR1A probe</b> 19351-L28591				10q23.2
266 *	Reference probe 07391-L30898	12q1			
274 ¥ ~	<b>LGI1 probe</b> 19294-L30901				10q23.33
280 *	<b>PTENP1 probe</b> 22040-L30965			<b>Exon 1</b>	
286	<b>PTEN probe</b> 17390-L14811		<b>Exon 2</b>		
292 *	Reference probe 16435-L30904	18q21			
299 ¥ ~	<b>RET probe</b> 18081-L30902				10q11.21
305 ~	<b>ANXA7 probe</b> 18380-L25185				10q22.2
312 *	Reference probe 13396-L30900	6q12			
319	<b>PTEN probe</b> 03639-L21321		<b>Exon 6</b>		
328	<b>PTEN probe</b> 19293-L25664		<b>Exon 3</b>		
337 ¥	<b>PTEN probe</b> 17396-L31245		<b>Exon 9</b>		
344	<b>PTEN probe</b> 18694-L24032		<b>Exon 3</b>		
352	Reference probe 05273-L25208	2p22			
359 ±	<b>PTEN probe</b> 17397-L25715		<b>Exon 8</b>		
369 ¥ ~	<b>HTRA1 probe</b> 08602-L30903				10q26.13
379	<b>PTEN probe</b> 03638-L24933		<b>Exon 5</b>		
391 *	Reference probe 08872-L30905	1p31			
400 *	<b>PTENP1 probe</b> 22042-L31246			<b>Exon 1</b>	
409	<b>PTEN probe</b> 13032-L22244		<b>Exon 5</b>		
418 *	Reference probe 11562-L12309	5q31			
427	Reference probe 08839-L22026	2p13			
436	<b>PTEN probe</b> 13692-L21061		<b>Exon 7</b>		
444	<b>PTEN probe</b> 17395-L21062		<b>Exon 9</b>		
454 ~	<b>ITIH5 probe</b> 17392-L21057				<b>10p14</b>
465	<b>PTEN probe</b> 17394-L21385		<b>Exon 1</b>		
475	<b>PTEN probe</b> 17386-L22174		<b>Exon 7</b>		
486	Reference probe 13594-L22376	19p13			
496 *	Reference probe 14894-L31209	15q21			

<sup>a</sup> See section Exon numbering on page 3 for more information.

\* New in version E1.

¥ Changed in version E1. Minor alteration, no change in sequence detected.

± SNP rs146326040 and rs562164491 (probe 196 nt and 359 nt, respectively) could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

– Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the *ITIH5*, *CELF2*, *ZNF25*, *PCDH15*, *ANXA7*, *LGI1*, *SUFU*, *HTRA1* flanking probes or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Please note: not all known SNVs are mentioned in the tables above. Single probe aberration(s) must be confirmed by another method.

**Table 2. P225 probes arranged according to chromosomal location**

Table 2a. *PTEN*

Length (nt)	SALSA MLPA probe	Gene/exon <sup>a</sup>	Ligation site / location	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
Centromeric flanking probes					
454 ~	17392-L21057	<i>ITIH5</i>	10p14	CGAGCAGAGTCA-TCTTGGATGGTG	3.6 Mb
252 ~	17393-L22030	<i>CELF2</i>	10p14	TCCCCCGGTCAT-GGTCGGAAAAGG	27.0 Mb
184 Δ ~	05760-L06666	<i>ZNF25</i>	10p11.1	ATGTTATTGTGG-AATTCACCAAGG	4.6 Mb
299 ~	18081-L30902	<i>RET</i>	10q11.21	TCCTCTACCTTA-ACCGGAGCCTGG	12.9 Mb
239 ~	08751-L22240	<i>PCDH15</i>	10q21.1	ACAACATGCTGA-TCAAAGGACTG	19.0 Mb
305 ~	18380-L25185	<i>ANXA7</i>	10q22.2	AGTCCCACCAGG-TGGAGCAGGCTT	13.8 Mb
259 ~	19351-L28591	<i>BMPR1A</i>	10q23.2	TACTTACATGCA-TGTGTTATTAAT	0.9 Mb
<b><i>PTEN</i> on 10q23.31. Ligation sites are according to NM_000314.8.</b>					
		<i>start codon</i>	846-848 (Exon 1)		
201	18254-L31140	Exon 1	14-13 reverse	GCTGCAGCTTCC-GAGAGGAGAGAA	0.8 kb
222	17387-L30897	Exon 1	781-782	CCTGCAGAAGAA-GCCCCGCCACCA	0.2 kb
465	17394-L21385	Exon 1	5 nt after Exon 1	TTGACCTGTATC-CATTTCTGCGGC	29.5 kb
286	17390-L14811	Exon 2	Intron 1 - 925	TAAAGTACTCAG-ATATTTATCCAA	0.1 kb
173	22000-L30838	Exon 2	998-997 reverse	CTTACTACATCA-TCAATATTGTTT	0.2 kb
196 ±	06729-L31036	Exon 2	217 nt after Exon 2 reverse	TATCACATAAGT-ACCTGATTATGT	31.0 kb
344	18694-L24032	Exon 3	145 nt before Exon 3	GGGGTATTTGTT-GGATTATTTATT	0.2 kb
328	19293-L25664	Exon 3	28 nt after Exon 3	ATTTGTATGCTT-GCAAATATCTTC	0.2 kb
178	17314-L20922	Exon 3	226 nt after Exon 3	TTGATCTGCTTT-AAATGACTTGGC	5.0 kb
142	21999-L30837	Exon 4	282 nt before Exon 4	AGCACCTGAATT-TACAGTACTCTG	2.3 kb
208	22001-L30839	Exon 4	1093-1094	CAAATTTAATTG-CAGAGGTAGGTA	0.1 kb
379 #	03638-L24933	Exon 5	1250-1251	GGTGTAAATGATA-TGTGCATATTTA	0.1 kb
409	13032-L22244	Exon 5	7 nt after Exon 5	AAAAGGTAAGTT-ATTTTTTGATGT	18.9 kb
155 #	13690-L15159	Exon 6	1389-1390	ATAGTACTCTGT-TAAAGAATCATC	0.1 kb
319	03639-L21321	Exon 6	1473-1472 reverse	CTTACTGCAAGT-TCCGCCACTGAA	5.7 kb
475 #	17386-L22174	Exon 7	1550-1551	ACACGACGGGAA-GACAAGTTCATG	0.1 kb
436	13692-L21061	Exon 7	4 nt after Exon 7	TAAAAAAGGTTT-GTACTTTACTTT	3.0 kb
359 # ±	17397-L25715	Exon 8	1787-1788	AATGACAAGGAA-TATCTAGTACTT	0.1 kb
162 # Δ	07685-L31034	Exon 8	45 nt after Exon 8	GACTTGTATGTA-TGTGATGTGTGT	4.1 kb
444	17395-L21062	Exon 9	6 nt before Exon 9	TAAATTTTCTTT-CTCTAGGTGAAG	0.3 kb
337 #	17396-L31245	Exon 9	2171-2170 reverse	AGAGAATTGTTT-CTATAACTGGTA	0.8 kb
214 #	07686-L15591	Exon 9	3003-3002 reverse	ACAGCATCTGAA-TTTTAGCACTGG	5.8 Mb
		<i>stop codon</i>	2055-2057 (Exon 9)		
Telomeric flanking probes					
274 ~	19294-L30901	<i>LGI1</i>	10q23.33	CTGAAATGGCTA-GTGAATGGCTT	8.8 Mb



168 « -	21051-L31033	<i>SUFU</i>	10q24.32	GAGGCCTGGTGA-GAAATGTGTGAT	19.9 Mb
369 -	08602-L30903	<i>HTRA1</i>	10q26.13	AATTGTTTCGCA-AGTAAAGAGAGC	

 Table 2b. *PTENP1*

Length (nt)	SALSA MLPA probe	<i>PTENP1</i> exon <sup>a</sup>	Ligation site NR_023917.1	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
400 #	22042-L31246	Exon 1	2109-2110	TTTCTCTTTTCG-TGACCAATCTTG	0.3 kb
280 #	22040-L30965	Exon 1	2423-2424	AGTTTGCAGTTA-GCTAAGAGAAGT	

<sup>a</sup> See section Exon numbering on page 3 for more information.

<sup>b</sup> Only partial probe sequences are shown. Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com). Please notify us of any mistakes: [info@mrcholland.com](mailto:info@mrcholland.com).

± SNP rs146326040 and rs562164491 (probe 196 nt and 359 nt, respectively) could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the *ITIH5*, *CELF2*, *ZNF25*, *PCDH15*, *ANXA7*, *LGI1*, *SUFU*, *HTRA1* flanking probes or reference probes are unlikely to be related to the condition tested.

# This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Please note: not all known SNVs are mentioned in the tables above. Single probe aberration(s) must be confirmed by another method.

Table 3. Reference probes arranged according to chromosomal location.

Length (nt)	SALSA MLPA probe	Gene	Chromosomal band (hg18)	Partial sequence (24 nt adjacent to ligation site)
391	08872-L30905	<i>LEPR</i>	1p31	AGCACATACTGT-TACGGTTCTGGC
427	08839-L22026	<i>DYSF</i>	2p13	TGCCATGAAGCT-GGTGAAGCCCTT
352	05273-L25208	<i>SPAST</i>	2p22	CGAGCCACAGCA-AAAAGAGCCCTC
148	14199-L23450	<i>EDAR</i>	2q13	GAGAGTTCTGTG-GGTGGAGAGAAG
229	00967-L31037	<i>SENP7</i>	3q12	CAGTCATTTGTT-TTCCATGGTTAG
130	19616-L26704	<i>ATP8A1</i>	4p13	CAGATTCTTCTT-CGAGGAGCTCAG
418	11562-L12309	<i>MYOT</i>	5q31	TCAGGATCTCAA-CAAGGAAGAGCA
190	10710-L31035	<i>PKHD1</i>	6p12	GGTTCCTGCTCT-TTCCAGTACCTC
312	13396-L30900	<i>EYS</i>	6q12	AAGGTTTGATGT-ACTCACCTACAA
246	08715-L30393	<i>PCSK5</i>	9q21	AAGCTGAGACCT-AGTTCCAGAGGG
266	07391-L30898	<i>COL2A1</i>	12q13	TGAACCTGGTGA-ACCTGGTGTCTC
496	14894-L31209	<i>SPG11</i>	15q21	GGACAATTCGCT-TTGGCCAGGAGG
292	16435-L30904	<i>MYO5B</i>	18q21	GAACAGCTCAAC-AACCAAATCCTG
486	13594-L22376	<i>CACNA1A</i>	19p13	ACTGGAGGAATG-GCAGCCCCTGGT
137	03797-L04594	<i>KCNJ6</i>	21q22	CTCGAAGCTCCT-ACATCACCAGTG

Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com).

### Related SALSA MLPA probemixes

ME001 Tumour suppressor mix 1 Tumour suppressor genes; including *KLLN/PTEN* shared promoter region.

P067 *PTCH1* Gorlin syndrome; gene included: *PTCH1*.

P081/P082 NF1 mix 1 & mix 2 Neurofibromatosis type 1; gene included: *NF1*.

P101 STK11	Peutz-Jeghers syndrome; gene included: <i>STK11</i> .
P158 JPS	Juvenile Polyposis; genes included: <i>BMPR1A</i> , <i>SMAD4</i> and <i>PTEN</i> .
P256 FLCN	Birt-Hogg-Dube syndrome; gene included: <i>FLCN</i> .
P472 SUFU	Familial Medulloblastoma and Meningioma; gene included: <i>SUFU</i> .

## References

- Eng C et al. (2003). PTEN: one gene, many syndromes. *Hum Mutat.* 22(3):183-98.
- Hömig-Hölzel C and Savola S. (2012). Multiplex ligation-dependent probe amplification (MLPA) in tumor diagnostics and prognostics. *Diagn Mol Pathol.* 21:189-206.
- Liu J et al. (2017). Decreased expression of pseudogene PTENP1 promotes malignant behaviours and is associated with the poor survival of patients with HNSCC. *Sci Rep.* 7:41179.
- Loffeld A et al. (2006). Epidermal naevus in Proteus syndrome showing loss of heterozygosity for an inherited PTEN mutation. *Br J Dermatol.* 154(6):1194-8.
- Orloff MS and Eng C. (2008). Genetic and phenotypic heterogeneity in the PTEN hamartoma tumour syndrome. *Oncogene.* 27(41):5387-97.
- Pilarski R et al. (2011). Predicting PTEN mutations: an evaluation of Cowden syndrome and Bannayan-Riley-Ruvalcaba syndrome clinical features. *J Med Genet.* 48:505-12.
- Poliseno L et al. (2010). A coding-independent function of gene and pseudogene mRNAs regulates tumour biology. *Nature.* 465:1033-8.
- Poliseno L et al. (2011). Deletion of PTENP1 pseudogene in human melanoma. *J Invest Dermatol.* 131:2497-2500.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Smith JM et al. (2002). Germline mutation of the tumour suppressor PTEN in Proteus syndrome. *J Med Genet.* 39(12):937-40.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.
- Zhou X et al. (2001). Association of germline mutation in the PTEN tumour suppressor gene and Proteus and Proteus-like syndromes. *Lancet.* 358(9277):210-1.

## Selected publications using SALSA MLPA Probemix P225 PTEN

- Ebrahimizadeh W et al. (2020). Design and development of a fully synthetic MLPA-based probe mix for detection of copy number alterations in prostate cancer formalin-fixed, paraffin-embedded tissue samples. *J Mol Diagn.* 22(10):1246-63.
- Ge NL et al. (2000). Expression of PTEN in PTEN-deficient multiple myeloma cells abolishes tumor growth in vivo. *Oncogene.* 19:4091-4095
- Hollander MC et al. (2011). PTEN loss in the continuum of common cancers, rare syndromes and mouse models. *Nat Rev Cancer.* 11(4):289-301.
- Jonosky JB et al. (2020). Large bowel obstruction from advanced prostate cancer in a 41-year-old: the role of genetic testing. *S Afr J Surg.* 58(3):163.
- Maiques O et al. (2014). FISH analysis of PTEN in endometrial carcinoma. Comparison with SNP arrays and MLPA. *Histopathology.* 65(3):371-88.
- Nathanson KL et al. (2013). Tumor genetic analyses of patients with metastatic melanoma treated with the BRAF inhibitor dabrafenib (GSK2118436). *Clin Cancer Res.* 19(17):4868-78.
- Sandell S et al. (2013). An intronic polymorphic deletion in the PTEN gene: implications for molecular diagnostic testing. *Br J Cancer.* 108(2):438-41.
- Wartenberg M et al. (2016). PTEN alterations of the stromal cells characterise an aggressive subpopulation of pancreatic cancer with enhanced metastatic potential. *Eur J Cancer.* 65:80-90.

- Yehia LE et al. (2020). The Clinical Spectrum of PTEN Mutations. *Annu Rev Med.* 71:103-116.
- Zhou XP et al. (2003). Germline PTEN promoter mutations and deletions in Cowden/Bannayan-Riley-Ruvalcaba syndrome result in aberrant PTEN protein and dysregulation of the phosphoinositol-3-kinase/Akt pathway. *Am J Hum Genet.* 73(2):404-11.

P225 product history	
Version	Modification
E1	All <i>KLLN</i> and <i>HhaI</i> digestion control probes have been removed, thus methylation detection is no longer possible in the shared promoter region of <i>PTEN/KLLN</i> genes by MS-MLPA with this probemix. In addition, <i>PTEN</i> exon 4 probes have been replaced and one probe has been added for exon 2; two <i>PTEN</i> -flanking and two <i>PTENP1</i> probes have been replaced; ten reference probes have been replaced and one has been added. Nine probes have a modification in length, not in the targeted sequence.
D2	One probe has a small change in length but no change in the sequence detected.
D1	One probe for <i>PTEN</i> exon 3 and one flanking probe for <i>PTEN</i> have been added and four flanking probes have been replaced. An <i>HhaI</i> digestion control probe has been included and several reference probes have been replaced.
C1	Several new target and reference probes have been added. P225 may also be used to detect methylation of <i>PTEN</i> and <i>KLLN</i> .
B3	The 88 and 96 nt control fragments have been replaced (QDX2).
B2	The number of <i>PTEN</i> probes has been increased to 25.
B1	Eight new <i>PTEN</i> probes have been added and three variable <i>PTEN</i> probes have been removed. In order to facilitate analysis of tumour DNA, several chromosome 10 probes have been included.
A1	Extra control fragments at 88, 96, 100 and 105 nt have been added.
A0	First release.

Implemented changes in the product description
<p>Version E1-04 – 16 June 2021 (04P)</p> <ul style="list-style-type: none"> <li>- Product description rewritten and adapted to a new template.</li> <li>- Various minor textual or layout changes.</li> <li>- Various minor textual changes in the Performance characteristics.</li> <li>- One point from the Interpretation of results section updated to state that copy number changes detected by reference probes or <i>ITIH5</i>, <i>CELF2</i>, <i>ZNF25</i>, <i>PCDH15</i>, <i>ANXA7</i>, <i>LGI1</i>, <i>SUFU</i>, <i>HTRA1</i> flanking probes are unlikely to have any relation to the condition tested for.</li> <li>- Ligation sites of the probes targeting the <i>PTEN</i> gene updated according to new version of the NM_ reference sequence.</li> <li>- Warning added to Table 1 and Table 2a for probes at 162 nt (07685-L31034) and 184 nt (05760-L06666) being more variable.</li> <li>- Warning added to Table 1 and Table 2a for probe at 168 nt (21051-L31033) located in or near a GC-rich region.</li> <li>- Warning for the flanking probes in Table 1 and Table 2a updated: Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the <i>ITIH5</i>, <i>CELF2</i>, <i>ZNF25</i>, <i>PCDH15</i>, <i>ANXA7</i>, <i>LGI1</i>, <i>SUFU</i>, <i>HTRA1</i> flanking probes or reference probes are unlikely to be related to the condition tested.</li> <li>- Updated all references in the Selected Reference section.</li> <li>- UK has been added to the list of countries in Europe that accept the CE mark.</li> </ul> <p>Version E1-03 – 29 April 2020 (04)</p> <ul style="list-style-type: none"> <li>- Product is now registered for IVD use in Israel.</li> <li>- Version E1-02 – 10 April 2019 (04)</li> <li>- Corrected <i>PTEN</i> exon number in Table 1 for <i>PTEN</i> probe 18254-L31140.</li> <li>- Changed the number of <i>PTEN/KLLN</i> probes in ME001 probemix to three in “Probemix P225 not suitable for methylation analysis” on page 3 and “Confirmation of results paragraph” on page 5.</li> </ul>

<p>Version E1-01 – 19 October 2018 (04)</p> <ul style="list-style-type: none"> <li>- Product description restructured and adapted to a new template.</li> <li>- Several references removed and added.</li> <li>- Warning added to Table 1 and 2a for SNPs that could influence probe signal: 196 nt probe 06729-L31036, and 359 nt probe 17397-L25715.</li> <li>- Warning added to Table 2a and 2b for probe specificity relying on a single nucleotide difference between target gene and related gene or pseudogene.</li> </ul> <p>Version 23 – 15 February 2018 (T08)</p> <ul style="list-style-type: none"> <li>- Warning added in Tables 1 and 2 concerning a single nucleotide difference between the PTEN gene and the PTEN pseudogene (PTENP1) for several probes.</li> <li>- Warning added in Tables 1 and 2 for a frequent SNP (rs12000677) in probe 17308-L25183</li> </ul> <p>Version 22 – 9 September 2017 (T08)</p> <ul style="list-style-type: none"> <li>- Warning added in Table 1, 172 nt probe 17388-L08261, 227 nt probe 17387-L26030, 364 nt probe 18110-L25928, and 465 nt probe 17394-L21385.</li> <li>- Small changes of S0750-L21493 probe length in Table 1 and 2 in order to better reflect the true length of the amplification product.</li> </ul> <p>Version 21 – 21 March 2017 (T08)</p> <ul style="list-style-type: none"> <li>- Various minor textual changes.</li> </ul> <p>Version 20 – 02 December 2016 (T08)</p> <ul style="list-style-type: none"> <li>- Warning regarding HhaI enzymes that are resistant to heat inactivation added under Methylation-specific MLPA section.</li> <li>- Various minor textual changes.</li> </ul> <p>Version 19 – 17 June 2016 (T08)</p> <ul style="list-style-type: none"> <li>- New references added on page 2.</li> <li>- Warning added in Table 1 and Table 2, 208 nt probe 17391-L21278.</li> </ul>
--

<b>More information:</b> <a href="http://www.mrcholland.com">www.mrcholland.com</a> ; <a href="http://www.mrcholland.eu">www.mrcholland.eu</a>	
	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	<a href="mailto:info@mrcholland.com">info@mrcholland.com</a> (information & technical questions) <a href="mailto:order@mrcholland.com">order@mrcholland.com</a> (orders)
Phone	+31 888 657 200

	EUROPE*  ISRAEL
	ALL OTHER COUNTRIES

\*comprising EU (candidate) member states, members of the European Free Trade Association (EFTA), and the UK. The product is for RUO in all other European countries.